i-walker
i-rollator
i-ambulateur

item #s:

RTL10555BL

RTL10555RD

RTL10555SL

TM
drive

www.drivemedical.com

parts list

Item # RTL10555SL

WEIGHT CAPACITY: 300 LBS (136 KG)

- | | |
|-----------------------------------|-----------------|
| 1. Handgrip | 7. Seat |
| 2. Adjustable tubing handgrip | 8. Bag |
| 3. Rear reflector on handgrip | 9. Front wheel |
| 4. Brake lever | 10. Rear wheel |
| 5. Interlock button for handgrips | 11. Cane Holder |
| 6. Backrest | |

introduction

The Walking aid is suitable for a single user with limited mobility. The I-Walker is ideal for indoor and outdoor use.

includes:

- | | | |
|------------|--------------|---------------------------|
| 1. Wheels | 3. Carry Bag | 5. Cane Holder |
| 2. Handles | 4. Backrest | 6. Operating Instructions |

assembly instructions

1. Remove contents from carton. (Figure 1)


Figure 1

2. Insert front wheel (A) into receptacle located on the front of the frame. A click can be heard when the wheel is securely in place. Repeat for the second front wheel. (Figure 2)


Figure 2

3. Insert back wheel (B) into back of frame by depressing push pin and sliding into the back receptacle, making sure the wheels face outwards from the frame. The push pin will snap into place when secure. Repeat for the second back wheel. (Figure 3)


Figure 3

4. Unfold the I-Walker by pressing down on both sides of the seat frame. (Figure 4)

Figure 4

5. Insert the handles (C) into the handle receptacles. Press in the red buttons located on the back of the I-Walker and slide the handles into the frame to secure the handles in place. (Figure 5)

Figure 5

6. Attach the bag (D) to the seat frame by rotating and removing the knobs located on the front of the seat frame and sliding the bag into the grooves located on the inside of the frame. Replace the knobs to secure in place. (Figure 6)

NOTE: Fold the I-Walker slightly to make it easier to install the bag.

Figure 6

7. Attach the backrest (E) by sliding the backrest hooks down into the receptacles located on the front of the frame. Secure in place with second thumb screw. (Figure 7)

Figure 7

8. Attach the cane holder (F) to the handle by sliding it down over the receptacle located on the front of the handle. (Figure 8)

Figure 8

height adjustment

9. Press the red release buttons (G) for handgrips, as shown in Figure 9, and pull the handgrip upwards at the same time. As the handle moves, release the button. This automatically locks in the next possible height. Repeat this process until the desired height is reached.

Figure 9

braking / locking

10. Pull both brake levers upwards to apply the brakes while the I-Walker is moving. (Figure 10) Always brake slowly and steady.

Figure 10

Press both brake levers down to apply the locks while the I-Walker is stationary. (Figure 11) Pull up the brake levers and release to release the locks.

Figure 11

The effect of the brakes can be affected by abrasion of the wheels.

Regularly test the brakes.

sitting

11. Position the I-Walker on an even and solid surface. Lock the brakes (parking position). The I-Walker must stand and must not slip. Sit straight up and do not lean backwards.

DO NOT MOVE ROLLATOR WHILE SEATED.

folding / transport / storage

12. To fold, pull seat upwards by the red handle on the seat plate as shown in Figure 12. The I-Walker can now be folded completely and secured by using the safety clasp (G) as shown in Figure 13.

Figure 12

Figure 13

safety instructions

- Use only as a walking aid!
- Use only on even and solid surfaces.
- Maximum user weight: 300 lbs.
- Regularly check the screws and fasteners on the I-Walker and retighten if necessary.
- Use bag only for transport of objects and do not overload. Maximum load: 11 lbs.
- Do not expose the product to temperatures above 104°F for an extended period of time, as this can cause the I-Walker parts to overheat and burn the user.
- Please make sure the Rollator is fully opened, seat is down and the folding bar under seat is locked in the down position.
- Be sure that both knobs are fully tightened and are at equal heights. Test Rollator for stability before applying full weight to hand grips.
- Test the brake lever before operating the Rollator. When wheels are locked, lock bar on rear wheels should prevent wheels from turning. When locks are released, the wheels should turn freely.
- Always lock the brake lever before sitting and do not use seat when unit is on an incline or uneven ground.
- Replace tires when significantly worn or damaged.
- **DO NOT MOVE ROLLATOR WHILE SEATED.** This is a walking aid only and is not to be used as a transportation device.

cleaning / maintenance

Clean regularly by hand with a damp cloth. Do not use aggressive detergents.

lifetime limited warranty

The frame of your Drive Medical product is warrantied to be free of defects in material and workmanship for the lifetime of the product for the original consumer purchaser. The brake handle and housing are covered for a period of 5 years from the date of purchase for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.

© 2016 Medical Depot, Inc. All rights reserved.

Drive is a trademark of Medical Depot, Inc.

Port Washington N.Y. 11050 USA Made in China

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© 2016 Medical Depot, Inc. Tous droits réservés.
Drive est une marque de commerce de Medical Depot, Inc.