

Delta™ Ultra- Light 1000, Semi-Electric Bed



FIGURE A



FIGURE B



FIGURE C



FIGURE D



INSET E

- The Delta 1000 is truly universal because the bed frame can be used with Drive's and most other manufacturers' bed ends, old and new
- The transition box mounted on the foot section is key to interchangeability (Figure A)
- The transition box allows you to change the rotation of the high/low shaft
- Headboard is taller than footboard to better fit home decor (Figure B)
- The wood grain panels are break and scratch resistant
- Foot section weighs 50% less than conventional foot sections to make deliveries safe and easy
- Once installed, the high/low shaft will stay with the bed to eliminate loss or forgetting on future installations (Figure C)
- Labels and color-coded springs assure accurate installation of side rails for patient safety (Figure D)
- Unique motor is completely self-contained to reduce weight and noise
- The motor assembly can be installed and removed with patient in bed
- Bed, side rails and mattress safety-tested by SGS, a nationally recognized testing facility
- Satisfies all applicable FDA entrapment guidelines
- The motor is UL approved
- In the event of a power failure, a 9 volt battery, located in the motor, can lower the head and foot sections nine times
- No crank needed
- New and improved hand pendant has large, easy-to-use controls (Figure E)
- Bed ships in two cartons
- Channel frame construction provides superior strength and reduced weight
- Zinc-coated spring deck
- Extension kit available (model #15030EXTKIT) extends length from 80" to 84"
- Patent Pending
- Weight Capacity 450 lbs

PRODUCTS	Description
15030	Ultra-Light Semi-Electric Bed, 1/cs
15030BV-FR	Ultra-Light Semi-Electric Bed with Full-Length Side Rails (item #15001ABV), 1/cs
15030BV-HR	Ultra-Light Semi-Electric Bed with Half Rails (item #15201BV), 1/cs
15030BV-PKG	Ultra-Light Semi-Electric Bed with Full-Length Side Rails (item #15001ABV) and 80" Innerspring Mattress (item #15006), 1/cs
15030BV-PKG-T	Ultra-Light Semi-Electric Bed with Full-Length Side Rails (item #15001ABV) and 80" Therapeutic Support Mattress (item #15019), 1/cs
15030BV-PKG-1	Ultra-Light Semi-Electric Bed with Half Rails (item #15201BV) and 80" Innerspring Mattress (item #15006), 1/cs
15030BV-PKG-1-T	Ultra-Light Semi-Electric Bed with Half Rails (item #15201BV) and 80" Therapeutic Support Mattress (item #15019), 1/cs
15030BV-PKG-2	Ultra-Light Semi-Electric Bed with Full-Length Side Rails (item #15001ABV) and 80" Foam Mattress (item #15007), 1/cs
15096-SEMI	Ultra-Light Semi-Electric Bed Reclamation Kit, 1/ea

SPECIFICATIONS	Description
Bed Height Minimum Without Casters	12.75"
Bed Height Maximum With Casters	24"
Unit Dimensions	36" (W) x 88" (L)
Weight	163 lbs.
Carton Shipping Weight	179 lbs.
Weight Limit	450 lbs Maximum weight limit includes weight of occupant, all accessories such as mattress, side/assist rails, trapeze, any personal affects, etc.
Warranty	(Original Purchaser) 5 Years on Motor and Hand Control, Limited Lifetime on Welds and Frame, 1 Year on All Other Parts and Components

ACCESSORIES	Description
15030HBC	Bedding in a Box - contains 1 pillow case, 1 fitted sheet, 1 cover sheet, and 1 blanket., 1/ea, 5/cs
15030HBL	Fitted Sheet - Contains 2 Fitted Sheets, 2 sheets per bx/ 10 bxs per cs

PARTS	Description
15544H	
15544F	
15541	
15548	
15544HPF	1/ea
15544FPF	1/ea
15544HGB	1/ea
15544FGB	1/ea
15044CNL	1/ea
15044CL	1/ea
15044CSET	(2 of ea)
15030P-1001	1/ea
15038MO	1/ea
15030P-1006	1/bx
15030-7RMM	1/ea
15030HC-G	1/ea
15030PN-1035	1/ea
15030PN-1036	1/ea
15544HP-1	1/ea
15044P	1/ea
15030UBS	1/ea
1503C	1/ea
15544C	1/ea
15005P1018	1/ea
15030P-1008	1/ea
15030MC	1/ea
15005P-1040	1/bx
15544GBD	1/ea
15544GBMB	1/ea
15030-FRLINK	1/ea
15030-MS	2/pkg
15544F-CH	1/ea
15030	1/cs