

**KARMAN HEALTHCARE**

19255 San Jose Ave.

City Of Industry, CA 91748

Tel: 626-581-2235 Fax: 626-581-2335

Toll Free: 800-805-2762



QR Reader Smartphone Info

**FLEXX WHEELCHAIR ORDER FORM**Submitting for ☐ QUOTE ☐ ORDER**Mark For:**

Date: _____

Dealer Acct#: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____

Dealer Phone: (____) _____

Confirmation Email: _____

ST: _____

ZIP: _____

FAX: (____) _____

Purchase Order #: _____

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

Ship To Phone: _____

ST: _____

ZIP: _____

The **HCPCS CODES** are not intended to be, nor should be considered billing or legal advice. They are provided based on code verification through PDAC or our interpretations of billing. The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, please verify the appropriate payer's coverage policy

FLEXX Wheelchair:

FLEXX Ultra lightweight Adjustable Wheelchair standard features: swing away removable footrest, fixed backrest, quick release 24"x1" solid poly-urethane rear wheels, 8"x1" solid front casters, adjustable front forks, weighs 28.5 lbs, flip back height adjustable armrest, weight capacity 286 lbs.

BASE FRAME				
✓	SKU#	HCPCS	Description	Price
☐	KM8522Q	K0005 pending	"FLEXX" Wheelchair Diamond Black Color	\$1899.00

SEAT WIDTH				
☐	16B		16" wide	Std
☐	18B		18" wide	Std
☐	20B	E2201	20" wide	\$300.00

SEAT DEPTH				
☐	D16		16" depth (not available for 20" seat width)	Std
☐	D18		18" depth	Std

FRONT SEAT TO FLOOR HEIGHT				
☐	FSH15	K0056	15" seat to floor height	\$100.00
☐	FSH16	K0056	16" seat to floor height	\$100.00
☐	FSH17		17" seat to floor height	Std
☐	FSH18		18" seat to floor height	Std
☐	FSH19		19" seat to floor height	Std

FRONT SEAT TO FLOOR HEIGHT cont.				
✓	SKU#	HCPCS	Description	Price
☐	FSH20		20" seat to floor height	Std
☐	FSH21		21" seat to floor height	Std

REAR SEAT TO FLOOR HEIGHT*Maximum 2" difference from front seat to floor height*				
☐	RSH15	K0056	15" seat to floor height	\$100.00
☐	RSH16	K0056	16" seat to floor height	\$100.00
☐	RSH17		17" seat to floor height	Std
☐	RSH18		18" seat to floor height	Std
☐	RSH19		19" seat to floor height	Std
☐	RSH20		20" seat to floor height	Std

SEAT UPHOLSTERY				
☐	US-8500		Non-Adjustable	Std
☐	OMITUS-8500		Omit Seat Upholstery	Std

FRONT CASTER SELECTION (pair)				
☐	C51LT		5"x1" composite hub poly urethane flat free	Std
☐	C5EPIC		5"x1" aluminum hub urethane narrow court caster	\$90.00
☐	C61LT		6"x1" composite hub poly urethane flat free	Std
☐	C6EPIC		6"x1" aluminum hub urethane narrow court caster	\$108.00
☐	C6EPICW		6"x1.4" aluminum hub wide soft roll caster	\$132.00
☐	C71LT		7"x1" composite hub poly urethane flat free	Std
☐	C81LT		8"x1" composite hub poly urethane flat free	Std

Caster Size/Front Seat to Floor Selection				
Front Seat to Floor	Caster Size			
	5"	6"	7"	8"
15"	X			
16"	X	X		
17"	X	X		
18"	X	X	X	
19"		X	X	X
20"				X
21"				X

ARM OPTIONS				
☐	NAARMS		No armrest, No arm pad	Std
☐	AR8500SAP		Swing away padded arms with round tube arm pads	N/A
☐	AR8500HASP	E0973	Flip back height adjustable single post arms with ergonomic arm pads	Std

BACK POST				
☐	FBH8500-17		17" w/ergonomic push handles, 8 degree bend	Std
☐	FBH8500-20		20" w/ergonomic push handles, 8 degree bend	Std
☐	ABH8500-1521		15"-21" range, w/ergonomic push handles, 8 degree bend	\$165.00

BACK UPHOLSTERY				
✓	SKU#	HCPCS	Description	Price
☐	UB-8500-FLEXX		Slip on padded back upholstery	Std
☐	OMITUB-8500		Omit back upholstery	Std
☐	ADJTENS8500		Adjustable Tension back upholstery	N/A

CENTER OF GRAVITY				
☐	COGADJ1		Rear wheel center of gravity adjustment: Stable	Std
☐	COGADJ0		Rear wheel center of gravity adjustment: Normal	Std
☐	COGADJ-1		Rear wheel center of gravity adjustment: Active	Std

REAR WHEELS (pair)				
☐	RW-8500Q		24x1 quick release flat free polyurethane wheel assem. with Ergonomic handrim	Std
☐	RW88		24x1 MAG 8 spoke flat free polyurethane wheel assembly	N/A
☐	SPINLX		24x1 Spinergy LX 12 spoke, Pneumatic Kenda Tire and Tube w/ anodized polished handrim	\$895.00
☐	SPINSPOX		24x1 Spinergy SPOX 18 spoke, Pneumatic Kenda Tire and Tube w/anodized polished handrim	\$695.00
☐	SPINWIRE		24x1 Spinergy WIRE 30 spoke, Pneumatic Kenda Tire and Tube w/anodized polished handrim	\$295.00

KENDA TIRE COLOR/SPINERGY SPOKE COLOR (Selection applies only to Spinergy Wheel)				
☐	KENDAYLW		Yellow Kenda tire	Std
☐	KENDARED		Red Kenda tire	Std
☐	KENDABLU		Blue Kenda tire	Std
☐	SPINSPYLW		Yellow Spinergy Spoke color	Std
☐	SPINSPRED		Red Spinergy Spoke color	Std
☐	SPINSPBLU		Blue Spinergy Spoke color	Std
☐	SPINSPBLK		Black Spinergy Spoke color	Std

BACKREST ANGLE				
☐	BAADJ-4		Back angle adjustment -4 degrees	Std
☐	BAADJ0		Back angle adjustment 0 degrees	Std
☐	BAADJ4		Back angle adjustment 4 degrees	Std
☐	BAADJ8		Back angle adjustment 8 degrees	Std
☐	BAADJ12		Back angle adjustment 12 degrees	Std
☐	BAADJ16		Back angle adjustment 16 degrees	Std

HANGAR OPTIONS				
☐	FR8500Q-16		16" 70 degree Swing in/out non tapered footrest (pair)	Std
☐	FR8500Q-18		18" 70 degree Swing in/out non tapered footrest (pair)	Std
☐	FR8500Q-20		20" 70 degree Swing in/out non tapered footrest (pair)	Std
☐	EL8500Q-16		16" Articulating elevating leg rest with calf support (pair)	\$298.00
☐	EL8500Q-18		18" Articulating elevating leg rest with calf support (pair)	\$298.00
☐	EL8500Q-20		20" Articulating elevating leg rest with calf support (pair)	\$298.00

SEAT CUSHIONS				
✓	SKU#	HCPCS	Description	Price
☐	CU-ERGO-16	E2601	16x17 seat memory foam 2" cushion	\$196.00
☐	CU-ERGO-18	E2601	18x17 seat memory foam 2" cushion	\$176.00
☐	CU-ERGO-20	E2601	20x17 seat memory foam 2" cushion	\$216.00
☐	CU-FT-16	E2601	16x16 seat foam 2" cushion	\$72.00
☐	CU-FT-18	E2601	18x16 seat foam 2" cushion	\$72.00
☐	CU-FT-20	E2601	20x16 seat foam 2" cushion	\$78.00
☐	CU-GFR-16	E2601	16x16 seat gel 2" cushion	\$98.00
☐	CU-GFR-18	E2601	18x16 seat gel cushion	\$98.00
☐	CU-GFR-20	E2601	20x16 seat gel 2" cushion	\$108.00

BACK CUSHIONS				
☐	CU-ERGO-16-back	E2611	16" wide memory foam 2" back cushion	\$156.00
☐	CU-ERGO-18-back	E2611	18" wide memory foam 2" back cushion	\$156.00
☐	CU-ERGO-20-back	E2611	20" wide memory foam 2" back cushion	\$176.00
☐	BKF-1617	E2611	16" wide back foam cushion	\$98.00
☐	BKF-1817	E2611	18" wide back foam cushion	\$98.00
☐	BKF-2019	E2611	20" wide back foam cushion	\$126.00

OPTIONS & ACCESSORIES				
☐	SB22	E0978	Clamp style seat belt 2 piece	\$19.00
☐	SB88	E0978	Wrap around seat belt 1 piece	N/A
☐	SB99-48	E0978	Auto Style seat belt 2 piece	\$24.00
☐	BT-115	E0961	Brake tip-Extending	Std
☐	HR-FLD-115	E0955 + E1028	16-18" width Foldable Headrest with hardware	\$327.00
☐	HR-FLD-115W	E0955 + E1028	20-22" width Foldable Headrest with hardware	\$357.00
☐	AT-115	E0971	Anti Tippers-pair	\$100.00
☐	LS16	K0038	16" width leg strap-calf support	Std
☐	LS18	K0038	18" width leg strap-calf support	Std
☐	LS20	K0038	20" width leg strap-calf support	Std
☐	BKR-EXT-S-105-16		16" width back rest extension	N/A
☐	BKR-EXT-S-105-18		18" width back rest extension	N/A
☐	PSG25F60-BLK	E1015	Frog Legs front fork suspension pair. Patient weight <160 lb	\$299.00
☐	PSG25F61-BLK	E1015	Frog Legs front fork suspension pair. Patient weight >160 lb	\$299.00
☐	CP2	E1399	Universal carry pouch-small	\$11.00
☐	CP4	E1399	Universal carry pouch-large	\$19.00
☐	CANO-115		UV protectant Canopy for 7/8" handle frame	N/A
☐	PB-115-16		Foldable rigidifying push bar for 16" width	\$88.00
☐	PB-115-18		Foldable rigidifying push bar for 16" width	\$88.00
☐	PB-115-20		Foldable rigidifying push bar for 16" width	\$88.00