

Full-Electric Bariatric Bed, 42"



FIGURE A

- Quiet, smooth operation
- Heavy-duty frame ensures strength and patient safety
- All steel construction
- Larger sleep surface than conventional bed
- Split-pan design with removable bed ends is easy to set up
- No tools necessary for set up
- Hand control (included) provides for multiple bed positioning for patients
- Low air loss/alternating pressure mattress available
- "T" rails available (Sold separately) (Figure A)
- Emergency manual crank included
- Weight Capacity 600 lbs

PRODUCTS	Description
15300	Bariatric Bed, 1/cs
15300BV-1HR	Bariatric Bed with 1 pair "T" Rails (item #15021BV), 1/cs
15300BV-PKG	Bariatric Bed Package (contains Bariatric Bed, 42" Bariatric Foam Mattress, item #15301, and 1 pair "T" Rails, item #15021BV, 1/cs

SPECIFICATIONS	Description
Bed Height	18.25" - 26.25"
Unit Dimensions	42" (W) x 88" (L)
Sleeping Surface	80"
Weight	231 lbs.
Carton Shipping Weight	276 lbs.
Weight Limit	600 lbs Maximum weight limit includes weight of occupant, all accessories such as mattress, side/assist rails, trapeze, any personal affects, etc.
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components

ACCESSORIES	Description
15301	1/cs

PARTS	Description
15307H	
15305	
15306	
15308	
15300P1005-1	1/ea
15300P1006-1	1/ea
15300P1051	1/ea
15300P1052	1/ea

PARTS	Description
15005PS1051	1/ea
15300MH	1/ea
15005P1025	1/ea
15300PH1014	1/ea
15005P1001	
15005P1028	1/ea
15005PF1028	1/ea
15005P1023	1/ea
15300HL	1/ea
15005HLS	1/ea
15300P1015L	1/ea
15300PF1015L	1/ea
15300EC	1/ea
15300CL	1/ea
15300MF	1/ea
15005P-1028	1/ea
15005JB	
15005JBN	
15005JBU	1/ea
15005HC	1/ea
15005HCN	1/ea
15005PC	1/ea
15005PF1025	1/ea
15300P1014	1/ea
15005P1016	1/ea
15300C	1/ea
15005S	1/ea
15005P1029	1/ea