

Delta™ Ultra-Light 1000 Full-Electric Low Bed



FIGURE A



FIGURE B



FIGURE C



FIGURE D



INSET E



FIGURE F

- The Delta 1000 is truly universal because the bed frame can be used with Drive's and most other manufacturers' bed ends, old and new
- The transition box mounted on the foot section is key to interchangeability (Figure A)
- The transition box allows you to change the rotation of the high/low shaft
- Headboard is taller than footboard to better fit home decor (Figure B)
- The wood grain panels are break and scratch resistant
- Once installed, the high/low shaft will stay with the bed to eliminate loss or forgetting on future installations (Figure C)
- Labels and color-coded springs assure accurate installation of side rails for patient safety (Figure D)
- Unique motor is completely self-contained to reduce weight and noise
- 9 ½" to 23 ½" deck height
- Easy to set up
- Top mounted Hi/Lo crank accessibility
- Height-adjustable motor comes pre-installed on the foot section
- Caster easily moves from standard to low bed position
- The motor assembly can be installed and removed with patient in bed
- Bed, side rails and mattress safety-tested by SGS, a nationally recognized testing facility Satisfies all applicable FDA entrapment guidelines
- The motor is UL approved
- In the event of a power failure, a 9 volt battery, located in the hand control, can lower the head and foot sections nine times
- New and improved hand pendant has large, easy-to-use controls (Figure E)
- Bed ships in two cartons
- Hand crank located on the top of the footboard prevents bending for the caregiver when in the low position (Figure F)
- Channel frame construction provides superior strength and reduced weight
- Zinc-coated spring deck
- Extension kit available (model #15030EXTKIT) extends length from 80" to 84"
- Patent Pending
- Weight Capacity 450 lbs

PRODUCTS

Description

PRODUCTS	Description
15235	Ultra-Light Full-Electric Low Bed
15235BV-FR	Ultra-Light Full-Electric Low Bed with Full-Length Side Rails (item #15001ABV), 1/cs
15235BV-HR	Ultra-Light Full-Electric Low Bed with Half Rails (item #15201BV), 1/cs
15235BV-PKG	Ultra-Light Full-Electric Low Bed with Full-Length Side Rails (item # 15001ABV) and 80" Innerspring Mattress (item #15006), 1/cs
15235BV-PKG-T	Ultra-Light Full-Electric Low Bed with Full-Length Side Rails (item #15001ABV) and 80" Therapeutic Support Mattress (item #15019), 1/cs
15235BV-PKG-1	Ultra-Light Full-Electric Low Bed with Half Rails (item #15201BV) and 80" Innerspring Mattress (item #15006), 1/cs

PRODUCTS

	Description
15235BV-PKG-1-T	Ultra-Light Full-Electric Low Bed with Half Rails (item #15201BV) and 80" Therapeutic Support Mattress (item #15019), 1/cs
15235BV-PKG-2	Ultra-Light Full-Electric Low Bed with Full-Length Side Rails (item #15001ABV) and 80" Foam Mattress (item #15007), 1/cs

SPECIFICATIONS

	Description
Bed Height Minimum Without Casters	9.5"
Bed Height Maximum With Casters	23.5"
Unit Dimensions	36" (W) x 88" (L)
Weight	175 lbs.
Carton Shipping Weight	191 lbs.
Weight Limit	450 lbs Maximum weight limit includes weight of occupant, all accessories such as mattress, side/assist rails, trapeze, any personal affects, etc.
Warranty	(Original Purchaser) 5 Years on Motor and Hand Control, Limited Lifetime on Welds and Frame, 1 Year on All Other Parts and Components

ACCESSORIES

	Description
15030HBC	Bedding in a Box - contains 1 pillow case, 1 fitted sheet, 1 cover sheet, and 1 blanket., 1/ea, 5/cs
15030HBL	Fitted Sheet - Contains 2 Fitted Sheets, 2 sheets per bx/ 10 bxs per cs

PARTS

	Description
15541	
15549	
15536	
15547H	
15547H-HP	
15547H-PF	
15544HGB	
15547F	
15547F-FP	
15547F-PF	
15544FGB	
15547H-F	
15547-CRANK	
15033MHLN-MGBA	
15033MHLN-GB	
15030UBS	
1503C	
15037MN	
15033HC-G	
15030PN-1035	
15030PN -1036	
15030-7RMM	
15030P-1006	
15030MC	
15030-MS	
15044CL	
15044CNL	
15044CSET	
15030P-1001	
15030-FRLINK	
15033MHLNC	
15033MELNC-N	
15033MB	
15033MBN-DELTA	
15544GBMB	