

## Prime Care Bed Model P903



- 600 lbs. weight capacity\*
- 7" super low height can help minimize falls (Figure A) while 30" high working height can facilitate staff & caregiver functions (Figure B)
- Integrated width with push button expansion from 36"-39"-42" provides extra measure of comfort and dignity (Figure I) 80", 84" and 88" sleep surface length helps comfortably accommodate the tallest residents
- 10 function Backlit Hand Pendant includes Auto Contour, Cardiac Comfort Chair (Figure D) and under bed night light provides an added measure of security at night
- Synchronized Hi/Lo control system keeps bed level for resident comfort and positioning
- Welded Steel construction, slat support deck is durable and easy to clean
- Roll at any height convenience with 3" heavy duty dual footprint casters plus steer control facilitate easy maneuverability (Figure H)
- Color coded, dual pedal floor lock system keeps bed securely in place when engaged (Figure G)
- Adjustable wall bumper helps protect walls and bed from damage (Figure E)
- Standard features include: adj. mattress retainers, heel lift, IV pole holders
- Optional: bed ends, Assist bars, rails, deluxe assist rails, trapeze, embedded staff control, battery backup, transporter
- Trendelenburg/Reverse Trendelenburg capable with staff control option (Figure C)
- Warranty: 15 yr. on frame, lifetime on welds, 3 yr. on electronics

PRODUCTS	Description
P903	Prime Care Bed Model P903, 1/ea

SPECIFICATIONS	Description
Sleep Surface	80", 84", 88"
Base Width	36"
Optional Width Expansion	36"- 39"-42"
Working Height	30"
Super Low Height	7"
Weight Capacity*	600 lbs
Bed Weight w/o options	222 lbs.
Warranty	Limited Lifetime on Welds; 15 Year on structural steel; 3 year on electrical & mechanical parts /components; 1 year on all other parts & components

\*Includes mattress and accessories

ACCESSORIES	Description
PLTCAB	Rotating Assist Bar, 1 pr
PLTCABL	Rotating Assist Bar, Left, 1/ea
PLTCABR	Rotating Assist Bar, Right, 1/ea
PLTCAR	Rotating Assist Rail, 1 pr
PLTCARL	Rotating 1/4 Assist Rail Left, 1/ea

**ACCESSORIES**
**Description**

PLTCARR	Rotating Assist Rail, Right, 1/ea
PLTCDR	Deluxe Rotating Rail, 1 pr
PLTCDRL	Deluxe Rotating Rail, Left, 1/ea
PLTCDRR	Deluxe Rotating Rail, Right, 1/ea
P903TRP	Transporter, 1/ea
P903TPZ	Trapeze, 1/ea
P903-BBK	Battery Backup w/Bracket, 1/ea

**PARTS**
**Description**

35BWBE-SO-SC03	35" Brownlee Woods Bed End w/SC03, Solar Oak, 1/ea
35BWBE-MW-SC03	35" Brownlee Woods Bed End w/SC03, Montana Walnut, 1/ea
35BWBE-FM-SC03	35" Brownlee Woods Bed End w/SC03, Figured Mahogany, 1/ea
35BWBE-WM-SC03	35" Brownlee Woods Bed End w/SC03, Windsor Mahogany, 1/ea
35MCBE-SO-SC03	35" Mill Creek Bed End w/SC03, Solar Oak, 1/ea
35MCBE-WC-SC03	35" Mill Creek Bed End w/SC03, Wild Cherry, 1/ea
35MCBE-MW-SC03	35" Mill Creek Bed End w/SC03, Montana Walnut, 1/ea
35MCBE-FM-SC03	35" Mill Creek Bed End w/SC03, Figured Mahogany, 1/ea
35MCBE-WM-SC03	35" Mill Creek Bed End w/SC03, Windsor Mahogany, 1/ea
35SHBE-WC-SC03	35" Smokey Hollow Bed End w/SC03, Wild Cherry, 1/ea
35SHBE-MW-SC03	35" Smokey Hollow Bed End w/SC03, Montana Walnut, 1/ea
35SHBE-FM-SC03	35" Smokey Hollow Bed End w/SC03, Figured Mahogany, 1/ea
35BWBE-B	35" BRNWOOD BED END,OSLO OAK, 1/ea
35BWBE-A	35" BRNWOOD BED END,OSLO CHERRY, 1/ea
35MCBE-B	35" MILL CRK BED END,OSLO OAK, 1/ea
35MCBE-A	35" MILL CRK BED END,OSLO CHERRY, 1/ea
42BWBE-B	42" BRNWOOD BED END,OSLO OAK,EA, 1/ea
42BWBE-A	42" BRNWOOD BED END,OSLO CHERRY, 1/ea
42MCBE-FM	42" MILL CRK BED END,OSLO FIG MAHOG, 1/ea
2MCBE-A-SC03	42" Mill Creek Bed End w/SC03, Cherry, 1/ea
42BWBE-SO-SC03	42" Brownlee Woods Bed End w/SC03, Solar Oak, 1/ea
42BWBE-WC-SC03	42" Brownlee Woods Bed End w/SC03, Wild Cherry, 1/ea
42BWBE-MW-SC03	42" Brownlee Woods Bed End w/SC03, Montana Walnut, 1/ea
42BWBE-FM-SC03	42" Brownlee Woods Bed End w/SC03, Figured Mahogany, 1/ea
42MCBE-SO-SC03	42" Mill Creek Bed End w/SC03, Solar Oak, 1/ea
42MCBE-MW-SC03	42" Mill Creek Bed End w/SC03, Montana Walnut, 1/ea
42MCBE-FM-SC03	42" Mill Creek Bed End w/SC03, Figured Mahogany, 1/ea
42MCBE-WM-SC03	42" Mill Creek Bed End w/SC03, Windsor Mahogany, 1/ea
42SHBE-SO-SC03	42" Smokey Hollow Bed End w/SC03, Solar Oak, 1/ea
42SHBE-WC-SC03	42" Smokey Hollow Bed End w/SC03, Wild Cherry, 1/ea
42SHBE-MW-SC03	42" Smokey Hollow Bed End w/SC03, Montana Walnut, 1/ea
42SHBE-FM-SC03	42" Smokey Hollow Bed End w/SC03, Figured Mahogany, 1/ea
42SHBE-WM-SC03	42" Smokey Hollow Bed End w/SC03, Windsor Mahogany, 1/ea
SP01-HB70-C101-00	
SP01-P903-16AE	
SP01-P903-15AE	
SP01-P903-10019E	
SP01-P903-10003E	
SP01-P903-10004E	
SP01-P903-10018E	
SP01-P010040209E-G	
SP01-P010040209E-R	
SP01-P903-10005E	
SP01-P903-14AE	
SP01-P903-13AE	