



Lightweight Expedition Aluminum Transport Chair With 12" Rear "Flat-Free" Wheels



FIGURE A



FIGURE B

- Weighs only 19 lbs
- Available in blue or red frames
- Folds flat for easy transport (Figure A)
- Handy lever on chair back makes it easy and convenient to fold back down (Figure B)
- Comes with seat belt
- Companion activated wheel locks
- 12" rear "flat-free" wheels

PRODUCTS	Description
EXP19LTBL	Blue, Expedition X-Light Blue TransChair, 1/cs
EXP19LTRD	Red, Expedition X-Light Red TransChair, 1/cs

SPECIFICATIONS	Description
Overall Width	22" (open) 10.75" (closed)
Overall height	38"
Seat-To-Floor Height	19"
Back Height	18"
Seat Dimensions	19" (W) x 16" (D)
Overall Length w/Riggings	38"
Weight	19 lbs.
Carton Shipping Weight	24 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

PARTS	Description
STDS871	Handgrip, 1/ea
STDS5D5501	Back Upholstery, 1/ea
STDS5D5502	Push Handle Blue, 1/ea
STDS5D5503	Push Handle Red, 1/ea
STDS5D5504	Back Upholstery Screws (8), 1/pkg
STDS5D0914	Folding Cane Hinge, 1/ea
STDS1070	Seat Rail Guide, 1/ea
STDS5D5505	Wheel Brake Left, 1/ea
STDS5D5506	Wheel Brake Right, 1/ea
STDS4020	tipping sleeve, 1/ea
STDS5D5507	Rear Axle Bearings, 1/pr
STDS5D5508	Rear Hub Cap, 1/ea
STDS5D5509	Rear Axle Bolt /sleeve, 1/ea
STDS5D5510	Rear Wheel, 1/ea

PARTS	Description
STDS5D5511	Cross frame Blue, 1/ea
STDS5D5513	Front Axle Bearings, 1/pr
STDS5D5514	Front Axle, 1/ea
STDS5D5515	Front Caster, 1/ea
STDS5D5516	Front Fork Blue, 1/ea
STDS5D5517	Front Fork Red, 1/ea
STDS5D5518	Fork Stem Bearing, 1/pr
STDS5D1069	Hanger Pin, 1/bag
STDS3J0929	Seat Upholstery 19", 1/ea
STDS5D5519	Seat Upholstery Screws, 1/pkg
STDS5D5520	Hand Brake/Cable, 1/ea
STDS3J091R	Footrest Release Lever Right, 1/ea
STDS3J091L	Footrest Release Lever Left, 1/ea
STDSFLVB	Arm Pad Full, 1/ea
STDS002001	arm pad bolts, 1/ea
EXPLTSF-RD	Footrest Pair Red, 1/pr
EXPLTSF-BL	Footrest Pair Blue
STDS5D1070	Seat Rail Guide, 1/ea