

Item #	Seat	Battery	UOM
SCOUTDST3	17.5" Folding Seat	20AH	1/ea




A



B




C



D



E



F

Scout DST 3-Wheel

Travel Scooter

- Dynamic Suspension Technology (DST) features front and rear suspension coils for the most comfortable ride on all terrain. (Figure D)
- Charge port located on the control panel and battery pack for easy charging.
- Headlight
- Interchangeable color panels in red and blue.
- Quick & easy disassembly.
- Armrests are padded, contoured, width and angle adjustable. (Figure F)
- Delta tiller with ergonomic throttle control. (Figure E)
- Armrest covers with reflectors.
- Large plastic basket with carry handle.
- Height adjustable swivel seat with fold down backrest.
- Anti-tip wheels.
- Extended footboard provides ample space for taller users. (Figure C)

specifications

Weight Capacity	325 lbs.	Freewheel Mode	Yes
Top Speed	4.25 mph	Front Wheels	8" x 2"
Estimated Range	15 miles	Rear Wheels	8" x 2"
Max Climbing Angle	6°	Anti-Tip Wheels	2" x 3/4"
Turning Radius	45"	Base Weight	60 lbs.
Ground Clearance	2.5"	Battery Weight (Pair)	30 lbs. (20 AH)
Dimensions	42.5" (L) x 21.25" (W) x 36" (H)	Seat Weight	23 lbs
Seat Dimensions	17.5" (W) x 16" (D)	Heaviest Piece Weight	Front Section: 34 lbs.
Back Height	14"	Total Weight	113 lbs.
Controller	CURTIS 45A	Warranty on Frame	Limited Lifetime
Motor	270W x 4700RPM	Warranty on Electronics, Controller, Drive-Train Components	14 months
Batteries	20AH	Warranty on Batteries	6 months
Charger	2A		
Brakes	Electromagnetic		