



Item #	Seat	UOM
2850HD-22	22" Captain Seat	1/ea
2850HD-24	24" Captain Seat	1/ea

accessories		
TRID-2A37	Swing-away Controller Arm - Trident	1/ea
SF8030	Cane/Crutch Holder	1/ea
SF8010	Oxygen Tank Holder	1/ea
TRID-2A37-LINX	Swing-Away Controller Arm	1/ea
LK3JELR	ELR Set	1/pr
SFELR	Front Rigging Hanger Bracket	1/ea



HCPCS  
**K0825**



## Trident HD

### Heavy-Duty Powerchair

- Shroudless design offers easy access to internal components and batteries
- Semi-Reclining captain seat with adjustable headrest is height and depth adjustable
- Armrests are padded, adjustable and removable
- Available with 22" or 24" width Captain's seat
- Quick-release lever allow for easy seat removal
- Flip-up footplate includes rubber bumper to prevent damage and is height, depth and angle adjustable
- Positioning belt included
- Programmable Controller
- Black, flat-free non-marking wheels
- Large, front-anti-tip wheels feature suspension for safety, smoother ride and enhanced obstacle handling
- Joystick is protected by metal casing
- Easy free-wheel operation
- Slide-out battery tray for easy access
- Optional swing-away controller arm sold separately
- Optional elevating leg rests sold separately

#### specifications

HCPCS	K0825	Charger	5A
Weight Capacity	450 lbs.	Brakes	Electromagnetic
Top Speed	4 mph	Caster Wheels	10" x 3"
Maximum Range	15 miles	Drive Wheels	14" x 3"
Turning Radius	29"	Base Weight	98 lbs.
Climbing Angle	6°	Battery Weight (pair)	66 lbs.
Ground Clearance	2.5"	Seat Weight	48 lbs.
Dimensions	48" (L) x 25.5" (W)	Total Weight	212 lbs.
Seat Dimensions	22" x 20"; 24" x 22"	Warranty on Frame	Limited Lifetime
Controller	PG VR2 70A	Warranty on Electronic Controller	24 Months
Motor	24V x 450W x 2700 rpm	Warranty on Drive Train Components	24 Months
Batteries (pair)	12V x 50AH	Warranty on Batteries	12 Months