

Go-Lite Bariatric Steel Rollator Padded Seat, Loop Locks



FIGURE A



FIGURE B

- 8" casters with soft-grip tires are ideal for indoor and outdoor use
- Soft padded oversized seat
- Comes with large 20.5" (W) x 4.5" (H) x 9" (D) basket that can be mounted under seat
- Tool-free removable padded backrest for comfort
- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure A)
- Handles easily adjustable in height with self-threading knob (Figure B)

PRODUCTS

Description

10215BL-1	Blue, E0149 / E0156, 1/cs
10215RD-1	Red, E0149+E0156, 1/cs

SPECIFICATIONS

Description

Handle Height	35.25" - 39.5"
Unit Dimensions	30.5" (W) x 26.75" (L)
Seat Dimensions	18" (W) x 22" (H) x 13" (D)
Weight	26 lbs.
Carton Shipping Weight	30.5 lbs.
Weight Limit	500 lbs.
Warranty	Limited Lifetime

PARTS

Description

10208GRP	1/ea
10208BRAKE	1/ea
9500021507	1/ea
95010215JHB	1/pr
10208K	1/ea
10215JK	1/pr
950300089901	1/ea
950000199901	1/ea
10208C	1/ea
9501021501B	1/ea
950000217-ER	1/ea
950000217-DL	1/ea
10215RA	1/ea
STDS002050	1/ea
10215RW	1/ea
9500004101	1/ea
950320139901	1/ea

PARTS	Description
10215F	1/ea
10215FA	1/ea
STDS1011	1/ea
9502P21503	1/ea
9501J21503	1/ea
10215S	1/ea
95108964225B	1/ea
10215W	1/ea
10215BSK	1/bx
950340059901	1/pr
951021820000	1/ea
10215SB	1/ea
9501A75000	1/ea