

## Nitro DLX Rollator



RETAIL BOX



INSET B

- Attractive, Euro-style design
- Handle height easily adjusts with unique push button
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- 10" front casters allow excellent steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Cross-brace design allows for side-to-side folding and added stability
- Durable and comfortable firm seat separates allowing unit to fold
- Handsome and removable deluxe zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight frame

### PRODUCTS

### Description

RTL10266CH-HS	Nitro, Firm Foam Seat, Champagne Frame, Champagne, 1/ea
RTL10266BL-HS	Nitro, Firm Foam Seat, Blue Frame, Blue, 1/ea

### SPECIFICATIONS

### Description

Handle Height	33.5" - 38.5"
Unit Dimensions	28" (L) x 24" (W)
Seat Dimensions	18" (W) x 10.75" (D)
Weight	22 lbs.
Seat Height	20"
Weight Capacity	300 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake Handle & Housing)	5 Years

### ACCESSORIES

### Description

892	Cane Holder, 1/ea
10266-CH	Nitro Cup Holder, 1/ea