

MedLite Beds are:

High-Quality

- Premium Linak® motors provide smooth reliable operation for years
- Recessed headspring supports and foot assembly reduce the chance of pinching and injuries

Light

- We have some of the industry's lightest beds, making transport, setup and storage much easier
- Our MedLite semi-electric bed weighs only 137 lbs (62 kgs)
- The 24V DC motors are also lightweight – and quiet

Strong

- An extra steel plate is welded to the middle of the frame ensuring strong and stable support
- Painted frames have three coats of enamel plus a clear coat for a strong, chip and scratch-resistant finish
- Durable plastic woodgrain end panels resist breakage
- Heavy-duty remotes resist breakage and moisture

Easy-To-Use

- The beds can be cleaned without having to remove the motor, and the electronics can be washed down with standard water pressure
- End panels are easy to wipe down and keep clean
- Smooth rolling, casters make the bed easy to move
- Convertible to a low bed using Med= Low Bed Conversion Kit – sold separately (see page 8)

MedLite Full-Electric Bed



MedLite Semi-Electric Bed

Among the lightest home care beds in the industry!



- ▶ Premium, quiet Linak motors provide dynamic movement without disturbing other family members
- ▶ Heavy-duty remote control
- ▶ Recessed headspring and foot assembly supports reduce the chance of pinching and injuries
- ▶ One of the longest warranties in the industry

Specifications

Assembled weight	145 lbs. (66 kg)
Height range	14.5" - 22.5" (37-57 cm)
Weight capacity	450 lbs. (204 kg)
Casters	Four (2 locking) 3" (8 cm)
Warranty	5 years for electrical and mechanical components; lifetime for the bed frame

Item No.	Description
MDR107003L	Full Electric Bed

"My dealership prefers Medline beds because they're a good value for the money. They're very durable and the long warranty gives us peace of mind. Thanks Medline!"

*Frank Lopez, AV Pulmonary Care,
Lancaster, CA*

- ▶ Features both automation and a traditional hand crank
- ▶ Premium, quiet Linak motors provide dynamic movement without disturbing other family members
- ▶ Weighs only 137 lbs (62 kg), one of the lightest beds in the industry
- ▶ Easy to clean without the need to remove the motor

Specifications

Assembled weight	137 lbs. (62 kg)
Height range	14.5" - 22.5" (37 -57 cm)
Weight capacity	450 lbs. (204 kg)
Casters	Four (2 locking) 3" (8 cm)
Warranty	5 years for electrical and mechanical components; lifetime for the bed frame

Item No.	Description
MDR107002L	Semi-Electric Bed

Standard on all Medline beds, the heavy-duty remote control resists breakage and moisture

