




S I G N A T U R E S E R I E S



**Pocketed Coil
Spring Seat
Standard!***



-   Hazelnut
-   Sterling
-   Buckskin*
-   Coffee Bean*

* upcharge applies

Imperial
PR404
shown in Coffee Bean



Imperial PR404

The Imperial with the industry's only articulating headrest and power lumbar system provides the ultimate in customized comfort. Recline and adjust your head, neck and shoulders for the best positioning using the powered adjustable headrest. The powered adjustable lumbar system offers support for the lower back and the pocketed coil seat system provides sacral comfort. A new designer chaise pad creates a mainstream look.



Imperial



Imperial	
Model Number	PR404
Size	Medium
Weight Capacity	375 Lb.
Recline Positions	3
Standard Fabrics	4
Fabric Options	ALL
Chaise Seat	YES
Coil Springs	YES
Back Type	Seam
Head Rest Cover	NO
Arm Covers	NO
Battery Back Up	YES
Overall Width	32.5"
Overall Height	43"
Width Between Arms	19"
Floor to Top of Seat	19"
Seat Depth	22"
Seat to Top of Back	28.5"
Seat to Top of Arm	6.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	34"
Back Angle Closed	110°
Back Angle Reclined	145°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	17"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Standard	YES
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	NO
Extra Pocket Option Available	NO
Heat and Massage Option Available	NO
Gel Memory Foam Option Available	NO
Power Pillow Option Available	Standard
AutoFlex Option Available	NO
Power Lumbar Option Available	Standard
Foot Rest Extension Option Available	NO
Straight Lift Option Available	YES
User Height Suggested	5'4" - 5'10"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



Headrest adjusts for personalized comfort



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
UL 73 and
CSA C22.2 No. 68

