

MAXICOMFORT® SERIES



MaxiComfort with
Twilight Positioning
PR514-MLA

MAXI
COMFORT®
with TWILIGHT



Hazelnut




Sterling




Brisa Coffee Bean*

* upcharge
applies


Coming
Summer 2018




Copper




Calypso



Shiraz



Brisa Buckskin*



Brisa Bridle*

* upcharge applies

PR514-MLA

Discover weightlessness in the MaxiComfort Cloud with Twilight positioning. Our patented three-motor technology offers an array of rejuvenating options to cradle the body and ease the mind. Achieve extreme Zero Gravity and TV positions while enjoying the world's finest rejuvenating seating experience!

 GOLDEN

PROUDLY MADE

BY AMERICANS
with Domestic & Imported Materials

PR514-MLA

MAXICOMFORT® SERIES

Model Number	PR-514-MLA
Size	Medium/Large
Weight Capacity	375 LB
Recline Positions	MaxiComfort
Standard Fabrics	8
Fabric Options	All
Chaise Seat	Yes
Coil Springs	No
Back Type	Biscuit
Head Rest Cover	Yes
Arm Covers	No
Battery Back Up	Yes
Overall Width	39.5"
Overall Height	46.5"
Width Between Arms	20"
Floor to Top of Seat	20.5"
Seat Depth	21"
Seat to Top of Back	31"
Seat to Top of Arm	7.5"
Seat Angle Raised	15°
Floor to Top of Foot Rest Reclined	25"
Floor to Top of Back Reclined	25"
Back Angle Closed	105°
Back Angle Reclined	170°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	24"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Removable Back Option Available	No
Smart Tek System Standard	Yes
Customizable *	No
Left Hand Control Option Available	Yes
Extra Pocket Option Available	Standard
Heat & Massage Option Available	No
Gel Memory Foam Option Available	No
Power Pillow Option Available	No
Auto Flex Option Available	No
Foot Rest Extension Option Available	No
Straight Lift Option Available	No
User Height (suggested)	5'7" -6'2"

* See Customize Chart for Dimension limits
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

Our luxurious lift and recline chairs are equipped with electronics to help assure a secure lift every time you need it.



TV Watching in Twilight



Lounge Reclined in Twilight



Zero Gravity in Twilight

Chair arm removed for demonstration purposes

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with
 UL 73 and
 CSA C22.2 No. 68



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com

SHT-PR514-040518