



Relaxer  
PR756L



Relaxer  
PR756MC



PR756L  
FABRICS



Hazelnut



Sterling



Copper



Calypso



Shiraz

PR756MC  
FABRICS



Hazelnut



Sterling



Copper



Calypso



Shiraz



Buckskin\*



Coffee Bean\*

\* upcharge applies

Relaxer  
PR756MC  
PR756L

The Relaxer is the original powered lift chair recliner featuring MaxiComfort® Ultimate Recline Technology standard! The Relaxer model is a true Golden customer favorite! Featuring the most luxurious four-pillow waterfall back design, and overstuffed, super plush rounded armrests, the Relaxer is available in two sizes – Medium and Large – for the absolute perfect fit. Every MaxiComfort chair comes standard with our exclusive AutoDrive™ 3.0 hand control to make reclining the chair into your favorite position a breeze!



# Relaxer



**Relaxer**  
PR-756MC      **Relaxer**  
PR-756L

<b>Model Number</b>	<b>PR-756MC</b>	<b>PR-756L</b>
<b>Size</b>	Medium	Large
<b>Weight Capacity</b>	375 LB	375 LB
<b>Recline Positions</b>	MaxiComfort	MaxiComfort
<b>Standard Fabrics</b>	7	5
<b>Fabric Options</b>	ALL	ALL
<b>Chaise Seat</b>	YES	YES
<b>Coil Springs</b>	NO	NO
<b>Back Type</b>	4 Pillow	4 Pillow
<b>Head Rest Cover</b>	YES	YES
<b>Arm Covers</b>	NO	NO
<b>Battery Back Up</b>	YES	YES
<b>Overall Width</b>	36.5"	39"
<b>Overall Height</b>	44"	46"
<b>Width Between Arms</b>	20"	22"
<b>Floor to Top of Seat</b>	19"	20.5"
<b>Seat Depth</b>	21"	22"
<b>Seat to Top of Back</b>	29"	30"
<b>Seat to Top of Arm</b>	8"	8"
<b>Seat Angle Raised</b>	15°	15°
<b>Floor to Top of Foot Rest Reclined</b>	23"	24"
<b>Floor to Top of Back Reclined</b>	21"	22"
<b>Back Angle Closed</b>	100°	100°
<b>Back Angle Reclined</b>	174°	174°
<b>Foot Rest Angle Closed</b>	90°	90°
<b>Foot Rest Angle Opened</b>	180°	180°
<b>Distance Required from Wall Reclined</b>	25"	27"
<b>Extension Reclined from Back of Seat to End of Foot Rest</b>	43"	44"
<b>Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism</b>	YES	YES
<b>Warranty on Electrical Parts</b>	3 Yrs.	3 Yrs.
<b>Pro Rated Years 4 Thru 7</b>	YES	YES
<b>Warranty on Mechanical Labor</b>	3 Yrs.	3 Yrs.
<b>Warranty on Electrical Labor</b>	1 Yr.	1 Yr.
<b>Removable Back Option Available</b>	NO	NO
<b>Smart Tek System Standard</b>	YES	YES
<b>Customizable</b>	NO	NO
<b>Left Hand Control Available</b>	YES	YES
<b>Extra Pocket Option Available</b>	YES	YES
<b>Heat and Massage Option Available</b>	YES	YES
<b>Foot Rest Extension Option Available</b>	YES	YES
<b>Straight Lift Option Available</b>	NO	NO
<b>User Height (suggested)</b>	5'4"-5'10"	5'11"-6'2"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness.  
All measurements taken on concrete floor using levels and metal rulers.



501 M shown

## Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

## Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

## Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



## Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



Literature is current at time of printing.  
Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68  
Listing # E112417  
FNAE 80-214A

