

**New  
Power Cloud!**



PR512MLA  
Medium/Large



PR510SME  
Small/Medium



PR-510MLA  
Medium/Large



**PR510-SME-  
FABRICS**



Hazelnut



Sterling



Copper



Calypso



Shiraz

**PR510-MLA &  
PR512-MLA  
FABRICS**



Hazelnut



Sterling



Copper



Calypso



Shiraz



Buckskin\*



Coffee Bean\*

\* upcharge applies

**MaxiComfort  
Cloud  
PR510 & PR512**

The Cloud is the world's most comfortable lift chair recliner which comes in two sizes, Small/Medium and Medium/Large. Featuring our custom bucket seat, this chair is the ultimate in relaxation and rejuvenation.

The newest edition to the Cloud family, the PowerCloud, also features our exclusive PowerPillow articulating headrest and a modern back design.



# MaxiComfort Cloud


M A X I C O M F O R T™ S E R I E S

	Cloud	Cloud	Power Cloud
<b>Model Number</b>	PR510-SME	PR510-MLA	PR512-MLA
Size	Small/Med	Med/Large	Med/Large
Weight Capacity	375 LBS	375 LBS	375 LBS
Recline Positions	MaxiComfort	MaxiComfort	MaxiComfort
Standard Fabrics	5	7	7
Fabric Options	ALL	ALL	ALL
Chaise Seat	YES	YES	YES
Coil Springs	NO	NO	NO
Back Type	Biscuit	Biscuit	Seam
Head Rest Cover	YES	YES	NO
Arm Covers	NO	NO	NO
Battery Back Up	YES	YES	YES
Overall Width	33.5"	38"	38"
Overall Height	41"	45"	45"
Width Between Arms	18.5"	20"	20"
Floor to Top of Seat	19"	20"	20"
Seat Depth	19.5"	21.5"	21"
Seat to Top of Back	27"	29"	30"
Seat to Top of Arm	7.5"	7"	7"
Seat Angle Raised	15°	15°	15°
Floor to Top of Foot Rest Reclined	23"	24.5"	24.5"
Floor to Top of Back Reclined	23"	24"	24.5"
Back Angle Closed	105°	105°	105°
Back Angle Reclined	170°	170°	170°
Foot Rest Angle Closed	90°	90°	90°
Foot Rest Angle Opened	180°	180°	180°
Distance Required from Wall Reclined	23"	24"	24"
Extension reclined from Back of Seat to end of Foot Rest	43"	44"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.
Removable Back Option Available	NO	NO	NO
Smart Tek System Standard	YES	YES	YES
Customizable	NO	NO	NO
Left Hand Control Option Available	YES	YES	YES
Extra Pocket Option Available	Standard	Standard	Standard
Heat & Massage Option Available	YES	YES	YES
Gel Memory Foam Option Available	NO	NO	NO
Power Pillow Option Available	NO	NO	Standard
AutoFlex Option Available	NO	YES	NO
Foot Rest Extension Option Available	NO	NO	NO
Straight Lift Option Available	NO	NO	NO
User Height (suggested)	5'1" - 5'6"	5'7" - 6'2"	5'7" - 6'2"

\* See Customize Chart for Dimension Limits  
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Literature is current at time of printing.  
 Golden Technologies reserves the right to make changes to the product or literature at any time.



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com

SHT-PR510+512-092616



501 M shown

## Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

## Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

## Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



## Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



## Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.