

LIGHTWEIGHT PORTABLE SCOOTER



**Comfort Spring System!**



Seat Color  
Black Vinyl

Body Colors  
Blue  
Red



3 wheel GB117S  
4 wheel GB147S

The next generation of the popular Buzzaround scooter series is now available with our comfort spring system! The Buzzaround XLS offers more standard features to make your riding experience better than ever. Features a wraparound Delta tiller, a new LED battery gauge and a new LED angle-adjustable headlight. The new Buzzaround XLS is also longer to offer more foot room, with a great turning radius. The Buzzaround XLS disassembles in seconds with just four easy steps making transporting it from one destination to another a breeze!





# LIGHTWEIGHT PORTABLE SCOOTER

Item	GB117S	GB147S
Description	XLS	XLS
Medicare Code	N/A	N/A
Weight Capacity	300 lbs.	300 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to <sup>1</sup>	8 miles	8 miles
Under Deck Clearance	3"	3"
Turning Radius	37"	47"
Type Batteries	2-12AH	2-12AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp Dynamic R-Series	50 Amp Dynamic R-Series

Dimensions		
Length	42.5"	41.5"
Width	21.5"	21.5"
Height (ground to top of back of seat)	33"-35"	33"-35"
Ground to Top of Seat	21"-23"	21"-23"
Ground to Top of Deck	4.75"	4.75"
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	33.5"	32"

Weight of Unit: (assembled) <sup>3</sup>		
Front Half	104 lbs.	112 lbs.
Rear Half	35 lbs.	43 lbs.
Seat w/Arms	28 lbs.	28 lbs.
Batteries	22 lbs.	22 lbs.
	19 lbs. pack	19 lbs. pack

Tires:		
Front	Solid	Solid
Rear	8"	8"
Rear Anti-Tip Wheels	8"	8"
	2"	2"

Standard Seat		
Seat Back Height - no headrest	Stadium Style	Stadium Style
Width x Depth	16"	16"
Color	18" x 16"	18" x 16"
	Black Vinyl	Black Vinyl

Optional Seat:		
Back Height - no headrest	Stadium Style	Stadium Style
Width x Depth	16"	16"
Color	20" x 18"	20" x 18"
	Black Vinyl	Black Vinyl

Warranty:		
Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics <sup>2</sup>	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

<sup>1</sup> Battery range tested at 200 lbs., and will vary due to rider weight, drive surface and terrain.  
<sup>2</sup> Electronics warranty excludes batteries.



**Comfort Spring System**

## Wireless disassembly in 4 easy steps!



**Step 1:** Pop off the seat!



**Step 2:** Pull out the battery pack!



**Step 4:** Pull forward on the drivetrain release lever to separate the front and rear frame sections!

**Step 3:** Fold down tiller and lock in place!



Disassembles into 5 pieces!

4 wheel model shown

## Features

- 300 lb. weight capacity
- Angle-adjustable LED headlight
- Tiller mounted basket
- Comfort Spring System

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

