

Stockton Youth Bed

Made in USA



Standard Features

- Head & Knee Elevation- Manual or Electric
- Hi-Lo System- Manual or Electric
- Foot End Control Box - Electric only
- EZ-Lift Sides with Fingertip Release
- 5" Innerspring Mattress
- Pan Mattress Surface
- 5" Heavy Duty Casters (1 Steer, 2 Brake, 1 Swivel)
- Red/Yellow/Blue Laminate Panels

Options

- IV Poles- Add up to 4
- Heavy Duty Bumper Pads
- 8" Pressure Reduction Mattress
- Patient Hand Control (Electric models only)
- Spring Link Mattress Surface
- Additional Colors Available for End Panels
- Additional Epoxy Colors Available

****Stockton Beds with Trendelenberg are available upon special request****

Articulation

- Head to 35 degrees
- Thigh/Foot to 15 degrees

Our Stockton Youth Beds are intended for use by an individual who has limited mobility (cannot sit up, high kneel or stand) and will not be at risk of falling out of bed. The side rails extend 15" above the surface of the 5" innerspring mattress, and 12" above the surface of the optional pressure reduction mattress.

HEMOCARE STOCKTON BEDS: MANUAL AND ELECTRIC SPECIFICATIONS

MODEL # SPECIFICATIONS	Electric N1030P	Manual N1035P	Electric N2070P	Manual N2085P
Overall Dimensions:	36" x 72"	36" x 72"	36" x 83"	36" x 83"
Mattress Size:	34" x 69"	34" x 69"	34" x 80"	34" x 80"
Floor to Top of Mattress Platform				
Highest Position	29"	29"	29"	29"
Lowest Position	19"	19"	19"	19"
Height of Sides	16.5"	16.5"	16.5"	16.5"
Overall Height				
Highest Position	50"	50"	50"	50"
Lowest Position	40"	40"	40"	40"



HARD

The nations #1 manufacturer of hospital cribs

LAMINATE COLOR OPTIONS



Red, Yellow & Blue
(Standard Colors)



White



Antique White



Red



Yellow



Green (Manual Beds Only)



Blue



Violet (Manual Beds Only)



Light Oak Wood Grain

Actual colors may vary slightly.

800-873-4273 ~ homecare@hardmfg.com