

MVP-502

Ergonomic V-Seat Recliner



- Patented Ergonomic/Anti-Sliding “V-Seating” prevents sliding, sores and relieves overall pressure on your body!
- Backrest reclines smoothly from 93-160 degrees.
- Rigidifying/Folding push-bar provides and strengthens back support
- Ergonomic Armrest allows the user to adjust the armrest position when reclined
- “Traveling-Axle” Leg-rest automatically stretches smoothly to reduce pressure/compression on abdomen.
- 3D Elevating Leg-rest support Ergonomically Designed to allow greater micro-adjustment: up/down, left/right or inward/outward
- Extending handbrakes equipped!
- Strap-Tension adjustable backrest & seat for ultimate comfort adjustments and stabilization for each user!
- Ergonomic Upholstery; Also pressure relieving 1.25” to 3” thick breathable material, offers superb comfort throughout the body
- 7” x 1” flat-free casters
- 22” x 1” flat-free rear wheels

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
MVP-502-16	16	18	9-11	20	19	52	24	14	47	36	250
MVP-502-18	18	18	9-11	20	19	52	26	15	47	36	250
MVP-502-20	20	18	9-11	20	19	52	28	16	47	36	250



KARMA HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335