

T-900 Transport Class

**TRANSPORTER
DELUXE WHEELCHAIR**

| SPECIFICATIONS | T-920W | T-922W |
|---------------------|----------|----------|
| Seat Width | 20 in. | 22 in. |
| Seat Depth | 18 in. | 18 in. |
| Seat Height | 21 in. | 21 in. |
| Back Height | 16 in. | 16 in. |
| Overall Height | 36 in. | 37 in. |
| Overall Open Width | 25 in. | 27 in. |
| Weight w/o Riggings | 40 lbs. | 44 lbs. |
| Weight Capacity | 450 lbs. | 450 lbs. |
| HCPCS Code | E1038 | E1038 |

(SUBJECT TO CHANGE WITHOUT NOTICE)



FIND THE **PERFECT** FIT FOR YOU...

- Padded detachable desk length armrests
- Dual Cross Bars
- 19" Seat with carry loops for transporting convenience
- 12½ rear solid polyurethane tires
- Swing-away footrests w/ heel loop
- Manual lock brakes
- Frame Color: Red
- High-grade flame retardant, breathable nylon upholstery