

Cruiser III Wheelchair

Lightweight, Dual-Axle Wheelchair



HCPCS
K0003

All the features of the original
Cruiser III and the New Built-in
Seat Rail System and
Extendable Upholstery

- Built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18" (Figure A - C)
- Weighs under 36 lbs (excluding front riggings)
- Carbon steel frame with silver vein finish
- Removable flip-back arms allow for easy transfer (Figure D)
- New frame style eliminates seat guides and allows for custom back inserts and accessories
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- Composite, Mag-style wheels are lightweight and maintenance free
- 8" front casters are adjustable in three positions
- Padded armrests provide additional comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings (Figure E)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Dual axle provides easy transition of seat height to hemi-level
- Comes with push-to-lock wheel locks


STEP 1



STEP 2



STEP 3



Item #	Seat Width	Product Overall Width Open	Armrest	Product Front Riggings	UOM
K316DDA-SF	16"	24"	Flip Back Detachable Desk Arm	Swing-away Footrest	1/cs
K316DDA-ELR	16"	24"	Flip Back Detachable Desk Arm	Elevating Leg rest	1/cs
K316DFA-SF	16"	24"	Flip Back Detachable Full Arm	Swing-away Footrest	1/cs
K316DFA-ELR	16"	24"	Flip Back Detachable Full Arm	Elevating Leg rest	1/cs
K316ADDA-ELR	16"	24"	Flip Back Detachable and Adjustable Height Desk Arm	Elevating Leg rest	1/cs
K316ADFA-ELR	16"	24"	Flip Back Detachable and Adjustable Height Full Arm	Elevating Leg rest	1/cs
K316ADFA-SF	16"	24"	Flip Back Detachable and Adjustable Height Full Arm	Swing-away Footrest	1/cs
K316ADDA-SF	16"	24"	Flip Back Detachable and Adjustable Height Desk Arm	Swing-away Footrest	1/cs
K318DDA-SF	18"	26"	Flip Back Detachable Desk Arm	Swing-away Footrest	1/cs
K318DDA-ELR	18"	26"	Flip Back Detachable Desk Arm	Elevating Leg rest	1/cs
K318DFA-ELR	18"	26"	Flip Back Detachable Full Arm	Elevating Leg rest	1/cs
K318DFA-SF	18"	26"	Flip Back Detachable Full Arm	Swing-away Footrest	1/cs
K318ADDA-SF	18"	26"	Flip Back Detachable and Adjustable Height Desk Arm	Swing-away Footrest	1/cs
K318ADFA-SF	18"	26"	Flip Back Detachable and Adjustable Height Full Arm	Swing-away Footrest	1/cs
K318ADFA-ELR	18"	26"	Flip Back Detachable and Adjustable Height Full Arm 26"	Elevating Leg rest	1/cs
K318ADDA-ELR	18"	26"	Flip Back Detachable and Adjustable Height Desk Arm	Elevating Leg rest	1/cs
K320DDA-SF	20"	28"	Flip Back Detachable Desk Arm	Swing-away Footrest	1/cs
K320DDA-ELR	20"	28"	Flip Back Detachable Desk Arm	Elevating Leg rest	1/cs
K320DFA-SF	20"	28"	Flip Back Detachable Full Arm	Swing-away Footrest	1/cs
K320DFA-ELR	20"	28"	Flip Back Detachable Full Arm	Elevating Leg rest	1/cs
K320ADDA-SF	20"	28"	Flip Back Detachable and Adjustable Height Desk Arm	Swing-away Footrest	1/cs
K320ADDA-ELR	20"	28"	Flip Back Detachable and Adjustable Height Desk Arm	Elevating Leg rest	1/cs
K320ADFA-ELR	20"	28"	Flip Back Detachable and Adjustable Height Full Arm	Elevating Leg rest	1/cs
K320ADFA-SF	20"	28"	Flip Back Detachable and Adjustable Height Full Arm	Swing-away Footrest	1/cs

specifications

Width Closed	12"	Product Weight (with 18" Seat)	36 lbs.
Height	35"	Product Weight (with 20" Seat)	38 lbs.
Seat--To-Floor Height	17.5" - 19.5"	Weight Limit (with 16" and 18" Seat)	300 lbs.
Back Height	16"	Weight Limit (with 20" seat)	350 lbs.
Seat Depth	16" - 18"	Carton Shipping Weight	43 lbs. or less
Length with Riggings	42"	Warranty	Limited Lifetime
Product Weight (with 16" Seat)	35 lbs.		

accessories

STDS801GT	Wheel Lock Extension	1/ea	STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS801	Wheel Lock Extension	1/ea	STDS851	Seat Belts, Hook and Loop Fastener Type Closure	1/ea
STDS1002	Standard 8" Caster	1/ea	STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS806	Anti-Fold Bar	1/ea	STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric	1/ea
STDS819	Anti-Tippers With Wheels	1/pr	STDS820	Telescoping IV Pole Attachment	1/ea
STDS835	Chart Carry Pocket	1/ea	STDS860	Universal Forward Stabilizer	1/pr
STDS831	Heel Loops	1/bx	STDS803	Universal Oxygen "E" Cylinder Holder	1/ea
WASR	Limb Support, Right	1/ea	STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea
STDS823	Overhead Anti-Theft Device	1/ea	STDS1034	Wheelchair Cane/Crutch Holder	1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks	1/ea	STDSSHK-K3	Super Hemi-Kit	1/ea
STDS834	Overhead Anti-Theft Device Single Pole	1/ea			
STDS3J24SF	Plastic Foot plate	1/pr			
LK3JELR	Plastic Foot plate	1/ea			