


Gravity 8 Deluxe Long Term Care Pressure Redistribution Mattress



INSET B

- This mattress includes a superior combination of multi-layered, latex free, die-cut and memory foam for therapeutic pressure redistribution
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution
- Mattress includes a 6" to 5" heel slope and a channel cut center layer of foam under the foot section to target heel ulcer prevention
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
- Stretch cover is fluid and tear resistant, low shear and easy to remove and clean
- High Density Foam Perimeter provides solid edge support and facilitates safe patient transfers
- Available with or without raised side rails
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS

PRODUCTS	Description
15886	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15870	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15884	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15888	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea
15885	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea

SPECIFICATIONS

SPECIFICATIONS	Description
Cover Material	Stretch, Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	3
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633
Weight Capacity	400 lbs.
Warranty	10 Year Limited Foam 1 Year Cover

ACCESSORIES

ACCESSORIES	Description
15870C	Replacement Cover for 15870, 1/ea
15884C	Replacement Cover for 15884, 1/ea
15885C	Replacement Cover for 15885, 1/ea
15886C	Replacement Cover for 15886, 1/ea
15888C	Replacement Cover for 15888, 1/ea
158709C	Replacement Cover