


Gravity 9 Premium Long Term Care Pressure Redistribution Mattress



INSET A

- The ultimate in the industry for multi-layered die-cut foam and true memory foam, this mattress provides therapeutic pressure redistribution and is latex free
- Mattress includes a tapered heel slope and channel cut center layer of foam under the entire lower leg to offload heels and prevent heel pressure ulcers
- Top layer is full-length and full width open-cell visco-elastic memory foam
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly
- Stretch cover is fluid and tear resistant, low shear and easy to clean
- High Density Foam Perimeter provides solid edge support and facilitates safe patient transfers
- Meets flammability standards CFR 1632 and CFR 1633
- Available with or without raised side rails

PRODUCTS

Description

15996	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15970	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15984	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15970-4280	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15970-4284	Multi-Layered/Multi-zoned Foam Mattress, E0184, 1/ea
15977	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea
15985	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea
15970-4280-RR	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea
15970-4284-RR	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out, E0184, 1/ea

SPECIFICATIONS

Description

Cover Material	Stretch, Fluid-Resistant Nylon Top Waterproof Vinyl Bottom
Mattress Layers of Foam	4
Side Rails	Yes
Dedicated Heel Section	Yes
Weight Capacity	450 lbs.
Warranty	36" W- 15 Years Limited Foam 1 Year Cover, 42" W- 7 Years Limited Foam 1 Year Cover

ACCESSORIES

Description

15970C	Replacement Cover for 15970, 1/ea
15977C	Replacement Cover for 15977, 1/ea
15984C	Replacement Cover for 15984, 1/ea
15985C	Replacement Cover for 15985, 1/ea
15996C	Replacement Cover for 15996, 1/ea