

Med-Aire Alternating Pressure Pump and Pad System



- Quiet pump alternately inflates and deflates 130 individual bubble air cells
- 4 LPM (liters per minute) pump produces consistent air flow and pressure over a 5 minute cycle time and is available fixed or variable (Fig.A)
- Variable pressure setting on pump allows for individualized comfort increasing patient compliance
- Durable, heavy gauge, premium vinyl bubble pad provides superior therapeutic treatment and is available with or without end flaps
- Pump has built-in brackets for mounting easily to bed

PRODUCTS	Description
14001E	Variable Pressure Pump and Pad, E0181, 1/ea
14001EFD	Variable Pressure Pump and Pad with End Flaps, E0181, 1/ea
14001E-220	Variable Pressure Pump (220V) and Pad, E0181, 1/ea
14002E	Fixed Pressure Pump and Pad, E0181, 1/ea
14002E-220	Fixed Pressure Pump (220V) and Pad, E0181, 1/ea
14001ECAN	14001E w/Hospital Grade Plug, 1/ea, E0181, 1/ea

SPECIFICATIONS	Description
Pump Airflow	4 LPM
Pump Cycle Time	5 Minutes
Pump Type	Compressor
Weight Capacity	300 lbs.
Warranty	1 Year Control Unit, 90 Day Pad

ACCESSORIES	Description
14003	Standard Replacement Pad (replaces Mason #2001); 78" (L) x 36" (W) x 2.5" (H), 1/ea
14003-EF	Replacement Pad with End Flaps (replaces Mason # 2001-FL); 123" (L) x 36" (W) x 2.5" (H), 1/ea
14001T	Beige Tubing for 14001 and 14001E, 1/ea
14005E	Variable Pressure Pump Only for 14001E (Replaces Mason # 2002-ADJ), 1/ea
14006E	Fixed Pressure Pump Only for 14002E (Replaces Mason # 2002-II), 1/ea