

Power Mobility Products

P320 Junior Super

Light Mini Powerbase Wheelchair



P320 Junior

Super Light Mini Powerbase Wheelchair

□ Technical Specs :

max.speed : 4.5mph
 range : up to 12.5mi.
 turning radius : 23"
 ground clearance : 2"
 motor: DC 24V, 120W
 controller : Dynamic Shark 40A / PG VR2 40A
 battery : 12V / 22A x2
 charger : 3A off-board
 gradient : 8°
 caster : 7 x 1 3/4" solid tires
 drive wheel : 9" foam filled tires
 brake : intelligent, regenerative
 electromagnetic brakes


□ Options :

Elevating leg rests. Oxygen tank holder.
 Cane & crutch holder. Walker holder.
 Foldable rear basket.

□ Weight : (lb)

base	battery/ea	seat	tl w/batt.	capacity
66	14	40	135	300

□ Dimensions : (in)



Medicare Codes:
 K0823 P320-2S C
 K0822 P320-2S S



□ Features :

- . Lightweight and easy to disassemble for handling, storage and transporting.
- . Ideal as personal secondary portable power wheelchair.
- . Programmable controller with a wide range of selections to suit the needs and demands of different users.
- . All standard power wheelchair features.
- . Excellent maneuverability.
- . Off-board charger.
- . Colors available: Red, Blue.