

Vision CF

P322A

Front-Wheel Drive

Powerbase Wheelchair



 **merits**

Vision CF

P322A

Front-Wheel Drive

Powerbase Wheelchair

- ◆ Weight capacity 300 lbs
- ◆ Detachable, height and width adjustable armrests
- ◆ Removable, replaceable seat and back covers
- ◆ 18" seat width, ergonomic seating design
- ◆ Semi-reclining back rest 90°-135°
- ◆ Simplified frame construction
- ◆ Compact, lightweight design
- ◆ Easy access to components
- ◆ 9" Drive wheels, 6" castors
- ◆ Excellent maneuverability
- ◆ Flip up footplate



Technical Specs

Max Speed	up to 4 mph
Range	up to 18 miles
Turning Radius	21"
Ground Clearance	2.4"
Motor	DC 24V, 130W
Controller	Dynamic LiNX
Battery	12V/U1(34AH) x 2pcs
Charger	4A off-board
Gradient	6°
Casters	6" foam filled tires
Drive Wheels	9" foam filled tires
Brakes	Intelligent, regenerative, electro-magnetic

Dimensions

Seat Width	18"
Seat Depth	16"
Seat Height	20" - 22"
Backrest Height	17"
Min. Footrest Extension	16.5"
Max. Footrest Extension	18.5"
Overall Length	33"
Overall Height	42"
Overall Width	24"
Total Weight	154 lbs
Capacity	300 lbs

www.meritsusa.com



730 N.E. 19th Place, Cape Coral, FL 33909

TEL: 1-900-963-7487, 1-239-772-0579

Sales fax: 1-239-574-2661