

Vision Sport

P326A

Mid-Wheel Drive

Powerbase

Wheelchair



 **merits**

1-800-963-7487 • Fax: 239-772-3252

www.meritsusa.com

P326A Vision Sport



Adjustable Armrests

Width and height Adjustment / Flip up Armrest



Easily disassembles



True six wheel design / Front & rear suspension
Helps climb over everyday obstacles smoothly



Excellent Stability



Lynx Joystick

Wireless Bluetooth programmability



Optional:

swing away joystick bracket

Dimensions	
Length	34"
Width	24"
Height	48.5"
Ground Clearance	2"
Seat Dimensions	
Width	20"
Depth	18"
Back Height	20"
Weight Capacity	300-lbs
Maximum Speed	5 mph
Range up to	18 miles
Tires	
Front/Rear Casters	6"
Drive Wheel	10"
Turning Radius	20"
Incline Capability	6°
Product Weight	
Total Weight without Battery	132-lbs
Total Weight with Battery	184-lbs
Weight of Heaviest Piece	90-lbs
Battery	
Type	U1 gel x 2
Weight	26-lbs x 2
Warranty	
Frame	5 year / Limited Lifetime
Electronics	1 year
Motor/Brakes	1 year
Batteries	1 year

- ❖ Mid-wheel drive design delivers outstanding maneuverability in tight, compact spaces and a smooth stable ride both indoors and outdoors even over rough terrain
- ❖ Full front suspension design offers a comfortable ride
- ❖ Low center of gravity design delivers a smooth and stable ride
- ❖ Height and width adjustable armrest, semi reclining backrest with height adjustable seat
- ❖ Full Suspension
- ❖ No rear tipping on footplate standing (up to 250-lbs)
- ❖ Dual in-line motors for enhanced efficiency, torque, range and performance
- ❖ Semi-Reclining Back
- ❖ Solid Caster Wheels & Foam Filled Drive Wheels
- ❖ 4A off-board Charger
- ❖ Lap Belt (Standard with chair)

