



COLOR & HYDROGRAPHIC OPTIONS



Radiant Red



Power Pink



Ocean Blue



Phantom Grey



Mossy Oak® Breakup® Country™



Sunburst Orange



Galactic Green



Titan Blue



Midnight Black



Carbon Fiber

FABRIC OPTIONS



Leatherette



Mesh

SPECIFICATIONS

Maximum User Weight	300 lbs	Seat-to-Floor Height	17.5", 18.5", 19.5"
Maximum Speed	7.5 mph	Seat Elevator	14"
Range	16 mi*	Tilt Options	0° – 50° Posterior
Base Width	25.5"	Recline Options	85° – 180° (power)
Base Length	43"	Legrest Angle	90° – 180°
Minimum Turning Radius	30"	Backrest Height	23" – 28" (by 1" incr.)
Ground Clearance	3"	Seat Depth	14" – 22" (by 1" incr.)
Recommended Battery Type	Group 24	Seat Width	17", 19", 21", 23"
Weight Including Batteries	421 lbs	Armrest Pad Lengths	10", 13", 16", 18"
Drive Electronics	R-Net 120A		

**Specifications may vary based upon configuration.*

**Actual driving range will vary based on driving and battery conditions.*

permobil

POWER WHEELCHAIRS

STAND WITHOUT
COMPROMISE



F5
Corpus VS

F5 Corpus VS

Adjustable Backrest Angle

Up to 17.5"
Seat-To-Floor
Height

We've combined all the features you love about the F5 Corpus with adjustable standing capabilities to create the F5 Corpus VS. Tailor your standing sequence with personalized features like adjustable backrest angles and articulating power footplates.

Experience the vast array of social and health benefits that come from standing. All of this with the freedom to drive at any point of the standing sequence.

NEW PROGRAMMABLE STAND SEQUENCE
ULTIMATE ADJUSTABILITY

Knee block and articulating power footplate for extra support.



LOWER & MORE COMFORTABLE RIDE



UNCOMPROMISED TILT & RECLINE

INTELLIGENT ANTI-TIPPERS
AUTOMATICALLY DEPLOY TO MAINTAIN
STABILITY DURING STANDING

