



USB  
Charger  
Port

LC-580i shown in  
Ultraleather® Charcoal fabric

MADE IN  
**USA**  
OF FOREIGN & DOMESTIC COMPONENTS

# Oasis

## COLLECTION

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
RECLINERS

The **Oasis** brings a new level of comfort to Pride® Power Lift Recliners. Overstuffed padding provides elegant comfort and design, while contoured memory foam seating and pocketed coils deliver individualized support and relaxation. Features include a unique quad pocket design for extra storage and an easy-to-use hand control with USB charger port. The Oasis is available in true infinite-position and 3-position.

US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |    

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —

# Oasis

COLLECTION

AMERICA'S #1  
POWER LIFT RECLINER

POWER  
LIFT  
RECLINERS



LC-580i shown in  
Saratoga Charcoal fabric

## OASIS FABRIC SELECTION

QS

### SARATOGA



CASHMERE



CHARCOAL



EGGPLANT



GODIVA



NAVY

### LEXIS STA-KLEEN®

sta-kleen.  
PERFORMANCE FABRIC



BLACK



MUSHROOM



BURGUNDY



CHESTNUT

### ULTRALEATHER™

Uf  
UltraLeather



BUFF



FUDGE



GARNET



PECAN



CHARCOAL

### ARIA<sup>2</sup> CRYPTON®

super fabric



ESPRESSO



COOL GREY



RED



LAZULI



SAND

Look for the Quick Ship  
Indicator for the shortest  
lead times available.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.  
\*Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



MADE IN USA



MADE IN USA



MADE IN USA

## SPECIFICATIONS

	Medium	Large	Medium/Large
	LC-580iM	LC-580iL	LC-380
	True Infinite-Position Zero Gravity, Trendelenburg Chaise Lounger 5'4" - 6'0"	True Infinite-Position, Zero Gravity, Trendelenburg Chaise Lounger 5'9" - 6'2"	3-Position Chaise Lounger 5'4" - 6'0"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	20"	21"	21"
Seat Depth	20"	21"	21"
Seat Width	20"	21"	24.5"
Top of Back to Seat	26"	29"	25"
Back Style	Blown Fiber Biscuit Back	Blown Fiber Biscuit Back	Blown Fiber Biscuit Back
Fully Padded Chaise	Yes	Yes	Yes
Heat & Massage	Optional	Optional	Optional
Seating Type	Pocketed Coil Springs w/ Memory Foam Topper	Pocketed Coil Springs w/ Memory Foam Topper	Pocketed Coil Springs w/ Memory Foam Topper
Distance From Wall	20"	21"	20"
Overall Chair Width	39"	39"	40"
Chair Weight	130 lbs.	132 lbs.	120 lbs.
Warranty	7-Year Pro-rated*	7-Year Pro-rated*	7-Year Pro-rated*

## BENEFITS

- Premium furniture grade laminate/hardwood frames for added strength and durability
- Contoured memory foam seating and pocketed coils offer individualized comfort and support
- Sinuous no-sag back springs combine durability with just the right amount of flexibility for personalized comfort
- Overstuffed biscuit back design provides body-conforming comfort
- Specially designed arms afford a luxurious feel and firm support
- Industry-exclusive quad pocket design offers extra storage – easy access top pockets provide convenient storage for your hand control and other devices
- Easy-to-use hand control with USB charger port  
(not available with Heat and Massage)
- Removable filling allows for custom lumbar support
- Headrest and armrest covers included

## OPTIONS\*

- Deluxe Heat & Massage (LC-580iM/L)
- Heat & Massage (LC-380 only)
- Residential cleaning kit
- Head pillow
- Lumbar pillow (Std. on LC-580iM/L)
- Protective pad
- Wireless remote (LC-380 only)

Additional lead time may  
apply for some options.

\*Additional charges apply for all options.

1. Batteries not included.  
2. Crypton® Fabric offers a five-year warranty  
on the moisture barrier. Does not cover  
punctures in the fabric.

US: 1-800-800-8586 • CAN: 1-888-570-1113  
pridemobility.com |

**Pride**  
Mobility Products Corp.  
— LIVE YOUR BEST —