

# Jazzy<sup>®</sup> 2.0

## SELECT 6

The Jazzy Select 6<sup>®</sup> 2.0 offers a truly unique ride experience with a sporty design and superior performance. Pride's patented Active-Trac<sup>®</sup> Suspension combines with Mid-Wheel 6<sup>®</sup> drive design for a smooth, comfortable and stable ride, both indoors and out. The depth-adjustable high-back seat folds flat for greater convenience, while offering maximum comfort and support.

**THIS FOOT PLATFORM OFFERS SUPERIOR ADJUSTABILITY FOR AMPLE FOOT ROOM.**



AMERICA'S FAVORITE  
POWER CHAIR

Learn more online.



[gojazzy.com](http://gojazzy.com)

**THIS CONVENIENT FOLDING SEAT IS DEPTH ADJUSTABLE FOR MAXIMUM COMFORT AND OFFERS REPLACEABLE BACK AND SEAT COVERS.**

# Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed<sup>1</sup></b>	Up to 4.3 mph
<b>Ground Clearance<sup>2</sup></b>	2.5" at footplate, 1.77" at motor bracket
<b>Turning Radius<sup>2</sup></b>	23.76"
<b>Overall Length<sup>2</sup></b>	43.27"
<b>Overall Width</b>	23.27"
<b>Drive Wheels</b>	10" flat free
<b>Front Wheels</b>	6" flat free
<b>Rear Wheels</b>	6" flat free
<b>Suspension</b>	Active-Trac® (ATX) with Mid-Wheel 6® Technology
<b>Drivetrain</b>	Two-motor, Mid-Wheel 6® drive
<b>Braking System</b>	"Intelligent Braking" electronic regenerative, disc park brake
<b>Standard Electronics</b>	40 amp, PG GC 3 controller
<b>Battery Charger</b>	Off-board, 3.5 amp
<b>Per-Charge Range<sup>2,4</sup></b>	Up to 19 miles
<b>Battery Requirements<sup>3,5</sup></b>	(2) 12 volt, deep cycle, U-1
<b>Battery Weight</b>	24.5 lbs. (U-1)
<b>Base Weight</b>	170 lbs. without batteries
<b>Standard Seat Weight</b>	50 lbs. (comfort, high-back)
<b>Weight of Heaviest Piece</b>	114 lbs.
<b>Max. Seat Dimensions</b>	20" x 20" (high-back)
<b>Warranty</b>	5-year limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Colors:




Sapphire  
Blue



Ruby  
Red



Pearl  
White



Black  
Onyx

\*Additional charges may apply.

# Features

- Active-Trac® Suspension combines with Mid-Wheel 6® Technology for a smooth, stable ride over uneven terrain
- Convenient, limited recline, high-back black vinyl folding seat is depth adjustable for maximum comfort and includes replaceable back and seat covers available in sizes: 16" x 16"-18", 18" x 18"-20" and 20" x 18"-20"
- Height-adjustable, flip-up foot platform offers ample foot room for maximum comfort
- Black, flat-free, non-marking tires saves your floors
- In-line joystick mount is equipped with performance-matched electronics for easy operation
- Equipped with 60" lap belt for security and peace of mind

## Optional Seating\*

- Limited recline, high-back black vinyl solid seat pan with depth adjustment available in sizes: 16" x 16"-18" or 18" x 18"-20"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 20" x 18" or 20" x 20"
- Power elevating seat (VR2 6-key required)

## Options

- Elevating leg rests
- Swing-away footrests
- Heel loops
- Swing-away joystick
- Residual limb support

## Accessories

- Cane/crutch holder
- Rear basket
- Cup holder
- Saddle bag
- Flag holder
- Walker holder
- Lap belt (50", 70", 80", 90")
- Weather cover
- Military patches
- XLR USB charger
- Oxygen tank holder

## Model

**Jazzy Select 6 2.0** | Code: K0822/K0823



Find a compatible lift for your vehicle online:  
[PRIDEMOBILITY.COM/LIFT-FINDER](https://PRIDEMOBILITY.COM/LIFT-FINDER)

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>5</sup> Battery weight may vary based on manufacturer.



(US) 800-800-8586 • (Canada) 888-570-1113 • [pridemobility.com](https://pridemobility.com) •



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.