

Jazzy[®]

ELITE HD

The Jazzy[®] Elite HD features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.



LARGE 14" DRIVE WHEELS



PERFORMANCE-MATCHED ELECTRONICS



COMPACT CASTERS FOR SUPERIOR INDOOR MANEUVERABILITY

AMERICA'S FAVORITE POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	450 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	3"
Turning Radius	24"
Overall Length	35.5"
Overall Width	25.25"
Drive Wheels	14" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, front-wheel drive
Braking System	Intelligent braking (electronic, regenerative disc brakes)
Standard Electronics	60A, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 9.8 miles
Battery Requirements	(2) 12 volt, deep cycle, NF-22
Battery Weight	37 lbs. each
Base Weight	99 lbs.
Standard Seat Weight	52 lbs. (high-back)
Standard Seat	20" x 18-20"
Max. Seat Dimensions	24" x 24"
Warranty	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Jazzy Red



Jazzy Blue

Features

- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 70" lap belt

Options

- Synergy® Seat
- Solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag
- Heel loops
- Lap belt (50", 60", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Residual limb support
- Military patches
- Swing-away joystick
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

Models

Jazzy Elite HD-2HD-SS	Code: K0824
Jazzy Elite HD-2HD-C	Code: K0825