

# Jazzy®

## 600 ES

The Jazzy® 600 ES features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an extraordinary level of performance, maneuverability and the tightest turning radius in the Jazzy line. The Jazzy 600 ES offers true all-around performance indoors and out making it an excellent power chair choice.



### LARGE 14" DRIVE WHEELS



SPHERICAL  
OMNI-CASTERS



ACTIVE-TRAC®  
ATX SUSPENSION

AMERICA'S FAVORITE  
POWER CHAIR

Learn more online.



[gojazzy.com](http://gojazzy.com)

## Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	2.75"
<b>Turning Radius</b>	20.5"
<b>Overall Length</b>	35.5"
<b>Overall Width</b>	24.5"
<b>Drive Wheels</b>	14"
<b>Front Wheels</b>	6" solid OMNI-casters
<b>Rear Wheels</b>	6" solid OMNI-casters
<b>Suspension</b>	Active-Trac® (ATX)
<b>Drivetrain</b>	Two-motor, in-line, mid-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	50A, Dynamic Shark controller
<b>Battery Charger</b>	Off-board, 5 amp
<b>Per-Charge Range</b>	Up to 17 miles
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, U-1
<b>Battery Weight</b>	25 lbs. each
<b>Base Weight</b>	124 lbs.
<b>Standard Seat Weight</b>	52 lbs. (high-back)
<b>Standard Seat</b>	18" x 18-20" and 20" x 18-20"
<b>Max. Seat Dimensions</b>	20" x 20"
<b>Warranty</b>	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## Standard Color



Cardinal Red

## Features

- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## Options

- Synergy® Seat
- Solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Heel loops
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Military patches
- Swing-away joystick
- XLR USB charger

## Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

[pridemobility.com/Lift-Finder](http://pridemobility.com/Lift-Finder)

## Models

<b>Jazzy 600 ES 2S-SS</b>	Code: K0822
<b>Jazzy 600 ES 2S-C</b>	Code: K0823