

# Jazzy<sup>®</sup>

## ELITE ES

The Jazzy<sup>®</sup> Elite ES delivers a reliable blend of power, performance and style for your active lifestyle. The in-line, front-wheel drive technology provides tight turns around corners for optimal indoor maneuverability in small spaces and large front wheels for extra absorption over obstacles.

### IN-LINE FRONT-WHEEL DRIVE TECHNOLOGY



**BLACK, FLAT-FREE  
NON-MARKING TIRES**



**COMFORT HIGH-BACK SEAT  
WITH HEADREST**



**AMERICA'S FAVORITE  
POWER CHAIR**

Learn more online.



[gojazzy.com](http://gojazzy.com)

## Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	0.75"
<b>Turning Radius</b>	23.5"
<b>Overall Length</b>	34.75"
<b>Overall Width</b>	22.75"
<b>Drive Wheels</b>	9" solid
<b>Front Wheels</b>	3" solid anti-tips
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, in-line, front-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	40 amp, PG GC 3 controller
<b>Battery Charger</b>	Off-board, 3 amp
<b>Per-Charge Range</b>	Up to 12 miles
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, U-1
<b>Battery Weight</b>	25 lbs. each
<b>Base Weight</b>	69 lbs.
<b>Standard Seat Weight</b>	43 lbs.
<b>Standard Seat</b>	18" x 18-20"
<b>Max. Seat Dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## Colors



Jazzy Red



Jazzy Blue

## Features

- Color-through shroud
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## Options

- Solid seat pan
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap Belt (50", 70", 80", 90")
- Elevating legrest
- Swingaway legrest
- Heel loops
- Cup holder
- Military patches
- Height-adjustable footrest bracket
- Swing-away joystick
- XLR USB charger

## Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

[pridemobility.com/Lift-Finder](http://pridemobility.com/Lift-Finder)

## Models

**Jazzy Elite ES-2S-SS** Code: K0822

**Jazzy Elite ES-2S-C** Code: K0823