

Jazzy[®]

SELECT 6

The Jazzy Select[®] 6 incorporates sporty design and performance features for a truly unique ride experience. Pride's patented Active-Trac[®] Suspension combined with six wheels on the ground and an innovative articulating beam allow for superior performance indoors and out, with added stability over uneven terrain and ramp transitions.

ACTIVE-TRAC[®] SUSPENSION



JAZZY ARMOR[®] AND
COLOR-THROUGH
SHROUD



COMFORT HIGH-BACK SEAT
WITH REPLACEABLE BACK
AND SEAT COVERS



AMERICA'S FAVORITE
POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.75"
Turning Radius	21.75"
Overall Length	35.5"
Overall Width	23.5"
Drive Wheels	10"
Front Wheels	6" solid casters
Rear Wheels	6" solid casters
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 15 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	25 lbs. each
Base Weight	100 lbs.
Standard Seat Weight	43 lbs. (comfort, high-back)
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Jazzy Red



Jazzy Blue

Features

- Jazzy® Armor® includes shroud and controller guards to protect against daily wear and tear
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Options

- Synergy® Seat
- Solid seat pan
- Power elevating seat (non-swivel)
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

Models

Jazzy Select 6-SS	Code: K0822
Jazzy Select 6-C	Code: K0823