

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com**Quantum 4Front™ Order Form with TRU-Balance 3 Positioning**

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

**INTRODUCTION**

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


Thank you for choosing Quantum!**PATIENT INFORMATION**

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

**1. BASE MODEL**

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 3MP-SS..... \$7,895

Part: 4FRONT 3MP-SS. HCPCS: K0861

Multiple actuator power option.

Quantum 4Front 3SP-SS \$7,895

Part: 4FRONT 3SP-SS. HCPCS: K0856

Single power option.

Quantum 4Front 3S-SS..... \$7,895

Part: 4FRONT 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Account Number: _____ Date: _____ Marked for: _____

2. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Patient's Weight: _____

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1700109-PRD1	 White Satin Part: DGN1700109-PWT2	 Grey Street Part: DGN1700109-PSL3	 Yellow Blaze Part: DGN1700109-PYL3
 Electric Blue Part: DGN1700109-PBU9	 Pretty in Pink Part: DGN1700109-PPK3	 Silver Bullet Part: DGN1700109-PSL7	 Military Green Part: DGN1700109-PGN8
 Back in Black Part: DGN1700109-PBK1	 Deep Purple Part: DGN1700109-PPR5	 Midnight Blue Part: DGN1700109-PBU10	To Be Determined Part: N/A
 Night Train Part: DGN1700109-PBK3	 Lime Light Part: DGN1700109-PGN7	 Orange Crush Part: DGN1700109-POR6	

3.1b. Hydrographic Shrouds

☐  Military Digital Camo \$250 Part: DGN1700109-WCM7	☐  Faded Glory \$250 Part: DGN1700109-WFL2	☐  Carbon Fiber \$250 Part: DGN1700109-WCB1
--	--	---

4. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color.

4.1. Drive Wheels

3"x14" Flat Free Black Tires with Accent Standard Part: PTOINDV3525	4"x14" Flat Free Knobby Black Tires \$275.50 ea. Part: WHL147601 x2
3"x14" Pneumatic Black Tires with Accent No Charge Part: WHL1704059 x2	4"x14" Pneumatic Knobby Black Tires \$275.50 ea. Part: WHL147599 x2

4.2. Drive Wheel Accent Color

Not available with knobby tires.

Back in Black Standard Part: DGN164245-PBK1 x2	Deep Purple No Charge Part: DGN164245-PPR5 x2	Yellow Blaze No Charge Part: DGN164245-PYL3 x2
Red Corvette No Charge Part: DGN164245-PRD1 x2	Lime Light No Charge Part: DGN164245-PGN7 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Electric Blue No Charge Part: DGN164245-PBU9 x2	Grey Street No Charge Part: DGN164245-PSL3 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
White Satin No Charge Part: DGN164245-PWH2 x2	Silver Bullet No Charge Part: DGN164245-PSL7 x2	Faded Glory \$50 ea. Part: DGN164245-WFL2 x2
Pretty in Pink No Charge Part: DGN164245-PPK3 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Carbon Fiber \$50 ea. Part: DGN164245-WCB1 x2
Night Train No Charge Part: DGN164245-PBK3 x2	Orange Crush No Charge Part: DGN164245-POR6 x2	To Be Determined Follow-up required.

5. FENDER LIGHTS

REQUIRED. Fender lights are standard and will be automatically included unless omitted or if base mounted lights are selected.

1. Fender Light Selection

Switch Operated Fender Lights Standard Part: DGN1700033	Omit Fender Lights No Charge Part: DGN1700032
---	---


2. Fender Light Switch Mount Type

Side Rail/Track Mount Bracket Standard Part: MEC1800866 Allows the egg switch to be mounted to seating side rails using a t-nut assembly. This may also be mounted to TRU-Balance 3 flip-back armrest and other seat or accessory bar rails. Not available with contoured seating.	Under Armrest Mounting Bracket \$50* Part: MEC1802726 *No charge with contoured seating. Allows the egg switch to mount to a joystick tube-style bracket mounted to any armrest. Includes joystick mounting bracket.
---	---

6. FULL LIGHT PACKAGES

Optional. Requires Q-Logic 3 Joystick for Full Light Package. Select one available package below.

Base Mounted, LED, Full Light Package \$995 Part: ELE1700035 Fender lights will be omitted with Base Mounted, LED, Full Light Package.	Seat Mounted, LED, Full Light Package \$995 Part: ELE131477
---	---



7. NO FRONT ANTI-TIP OPTION

Optional. With this option, maximum speed will be 5 mph and 2.7 mph when elevated, maximum weight will be 200 lbs., and seat towers will be fixed in lowest position.

No Front Anti-Tips No Charge Part: PTOINDV3851
--

8. JOYSTICK/SCIM

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

8.1a. Omit Joystick

No Joystick, Setup for Q-Logic 3 Electronics
Part: Configuration selection.

8.1b. Q-Logic 3 Non-Expandable Capacity Controls

Q-Logic 3 non-expandable capacity allows control up to one actuator.

Q-Logic 3 Non-Expandable Capacity Joystick Standard
Part: CTL167724-NE

Q-Logic 3 Non-Expandable Capacity Joystick No Charge
for Full Light Package
Part: CTL167725-NE

8.1c. Q-Logic 3 Expandable Capacity Controls

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

1. Q-Logic 3 Expandable Capacity Joysticks

Q-Logic 3 Expandable Capacity Joystick\$700
Part: CTL167724

Q-Logic 3 Expandable Capacity Joystick\$700
for Full Light Package
Part: CTL167725
Required with full light packages.

**2. Q-Logic 3 Expandable Capacity Required Components**

Q-Logic 3 expandable capacity controls require the controller module and expandable harness below. Customer Service will include these items in a quote when expandable controls or specialty electronics are selected on this form.

Q-Logic 3 Expandable Capacity Controller Module..... \$695
Part: See XRef. HCPCS: E2377

Q-Logic 3 Expandable Capacity Controller Harness..... \$425
Part: HAR144005. HCPCS: E2313

8.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN124307

Match Joystick Armor to Base Color \$50
Part: PTOINDV3421

Non-matching Joystick Armor\$50
Part: PTOINDV3422
Select color below

Part below will be selected with system configurator.

Red Corvette
Part: DGN124336-PRD1

White Satin
Part: DGN124336-PWH2

Grey Street
Part: DGN124336-PSL3

Orange Crush
Part: DGN124336-POR6

Military Digital Camo
Part: DGN124336-WCM7

Electric Blue
Part: DGN124336-PBU9

Pretty in Pink
Part: DGN124336-PPK3

Silver Bullet
Part: DGN124336-PSL7

Yellow Blaze
Part: DGN124336-PYL3

Faded Glory
Part: DGN124336-WFL2

Back in Black
Part: DGN124336-PBK1

Deep Purple
Part: DGN124336-PPR5

Midnight Blue
Part: DGN124336-PBU10

Military Green
Part: DGN124336-PGN8

Carbon Fiber
Part: DGN124336-WCB1

Night Train
Part: DGN124336-PBK3

Lime Light
Part: DGN124336-PGN7

8.3. Q-Logic SCIM Module

Optional.

SCIM Module Kit for Q-Logic 3..... \$525
Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 Expandable Capacity Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

**9. JOYSTICK HANDLES**

Optional. The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$85
Part: ACCASMB1860. HCPCS: E2323

Extended I-Handle (G) \$85
Part: ACCASMB1857. HCPCS: E2323

Chin Cup (C) \$85
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

Bodypoint 3" Goal Post Handle \$99
Part: PLSKNCP1057. HCPCS: E2323

T-Bar (D) \$120
Part: ACCASMB1858. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$99
Part: PLSKNCP1058. HCPCS: E2323

Medium Ball Handle (E) \$65
Part: ACCASMB1859. HCPCS: E2323


Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308

Soft Ball Handle (F) \$65
Part: ACCASMB1856. HCPCS: E2323



10. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.


Omit Joystick Mounting Bracket

10.1a. Fixed, In-line Joystick BracketFixed, In-line, Bracket Standard
Q-Logic Part: FRMASMB10617Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount \$245
Part: ACC142524. HCPCS: E1028Swing-away, 1" Offset, Left Side Mount \$245
Part: ACC142742. HCPCS: E1028Swing-away, Inline, Right Side Mount \$245
Part: ACC142525. HCPCS: E1028Swing-away, 1" Offset, Right Side Mount \$245
Part: ACC142743. HCPCS: E1028**10.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$350
Q-Logic Part: MEC138860. HCPCS: E1028
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Stealth Swing-Away/Flip Down Bracket \$294.25
QL3 SAJ: ST-TWBM480CJ-Q. HCPCS: E1028
QL3 Joystick Part: ST-TWBM480QL-Q. HCPCS: E1028
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.Multi-axis Mount \$100
Q-Logic Part: ACC131216. HCPCS: E1028
Allows height adjustment of joystick.**Specialty Mount Side**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**11. ELECTRONICS ACCESSORIES****Optional.**USB Mobile Device Charger \$175
Part: ELE124122
Can charge 1 device. See multiplier harness section for secondary multiplier.Q-Logic PC Programming Disk with Cable \$750
Part: ELEASMB5216**12. TRU-BALANCE 3® SEATING - QUANTUM INTELLIGENT SYSTEMS™****12.1. TRU-Balance 3 Configuration****REQUIRED.** Select one seating system. Be sure to complete the seat dimensions, seat to floor, and back selections.TRU-Balance 3 Static Seat \$2,295
Part: PTOINDV3442
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.TRU-Balance 3 Power Recline \$7,595
Part: PTOINDV3444. HCPCS: E1004
175° of power recline.TRU-Balance 3 Power Tilt \$6,995
Part: PTOINDV3443. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.TRU-Balance 3 Power Tilt & Recline \$11,895
Part: PTOINDV3446. HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.**12.2. 10" Power Lift or Manual Recline****Optional.** Select one option below.TRU-Balance 3 Power 10" Lift \$2,600
Part: PTOINDV3775. HCPCS: E2300
10" of scissor lift. 22" W x 22" D maximum seat size.Attendant Manual Recline \$500
Part: PTOINDV3501. HCPCS: E1226
Available with Static Seat or Power Tilt. Not available with Power Lift. Will require extended lead time when selected with Articulating Foot Platform or Fender Lights.**12.3. Occupied Transit****Optional.** Available in seat dimensions of 16" to 20" wide and 16"-20" deep. Requires TRU-Balance 3 armrests and 60" Sure-lok lap belt. Only available with swing-away legrests and center mount foot platforms.Occupied Transit \$275
Part: KIT152047To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>**12.4. Seat Dimensions - Width****REQUIRED with seat selection.** 15" width is only available with non-matching back width. Minimum seat width is 12". Maximum seat width is 22".

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

TRU-Balance 3 seating section continued on following page.

12. TRU-BALANCE 3 SEATING CONTINUED**12.5. Seat Dimensions - Width Package****REQUIRED with seat selection.** Width and depth package part: See XRef: TB3SPACK**12"-16" Width Package**

Not available with Power Lift.

14"-16" Width Package

Only available with Power Lift.

16"-20" Width Package**20"-24" Width Package**

Maximum width is 22".

12.6. Seat Dimensions - Depth**REQUIRED with seat selection.** Minimum seat depth is 12". Maximum seat depth is 22".**12" D**
(#1002)**13" D**
(#1003)**14" D**
(#1004)**15" D**
(#1005)**16" D**
(#1006)**17" D**
(#1007)**18" D**
(#1008)**19" D**
(#1009)**20" D**
(#1010)**21" D**
(#1011)**22" D**
(#1012)**12.7. Seat Dimensions - Depth Package****REQUIRED with seat selection.** Width and depth package part: See XRef: TB3SPACK**12"-16" Depth Package**

Not available with Power Lift.

14"-16" Depth Package

Only available with Power Lift.

16"-20" Depth Package**20"-24" Depth Package**

Maximum depth is 22".

12.8. Seat to Floor Height**REQUIRED.** Select a seat to floor height compatible with the seat functions previously selected. Seat to floor height stated below is based on unoccupied height. Actual height may vary depending on weight of patient.**Seat to Floor - Static Seat or Power Recline Only****15" Seat to Floor**

Part: PTOINDV3847

16" Seat to Floor

Part: PTOINDV3848

17" Seat to Floor

Part: PTOINDV3840

Seat to Floor - Power Tilt Only, Power Tilt and Recline, or any Power Lift**17.75" Seat to Floor**

Part: PTOINDV3649

18.75" Seat to Floor

Part: PTOINDV3650

19.75" Seat to Floor

Part: PTOINDV3651

12.9. Back Type Selection**REQUIRED.** Match choice in this section with the options in the following back cushion configuration section.**TRU-Balance 3 Back Assembly (Static Back)**

Part: Configuration selection

Select a back cushion configuration in the subsequent section.

Back Cane Option

Part: Configuration selection

Available only with Static Seat and Power Tilt selections. Not available with ADI Q-Backs. Select cane type, back cane height, and back cushion configuration below.

Flip Forward Back..... \$150

Part without Lift: MEC131957, Part with Lift: MEC148117

Available with Power Tilt and Lift. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-15" width or with attendant manual recline. Select a back cushion type in the subsequent section.

12.10. Back Cane Selection**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type****Angle Adjustable..... No Charge**

Part: MEC140266

1" Endo..... \$35

Part: MEC141455

Flip Forward..... \$25

Part: MEC140272

2" Endo..... \$35

Part: MEC141456

2. Back Cane Height**16"H Cane**

Part: FRMASMB15834

20"H Cane

Part: FRMASMB15838

24"H Cane

Part: FRMASMB15842

17"H Cane

Part: FRMASMB15835

21"H Cane

Part: FRMASMB15839

25"H Cane

Part: FRMASMB15843

12.11. Back Cushion Configuration - Type**REQUIRED.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**Sport Back..... No Charge**

Part: See XRef

Complete following width package and size.

TRU-Comfort 2 Back, 20" Width or less..... \$675

Part: See XRef. HCPCS: E2620

Complete following width package and size.

ADI Q-Back, Active Contour..... \$525

Part: See XRef. HCPCS: E2613

Complete following width package and size.

TRU-Comfort 2 Back, 22" Width or greater..... \$775

Part: See XRef. HCPCS: E2621

Complete following width package and size.

ADI Q-Back, Deep Contour..... \$595

Part: See Ref. HCPCS: E2615

Complete following width package and size.

Omit Back Pan/Cushion..... No Charge

Part: PTOINDV1287

Select a following width package.

ADI Q-Back, Extra Deep Contour..... \$650

Part: See XRef. HCPCS: E2620

Complete following width package and size.

After-market Back Interface..... \$200

Part 16"-20"W: SET131920. Part 20"-24"W: SET131921

20" Wide Extra Deep Contour is not compatible with Gel and Waterfall Armpads, Endo/Meso Downposts, or 2 Post Arms unless wider arm bars are installed on the chair.

Enter the set width below, range of 16"-24". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs. A Special Order is required on back widths less than 16".

TRU-Balance 3 seating section continued on following page.

12. TRU-BALANCE 3 SEATING CONTINUED**12.12. Back Cushion Configuration - Back Width****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections

Part: PTOINDV3468

10"W-14"W Package

Part: Configuration selection, back cushions

16"W-20"W Package

Part: Configuration selection, back cushions

20"W-24"W Package

Part: Configuration selection, back cushions

12.13. Back Cushion Configuration - Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.12" W
(#3562)14" W
(#3563)16" W
(#3564)17" W
(#3565)18" W
(#3566)19" W
(#3567)20" W
(#3568)21" W
(#3569)22" W
(#3570)24" W
(#3571)**12.14. Back Cushion Configuration - Overall Back Height****Static Seat and Power Tilt available height range: 14"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.**

14" Height

Part: SETHGHT1027

22" Height

Part: SETHGHT1035

15" Height

Part: SETHGHT1028

23" Height

Part: SETHGHT1036

16" Height

Part: SETHGHT1029

24" Height

Part: SETHGHT1037

17" Height

Part: SETHGHT1030

25" Height

Part: SETHGHT1038

18" Height

Part: SETHGHT1031

26" Height

Part: SETHGHT1039

19" Height

Part: SETHGHT1032

27" Height

Part: SETHGHT1040

20" Height

Part: SETHGHT1033

28" Height

Part: SETHGHT1041

21" Height

Part: SETHGHT1034

	Tilt or Static Width												Recline Width											
	12"	14"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	16"	17"	18"	19"	20"	21"	22"	24"				
Overall Height	14"	ST	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
	15"	ST	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
	16"	S	ST	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-				
	17"	S	ST	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-				
	18"	-	-	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-				
	19"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	S	ST	ST	ST	ST	ST	S	ST	S				
	20"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	S	ST	ST	ST	ST	ST	S	ST	S				
	21"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	-	-	ST	ST	ST	ST	ST	S	ST	ST			
	22"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	23"	-	-	ST	ST	ST	ST	ST	S	ST	ST	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	24"	-	-	S	S	ST	S	ST	S	S	S	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	25"	-	-	S	S	ST	S	ST	S	S	S	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	26"	-	-	-	-	-	-	-	-	-	-	-	-	ST	ST	ST	ST	ST	S	ST	ST			
	27"	-	-	-	-	-	-	-	-	-	-	-	-	S	S	ST	S	ST	S	S	S			
	28"	-	-	-	-	-	-	-	-	-	-	-	-	S	S	ST	S	ST	S	S	S			

S

Sport Back

ST

Sport or TRU-Comfort 2 Back

SQT

Sport, Q-Back, or TRU-Comfort 2 Back

Sport, Q-Back, and TRU-Comfort 2 Back Availability Chart

S Sport Back ST Sport or TRU-Comfort 2 Back

SQT Sport, Q-Back, or TRU-Comfort 2 Back

Sport, Q-Back, and TRU-Comfort 2 Back Availability Chart

12.15. TRU-Balance 3 Back Shroud**REQUIRED.** The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.Red Corvette
Part: DGN148034White Satin
Part: DGN148075Grey Street
Part: DGN148074Orange Crush
Part: DGN148038Military Digital Camo
Part: DGN1707171Electric Blue
Part: DGN148076Pretty in Pink
Part: DGN148045Silver Bullet
Part: DGN148047Yellow Blaze
Part: DGN148077Faded Glory
Part: DGN1707173Back in Black
Part: DGN148032Deep Purple
Part: DGN148046Midnight Blue
Part: DGN148035Military Green
Part: DGN148037Carbon Fiber
Part: DGN1707169Night Train
Part: DGN1706558Lime Light
Part: DGN148036**12.16. TRU-Balance 3 Power Lift Side Shrouds****REQUIRED with TRU-Balance 3 10" Power Lift selection.** The power lift side shrouds are located at the base of the TRU-Balance 3 system and is available in matching shroud colors. Reference colors shown on page 2.Red Corvette
Part: DGN1704532-PRD1White Satin
Part: DGN1704532-PWH2Grey Street
Part: DGN1704532-PSL3Orange Crush
Part: DGN1704532-POR6Military Digital Camo
Part: DGN1704532-WCM7Electric Blue
Part: DGN1704532-PBU9Pretty in Pink
Part: DGN1704532-PPK3Silver Bullet
Part: DGN1704532-PSL7Yellow Blaze
Part: DGN1704532-PYL3Faded Glory
Part: DGN1704532-WFL2Back in Black
Part: DGN1704532-PBK1Deep Purple
Part: DGN1704532-PPR5Midnight Blue
Part: DGN1704532-PBU10Military Green
Part: DGN1704532-PGN8Carbon Fiber
Part: DGN1704532-WCB1Night Train
Part: DGN1704532-PBK3Lime Light
Part: DGN1704532-PGN7**12.17. Additional TRU-Balance 3 Back Options****Optional.**TRU-Comfort 2 Lateral Positioning Wedge, pair \$30
Part: SET123465

Available with TRU-Comfort 2 backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

See Comfort Company Back order form

Privacy Flap, 9.75" \$35
Part: DWR1044E456B

Standard on recline seating systems.

Privacy Flap, 12.75" \$35
Part: DWR1044E456A

Standard on recline seating systems.

TRU-Balance 3 seating section continued on following page.

13. TRU-BALANCE 3 ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**13.1a. TRU-Balance 3 Flip-Back Armrests**

Choose a left and right TRU-Balance 3 armrest and set height in the following section. An optional downpost is also available. Different left (L) and right (R) armrest height packages are allowed. Not available with Back Cane Option.

☐ L ☐ R TRU-Balance 3 Flip-Back, 8"-12"H Standard
 Part: See XRef

☐ L ☐ R TRU-Balance 3 Flip-Back, 12"-16"H Standard
 Part: See XRef
Armrest Support

TRU-Balance 3 Armrest Downpost (Pair)..... \$87.50 ea.

Part: See XRef

Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available on 14" Width seat or back.

13.1b. Two-Post Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

2-Post, Flip Back, Height Adjustable, and Removable Armrests (9"-13" Height)No Charge

Part: See XRef

Not available with a depth less than 16". Required with flip forward back and lap trays.

13.1c. Heavy Duty, Removable Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$75 ea.
Heavy Duty (Sold as pair)

Part: See XRef. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable, \$75 ea.
Heavy Duty (Sold as pair)

Part: See XRef. HCPCS w/ Static Seat: E0973

13.2. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1. Left Armrest Height**8" H
PTOINDV353510" H
PTOINDV353712" H
PTOINDV353914" H
PTOINDV354116" H
PTOINDV35439" H
PTOINDV353611" H
PTOINDV353813" H
PTOINDV354015" H
PTOINDV3542**2. Right Armrest Height**8" H
PTOINDV345010" H
PTOINDV345212" H
PTOINDV345414" H
PTOINDV345616" H
PTOINDV34589" H
PTOINDV345111" H
PTOINDV345313" H
PTOINDV345515" H
PTOINDV3457**13.3. Armrest Size****REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.
☐ L ☐ R 14" Full Length
 Part: Configuration selection

☐ L ☐ R 10" Desk Length
 Part: Configuration selection
14. TRU-BALANCE 3 ARMPADS**REQUIRED.** Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.No Armpad Selection and/or See Positioning Components Page or Order FormNo Charge
Part: Configuration selection**14.1. Foam Armpads**
☐ L ☐ R Straight, Armpad Standard
 Part: See XRef

☐ L ☐ R Waterfall, Armpad No Charge
 Part: See XRef
14.1b. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for TRU-Balance 3 ArmrestsFalcon Gel, Full Length Left and Right \$185
Part: ARM136417Falcon Gel, Desk Length Left, Full Length Right \$162.50
Part: ARM136418Falcon Gel, Full Length Left, Desk Length Right \$162.50
Part: ARM136426Falcon Gel, Desk Length Left and Right \$140
Part: ARM136416**b. Falcon Gel Armpads for Tube Style Armrests**Falcon Gel, Full Length Left and Right \$296.40
Part: ARM136438Falcon Gel, Desk Length Left, Full Length Right \$274.22
Part: ARM136439Falcon Gel, Full Length Left, Desk Length Right \$274.22
Part: ARM136442Falcon Gel, Desk Length Left and Right \$252.40
Part: ARM136437

15. TRU-BALANCE 3 JOYSTICK BRACKET RECEIVER

REQUIRED.

Standard Fixed Receiver Standard
Part: ACC123525

Stealth Height Adjustable Receiver \$214
Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108
Track Arm Part: ST-ARM250T-Q. HCPCS: K0108



16. TRU-BALANCE 3 HEADRESTS

Optional. Select a headrest pad and bracket.

16.1. Headrest Pads

Stealth 8" Comfort Plus (CP450)..... \$151.94
Part: POS148484. HCPCS: E0955

Stealth 10" Comfort Plus (CP250)(A) \$151.94
Part: SETHEAD1057. HCPCS: E0955

Stealth 14" Comfort Plus (CP150) \$151.94
Part: INDPART2327. HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H) \$151.94

(CPN350)


Part: SETHEAD1069. HCPCS: E0955

Stealth Small Contoured (9.75"Wx3"H)..... \$64.20
(STL-OTSM250)

Part: SETHEAD1052. HCPCS: E0955

Stealth Large Contoured (11"Wx6"H) (B) \$80.25
(STL-OTLG150)

Part: SETHEAD1053. HCPCS: E0955



16.2. Headrest Brackets

Stealth's Worlds Best Removable \$220.50
Headrest Bracket (TWB485-1) (C)
Part: FRMASMB15999. HCPCS: E1028

Stealth Multi-Axis Removable Headrest \$160.50
Bracket (HMO475P) (D)
Part: HDWASMB1084. HCPCS: E1028

17. LAP BELTS

REQUIRED. Select a lap belt below.

17.1a. Standard Lap Belts

38" Lap Belt
Part: ACCBELT1006

50" Lap Belt
Part: ACCBELT1000

60" Lap Belt
Part: ACCBELT1002

70" Lap Belt
Part: ACCBELT1003

80" Lap Belt
Part: ACCBELT1004

90" Lap Belt
Part: ACCBELT1001

17.1b. Optional Lap Belts

i-Fit 60" Single Pull, Push Button, Medium \$86
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

i-Fit 60" Single Pull, Push Button, Extra Large \$86
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

i-Fit 60" Single Pull, Push Button, Large \$86
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

18. TRU-BALANCE 3 MANUFACTURER INSTALLED ACCESSORIES

Optional.

Trilogy Vent Tray \$400
Part: ACC136791. HCPCS: E1029

Push Handles (pair) \$175
Part: ACC147772

Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.

Compact Vent Tray \$1,395
Part: ACC133889. HCPCS: E1029

Accessory Bar \$100
Part: MEC125098

Standard with the Trilogy vent tray, and oxygen holder. May be added to mount additional accessory to back of unit. Not available with back cane option.



19. TRU-BALANCE 3 LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)
Part: Configuration selection

Articulating Foot Platform (Requires 1 Actuator Control)
Part: Configuration selection

Combined Legrests (Requires 1 actuator control)
Part: Configuration selection

20. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options.**20.1a. Tilt Only through Toggle**Single Actuator Control, Through Toggle (Power Tilt Only) **No Charge**

Part: See XRef

Only available with Power Tilt only. Select toggle side. Will require extended lead time when selected with Fender Lights.

Left Side Toggle

Part: PTOINDV1838

Right Side Toggle

Part: PTOINDV1839

20.1b. Q-Logic 3 Controls through JoystickSingle Actuator Through Q-Logic 3 \$1,895
Controller and Interface (AM1)

Part: See XRef. HCPCS: E2310

Multiple Actuator Through Q-Logic 3 \$2,895
Controller and Interface (AAM)

Part: See XRef. HCPCS: E2311

Single Actuator Through Q-Logic 3 \$2,895
Controller and Interface (AAM)

Part: See XRef. HCPCS: E2310

Optional Switch Control

Not available with AM1 electronics. Switch will be mounted to same side of joystick.

4 Way Toggle

Part: ELEASMB5925

4 Way Button Toggle

Part: ELEASMB5926

20.1c. Q-Logic 3 Controls with iAccess Switch Control

iAccess control is not available with combined or independent power legrests.

Single Actuator Control, Through Flat Button \$2,295
iAccess and Q-Logic 3 Controller and Interface (AM1)**COMING SOON.** Part: See XRef. HCPCS: E2310Single Actuator Control, Through Flat Button \$3,295
iAccess and Q-Logic 3 Controller and Interface (AAM)

Part: See XRef. HCPCS: E2310

Multiple Actuator Control Through Flat Button \$3,295
iAccess and Q-Logic 3 Controller and Interface (AAM)

Part: See XRef. HCPCS: E2311

Single Actuator Control, Through Toggle Button \$2,295
iAccess and Q-Logic 3 Controller and Interface (AM1)**COMING SOON.** Part: See XRef. HCPCS: E2310Single Actuator Control, Through Toggle Button \$3,295
iAccess and Q-Logic 3 Controller and Interface (AAM)

Part: See XRef. HCPCS: E2310

Multiple Actuator Control Through Toggle Button \$3,295
iAccess and Q-Logic 3 Controller and Interface (AAM)

Part: See XRef. HCPCS: E2311

iAccess Mount Location

Left Side iAccess

Part: PTOINDV3678

Right Side iAccess

Part: PTOINDV3679

20.1d. Q-Logic 3 Controls with 4 Function Push ButtonSingle Actuator Control Through 4 Function \$2,895
Push Button and Q-Logic 3 Controller and Interface (AAM)

Only available with TB3 armrests. Select a side to mount.

COMING SOON. Part: See XRef. HCPCS: E2310Multiple Actuator Control Through 4 Function \$2,895
Push Button and Q-Logic 3 Controller and Interface (AAM)

Only available with TB3 armrests. Select a side to mount.

COMING SOON. Part: See XRef. HCPCS: E2311**4 Function Push Button Mount Location**

Left Side Mount

Part: PTOINDV3678

Right Side Mount

Part: PTOINDV3679

21. TRU-BALANCE 3 ATTENDANT CONTROL**Optional.** *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.Attendant Control **\$500***

Part: ELE167862. HCPCS: E2331

22. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**22.1. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Power Articulating Foot Platform \$2,995

Part: PTOINDV3754. HCPCS: E1012

Manually Adjustable Center Mount Foot Platform... \$250

Part: RIG163640

2. Lower Extension and Footplates

8"-12" Extension, 5"x6" Footplates

Part: RIG152577

Not available with power articulating platform.

8"-12" Extension, 5"x8" Footplates

Part: RIG152578

Not available with power articulating platform.

10"-16" Extension, 5"x6" Footplates

Part: RIG150206

10"-16" Extension, 5"x8" Footplates with
Width Extension

Part: RIG150203

10"-16" Extension, 5"x12" Footplates with
Width Extension

Part: RIG150212

13"-19" Extension, 5"x6" Footplates

Part: RIG151505

13"-19" Extension, 5"x8" Footplates with
Width Extension

Part: RIG151506

13"-19" Extension, 5"x12" Footplates with
Width Extension

Part: RIG151508

Power Articulating
Foot Platform

Legrests section continued on following pages.

22. LEGRESTS CONTINUED

Center Mount Foot Platform Compatibility Matrix for TRU-Balance® 3	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
10"-16" Length Extension	!	✓	✓	✓	✓	✓	✓	✓	✓
13"-19" Length Extension	!	!	!	!	!	!	!	!	✓

✓ No interference issues at this seat-to-floor height.

! May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H \$150*
Part: RIGASMB7120018

Small Calf Panel 18"Wx8"H \$150
Part: FRMASMB7110012

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$150*
Part: RIGASMB7120019

Large Calf Panel 24"Wx8"H \$200
Part: RIG123824

Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

4" Outside Edge (pair) \$50
Part: RIG152510

Only available with 6"D footplates and 8"D footplate extension.

6" Outside Edge (pair) \$50
Part: RIG152516

Available with 8"D and 11"D footplates and 11"D footplate extension.

2" Bell-Shaped Length Extension \$50
Part: RIG152519

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

5. Additional Accessories

Reversed Dartex Calf Pad Cover (pair) \$150
Small part: RIG133717. Large part: RIG133719

Wheel Bumpers (pair) \$12.50 ea.
Part: WHLASMB1578 x2

Reversed Dartex Panel Cover (pair) \$150
Small part: RIG134673. Large part: RIG134744

Footplate Connector \$15
Part: RIGASMB7120023

22.2. 70° Power Articulating, Swing Away Legrests

4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Articulating Legrest Upper Assembly

P L R 70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 14.5"-18" Extension
Pair Part: FRMASMB9073. Left Part: FRMASMB9074.
Right Part: FRMASMB9075. HCPCS: E1010

P L R 70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 11"-13" Extension
Pair Part: FRMASMB11444. Left Part: FRMASMB11366.
Right Part: FRMASMB11367. HCPCS: E1010

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Angle Adjustable, 11"-13"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	!	!	!	✓

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 4"Wx6"D
Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566.

P L R 5"Wx8"D
Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.

P L R 6"Wx11"D¹
Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570.

P L R 5"Wx6"D
Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567.

P L R 6"Wx8"D¹
Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569.

P L R 8"Wx11"D²
Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571.

3. Power Articulating, Swing-Away Legrest Accessories

Calf Pad (pair) 6"Wx6"H No Charge
Part: MEC153735

Falcon Gel Padded Knee Adductor Buttons (KBE)(pair) \$125
Part: INDPART2242. HCPCS: E0953[^]

[^]Please verify the appropriate billing code to be used for this item with all third-party payors.

Legrests section continued on following pages.

22. LEGRESTS CONTINUED**22.3. 60°-90° Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Swing-Away Legrest Upper

60°-90° Swing-Away Legrest.....\$81 ea. or N/C*
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559

*No charge with Power Positioning. Not available with 11"-13" extensions.

60°-90° Swing-Away Legrest, Short.....\$81 ea. or N/C*
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.

*No charge with Power Positioning. For use only with 11"-13" lower extensions.

2. Lower Extension

11"-13" Length
 Part: See XRef
 N/A with Multi-axis.

13"-16.5" Length
 Part: See XRef
 N/A with Multi-Axis.

11"-16.5" Length
 Part: See XRef
 N/A with Angle Adjustable.

15.5"-19.5" Length
 Part: See XRef

3a. Extruded Footplate(s)

7.25"Wx6"D Extruded Footplate..... Standard
 Part pair: FRMASMB7847. Left: FRMASMB7861
 Right: FRMASMB7862

8.5"Wx6"D Extruded Footplate..... Standard
 Part pair: FRMASMB7848. Left: FRMASMB7863
 Right: FRMASMB7864
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

4"Wx6"D.....\$85 ea.
 Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
 HCPCS: K0040

6"Wx8"D¹.....\$85 ea.
 Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
 HCPCS: K0040

5"Wx6"D.....\$85 ea.
 Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
 HCPCS: K0040

6"Wx11"D¹.....\$85 ea.
 Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
 HCPCS: K0040

5"Wx8"D.....\$85 ea.
 Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562
 HCPCS: K0040

8"Wx11"D².....\$85 ea.
 Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565
 HCPCS: K0040

3c. Multi-Axis Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

4"Wx6"D.....\$247.50 ea.
 Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
 HCPCS: K0040

6"Wx8"D¹.....\$247.50 ea.
 Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
 HCPCS: K0040

5"Wx6"D.....\$247.50 ea.
 Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
 HCPCS: K0040

6"Wx11"D¹.....\$247.50 ea.
 Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
 HCPCS: K0040

5"Wx8"D.....\$247.50 ea.
 Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554
 HCPCS: K0040

8"Wx11"D².....\$247.50 ea.
 Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559
 HCPCS: K0040

4. 60°-90° Swing-Away Legrest Accessories

Heel Loops (pair).....\$50
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi-axis footplates.

Wheel Bumper.....\$12.50 ea.
 Part: See XRef

Gel Neoprene Footrest Sleeve (pair).....\$80
 Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support.....\$230 ea.
 Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover.....\$155 ea.
 Part: INDPART2823. HCPCS: K0108

QTY

Legrests section continued on following page.

22. LEGRESTS CONTINUED**22.4. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$162.50 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2. Lower Extension

☐ ☐ ☐ 11"-14" Length ☐ ☐ ☐ 14"-17" Length
 Pair: FRMASMB4755 Pair: FRMASMB4756
 Left: RIG138026 Left: RIG138033
 Right: RIG138030 Right: RIG138035
 For reinforced only For reinforced only

☐ ☐ ☐ 12"-14.5" Length ☐ ☐ ☐ 14.5"-18" Length
 Part: See XRef Part: See XRef
 For angle adjustable only For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-14.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	✓	!	✓	✓

☒ This rigging configuration causes no interference issues at this seat-to-floor height.
☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s) with Heel Loops

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 4"Wx6"D \$85 ea.
 Pair parts: FRMASMB8812 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx8"D¹ \$85 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040

☐ ☐ ☐ 5"Wx6"D \$85 ea.
 Pair parts: FRMASMB8811 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx11"D¹ \$85 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040

☐ ☐ ☐ 5"Wx8"D \$85 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040

☐ ☐ ☐ 8"Wx11"D² \$85 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040

4. Elevating Legrest Accessories

Heel Loops (pair) \$50
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable footplates.

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$125
 Part: INDPART2242. HCPCS: E0953[^]
[^]Please verify the appropriate billing code to be used for this item with all third-party payors.

Wheel Bumper \$12.50 ea.
 Part: See XRef
 Available on angle adjustable footplates.

23. SET LEGREST LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

24. QUICK STEALTH CUSHION CONFIGURATION

Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the complete [Stealth Cushion Order Form](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

24.1. Cushion Dimensions**Cushion Width**

14" W 16" W 17" W 18" W 19" W 20" W 22" W

Cushion Depth


14" D 16" D 17" D 18" D 19" D 20" D 22" D

24.2. Cushion Model Selection**SIMPLICITY G CUSHION**

General use with 2-way stretch, breathable cover.

Part series: ST-SIMGXXXX

E2601, 16"W to 20" W \$110
E2602, 22"W \$165

**SOLUTION SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SOLSPPXXXX

E2607, 14"W to 20" W \$375
E2608, 22"W \$430

**TRU-COMFORT 2 SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover.

Part series: ST-TC2SPPXXXX


E2607, 16"W to 20" W \$400
E2608, 22"W \$500

**SPECTRUM GEL SPP CUSHION**

Positioning with 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCGSPPXXXX


E2607, 14"W to 20" W \$460
E2608, 22"W \$530

**SPECTRUM FOAM SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover.

Part series: ST-SPCFSPXXXX

E2607, 16"W to 20" W \$460
E2608, 22"W \$530

**24.3. Additional Fluid Proof Liners**

Fluid Proof Urethane Liner \$50-\$70

Part: URL1709070-XXXX

Not available with Simplicity and TRU-Comfort 2.

Fluid Proof Polyurethane Liner Coming soon

Not available with Simplicity.

Liner Installation

Install Liner on Cushion

25. PROVIDER INSTALLED ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Cup Holder No Charge

Part: ACC125003

For use with TRU-Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest) No Charge

Part: ACCASMB2634

For use with tube-style armrests.

Backpack Holder (Pair) (A) \$50

Part: ACC147774

Not available with back cane option.

Hydration System (B) \$200

Part: ACC125010

Glove Box (C) \$125

Part: ACC125011

Transfer Bar, Quick Release (Pair) (D) \$200

Part: ACC147773. HCPCS: E0705

Front Mounted Cell Phone/Media Holder (E) \$50

Part: ACC155835

Large, Black, Weather Cover \$168

Part: ACCCOVER1021

Reflector Kit \$50

Part: ACC124670

Personal Item Hook \$50

Part: ACC125008

Not available with back cane option.

Rear View Mirror \$50

Part: ACC155833

Oxygen Holder \$200


Part: ACC125006. HCPCS: E2208

Clothing Guard (Pair) \$200

Part: ACC125007

Quantum Backpack \$75

Part: ACCINDV1102

TRU-Balance 3
Accessories**26. BATTERIES****REQUIRED.** The powerchair requires 2 flat top batteries and will be installed if ordered.

Group 24, AGM Cell Battery (MK) \$410 ea.

Part: BAT1712207. HCPCS: E2363

No batteries. Provider will assume responsibility of batteries.

Part: Configuration selection

Group 34, Gel Cell Battery (MK) \$395 ea.

Part: BATGEL1006. HCPCS: E2359

27. BRANCH PATCHES**Optional.**Air Force Patch \$20
Part: ACC140240Army Patch \$20
Part: ACC140241Coast Guard Patch \$20
Part: ACC140239Marines Patch \$20
Part: ACC140243Navy Patch \$20
Part: ACC140244Purple Heart Patch \$20
Part: ACC140239**28. MANUFACTURER INSTALLED POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**28.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, Rail Mounts, pair (TWBTLTC) \$472.50
Part: INDPART2677. HCPCS: E1028**2. Pads**Curved Pads, 3"x5", pair (TWBL-3X5) \$210
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6) \$210
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**28.2. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide BracketFlip Down Bracket (TWBTLADD-HW) \$236.25 ea.
Part: INDPART2804. HCPCS: E1028Removable Bracket (LP430) \$262.50 ea.
Part: INDPART2805. HCPCS: E1028Removable Bracket with 1" Offset \$400 ea.
(LP430 & LHW-121)
Part: INDPART2854. HCPCS: E1028Flip Down with LEFT Full Surface Contact \$236.25 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028Flip Down with RIGHT Full Surface Contact \$236.25 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028**2. Pelvic/Thigh Guide Adapter**2" Height Adapter Plate (LHW-122) \$26.25
Part: INDPART2821. HCPCS: K0108
Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.2 1/4" Adapter Plate (LHW-129) \$53.50
Part: POS140430. HCPCS: K0108
Only available with the flip down full surface bracket.**3. Pelvic/Thigh Guide Pad**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4) \$75 ea.
Part: INDPART2801. HCPCS: E0953^Medium Pad, 4"x8" (TWBADD-4X8) \$75 ea.
Part: INDPART2802. HCPCS: E0953^Long Pad, 4"x12" (TWBADD-4X12) \$75 ea.
Part: INDPART2803. HCPCS: E0953^**4. Pelvic/Thigh Guide Installation**Install on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719**29. PROVIDER INSTALLED POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**29.1. Stealth Foot Boxes**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$498.75
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954^Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954^Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$498.75
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954^

Provider installed positioning components continued on following page.

29. PROVIDER INSTALLED POSITIONING COMPONENTS CONTINUED**29.2. Stealth Padded Shoe Holder**

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
(SPLE111) (pair)
Part: INDPART2783

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266
(SPLE114) (pair)
Part: INDPART2785

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$266
(SPLE112) (pair)
Part: INDPART2784


29.3. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support \$183.75 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support \$183.75 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

Left, Fixed Elbow Support Bracket \$183.75 ea.
TRU-Balance 3 Arms (ARMS-ESFTB3-L)
Part: POS133465

Right, Fixed Elbow Support Bracket \$183.75 ea.
TRU-Balance 3 Arms (ARMS-ESFTB3-R)
Part: POS133466



Elbow Supports

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$52.50 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$52.50 ea.
Part: POS124067. HCPCS: K0108

29.4. Stealth i-Fit Chest Harness

Upper Classic, Medium (IFPS1103) \$142 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103) \$142 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960



Contoured and Classic Chest Harness

29.5. Stealth i-Fit Chest Straps

Chest Strap, Medium w/ Belt Anchor \$90 ea.
(IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$90 ea.
(IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960



Chest Strap

29.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

Facial Lateral for Comfort Series \$294.25 ea.
Right Side (CPS9-R-HW)
Part: INDPART2647. HCPCS: E1028

Facial Laterals for Comfort Series \$347.75 ea.
Left and Right Side (CPS9-B-HW)
Part: INDPART2649. HCPCS: E1028

Facial Lateral for Comfort Series \$294.25 ea.
Left Side (CPS9-L-HW)
Part: INDPART2648. HCPCS: E1028



Left and Right Facial Lateral with Pads

Facial Lateral Pads and Accessories

1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$80.25
Part: POS150532
Required with facial laterals to mount to comfort headrest brackets.

Facial Pad, 4"x2" (900-Q) \$53.50 ea.
Part: ST-900-Q. HCPCS: K0108


Facial Pad, 6"x2" (920-Q) \$53.50 ea.
Part: ST-920-Q. HCPCS: K0108

Egg Switch Mount (ESM) (A) \$52.50 ea.
Part: INDPART2815

Swing-Away Joystick Mount, 18" Rod (C) \$250 ea.
Part: ST-SUS18-Q

Figure 8 Bracket with 13" Rod \$267.50 ea.
and Egg Switch Mount (B)
Part: ST-F8675-13-Q

Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

30. SPECIALTY CONTROLS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**30.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,720
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard
Part: DGN137728**Match Joystick Shroud to Base Color** \$50
Part: See color below**Non-matching Joystick Shroud**..... \$50
Select color below**Red Corvette**
Part: DGN137729-PRD1**White Satin**
Part: DGN137729-PWH2**Grey Street**
Part: DGN137729-PSL3**Orange Crush**
Part: DGN137729-POR6**Military Digital Camo**
Part: DGN137729-WCM7**Electric Blue**
Part: DGN137729-PBU9**Pretty in Pink**
Part: DGN137729-PPK3**Silver Bullet**
Part: DGN137729-PSL7**Yellow Blaze**
Part: DGN137729-PYL3**Faded Glory**
Part: DGN137729-WFL2**Back in Black**
Part: DGN137729-PBK1**Deep Purple**
Part: DGN137729-PPR5**Midnight Blue**
Part: DGN137729-PBU10**Military Green**
Part: DGN137729-PGN8**Carbon Fiber**
Part: DGN137729-WCB1**Night Train**
Part: DGN137729-PBK3**Lime Light**
Part: DGN137729-PGN7

Enhanced Display

30.2. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ**.....\$700
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick..



EX SAJ

30.3. Head Array**Adult, i-Drive Tri-Array Series, 3-Switch** \$4,414
Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q, HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461
(IDHBT300-1)

Part: ST-IDHBT300-1-Q, HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

30.4. Sip and Puff Drive**Module****Sip and Puff Module for Q-Logic 3**.....\$1,895
Part: CTL168056

Includes module and mounting hardware.



Sip & Puff Module

Straw and Mount**Stealth Sip and Puff for Comfort Series,**\$294.25
Left Mount (SP625-CL)

Part: ST-SP625-CL-Q

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Comfort Series,\$294.25
Right Mount (SP622-CR)

Part: ST-SP625-CR-Q

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive,\$294.25
Head Array (SP625-U)

Part: ST-SP625-U-Q

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Stealth Sip & Puff

30.5. Sip and Puff Drive**Red Egg Switch (ESRED)**..... \$80.85
Part: SWTMCRO1077**Blue Egg Switch (ESBLU)** \$80.85
Part: SWTMCRO1079**Black Egg Switch (ESBLK)** \$80.85
Part: SWTMCRO1081**Green Egg Switch (ESGRN)** \$80.85
Part: SWTMCRO1078**Yellow Egg Switch (ESYEL)** \$80.85
Part: SWTMCRO1080

Egg Switch

31. TRU-BALANCE 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA)..... No Charge

Part: Configuration selection

Secondary Multiplier (SM)..... \$221.21

Part: HAR144004

Special Order Configuration (SOF)

Part: Configuration selection

Arm Mounts		Back Mounts			
		1 Arm	2 Arm	3 Arm	4+ Arm
Arm Mounts	Joystick (Hand Control)	NA	SM	SM	SOF
	Enhanced Display	NA	SM	SM	SOF
	iAccess	NA	SM	SM	SOF
	USB Charger	SM	SM	SM	SOF
	SAJ	SOF	SOF	SOF	SOF
Back Mounts	Attendant Control	NA	SM	SOF	SOF
	SCIM	NA	SM	SOF	SOF
	Sip n Puff	SM	SM	SOF	SOF
	ECU	SOF	SOF	SOF	SOF
Power Positioning Chart					
		Static Seat Chart			
		1 Arm	2 Arm	3 Arm	4+ Arm
0 Back	NA	SM	SOF	SOF	SOF
	NA	SM	SOF	SOF	SOF
	SM	SM	SOF	SOF	SOF
	SOF	SOF	SOF	SOF	SOF
Q-Logic 3 Multiplier Charts					

32. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2018 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.



[youtube.com/
QuantumRehab](https://www.youtube.com/QuantumRehab)



[facebook.com/
QuantumRehab](https://www.facebook.com/QuantumRehab)



[twitter.com/
quantumrehab](https://twitter.com/quantumrehab)



[instagram.com/
quantum_rehab](https://www.instagram.com/quantum_rehab)