

EDGE 3

EDGE 3 shown in Blue Jean
with iLevel® power seat elevation
and Q-Logic 3 with Bluetooth®

AVAILABLE WITH
12" iLEVEL® POWER
SEAT ELEVATION AT
4.5 MPH



PROUDLY MADE IN
AMERICA 
OF FOREIGN AND DOMESTIC PARTS



quantumrehab.com


QUANTUM®

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 6 mph Up to 4.5 mph w/iLevel activated
GROUND CLEARANCE	2.7" (battery tray)
TURNING RADIUS ²	20.5"
OVERALL LENGTH ²	32.23" without front riggings
BASE WIDTH	24.4"
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 10-22" Depth: 10-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-22" Depth: 12-22"
TRU-Balance® Flex Seating	Width: 15-22" Depth: 12-20"
SEAT-TO-FLOOR HEIGHTS ¹⁰	
Synergy®/TRU-Balance® 3 Static Seat ⁵	16.5-18.5"
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5-19.5"
TRU-Balance® Flex Seating	16.5" - 19.5"
MANUAL RECLINE ⁶	Yes
BATTERY SIZE ³	NF-22 (2)
BATTERY WEIGHT ⁷	38 lbs. each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6 [®]
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
 8) All Quantum Rehab® units can be ordered with specialty controls.
 9) Call for availability
 10) 16.5" seat-to-floor height is not available with power tilt or iLevel®. 19.5" seat-to-floor height is not available with static seat or power recline.

Features


- Upgraded SRS (Smooth Ride Suspension)
- 4.5 mph with iLevel®
- Fender lights and USB charger included
- Bold new colors
- Made in the USA of foreign and domestic parts



7 New Colors*



HYDROGRAPHICS



* Visit quantumrehab.com to see all 19 available colors.