

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com**Q6 Edge® 3 Order Form with Contoured, Synergy®, and TRU-Balance® Flex Seating**

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



TRU-Balance® 3 Power Lift Seating Shown

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


Thank you for choosing Quantum!**PATIENT INFORMATION**

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

**1. BASE MODEL**

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

☐ Quantum Q6 Edge 3 3S-SS \$7,895

Part: Q6 EDGE 3 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

☐ Quantum Q6 Edge 3 3SP-SS \$7,895

Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856

Single power option for power elevating seat.

☐ Quantum Q6 Edge 3 3S-C \$7,895

Part: Q6 EDGE 3 3S-C. HCPCS: K0849

Contoured seat (captain) standard.

☐ Quantum Q6 Edge 3 3SP-C \$7,895

Part: Q6 EDGE 3 3SP-C. HCPCS: K0857

Contoured seat (captain), single power option for 12" lift.

Account Number: _____ Date: _____ Marked for: _____

2. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Patient's Weight: _____

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1801101-PRD1	 White Satin Part: DGN1801101-PWH2	 Electric Blue Part: DGN1801101-PBU9	 Military Green Part: DGN1801101-PGN8
 Cherry Bomb Part: DGN1801101-PRD4	 Deep Purple Part: DGN1801101-PPR5	 Blue Jean Part: DGN1801101-PBU11	 Sugar Plum Part: DGN1801101-PPR7
 Back in Black Part: DGN1801101-PBK1	 Raspberry Beret Part: DGN1801101-PPK4	 Lemon Crush Part: DGN1801101-PYL4	To Be Determined Follow-up required.
 Night Train Part: DGN1801101-PBK3	 Evolution Orange Part: DGN1801101-POR7	 Green Machine Part: DGN1801101-PGN11	
 Grey Street Part: DGN1801101-PSL3	 Midnight Blue Part: DGN1801101-PBU10		

3.1b. Hydrographic Shrouds

☐  Military Digital Camo \$250 Part: DGN153201-WCM7	☐  Faded Glory \$250 Part: DGN153201-WFL2	☐  Carbon Fiber \$250 Part: DGN153201-WCB1
---	---	--

4. DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard Part: MOT147433	Accu-Trac Motor Technology \$1,295 Part: MOT147881. HCPCS: K0108 Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.
---	---

5. BASE MOUNTED FOOT PLATFORM

REQUIRED. May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard Part: FRMASMB16077	Omit base mounted foot platform No Charge Part: PTOINDV3674
Bell Shaped, Base Mounted Foot Platform \$235 Part: RIGASMB7120041	

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires Standard Part: PTOINDV3525 (pair)	3"x14" Pneumatic Black Tires No Charge Part: WHL1713925 x2
4"x14" Knobby, Flat Free Black Tires \$275.50 ea. Part: WHL147601 x2	4"x14" Knobby, Pneumatic Black Tires \$275.50 ea. Part: WHL147599 x2

6.2. Drive Wheel Accent Color

Not available with knobby tires.

Back in Black Standard Part: DGN164245-PBK1 x2	Raspberry Beret No Charge Part: DGN164245-PPK4 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Red Corvette No Charge Part: DGN164245-PRD1 x2	Evolution Orange No Charge Part: DGN164245-POR7 x2	Sugar Plum No Charge Part: DGN164245-PPR7 x2
Cherry Bomb No Charge Part: DGN164245-PRD4 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
Night Train No Charge Part: DGN164245-PBK3 x2	Electric Blue No Charge Part: DGN164245-PBU9 x2	Faded Glory \$50 ea. Part: DGN164245-WFL2 x2
Grey Street No Charge Part: DGN164245-PSL3 x2	Blue Jean No Charge Part: DGN164245-PBU11 x2	Carbon Fiber \$50 ea. Part: DGN164245-WCB1 x2
White Satin No Charge Part: DGN164245-PWH2 x2	Lemon Crush No Charge Part: DGN164245-PYL4 x2	To Be Determined Follow-up required.
Deep Purple No Charge Part: DGN164245-PPR5 x2	Green Machine No Charge Part: DGN164245-PGN11 x2	

7. SEAT TOWERS/POWER ADJUSTABLE SEAT HEIGHT**REQUIRED.** Select one option below.**7.1a. Power Base Seat Towers**

Standard Seat Towers Standard
 Part: Configuration selection, no part number

Extended Seat Towers (increased seat to floor height) ... No Charge
 Part: MEC120301

7.1b. Power Adjustable Seat Height

The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit.

12" Power Lift for Captain Seat with iLevel.....\$2,600
 Part: PTOINDV3857. HCPCS: E2300
 For use with K0857 base code. Must select Power Lift Captain Seat in contoured seat section. Includes armrests and lap belt. Available with swing away legrests or manually adjustable center mount foot platform.

Power Elevating Seat with 7" Stroke \$1,895
 Contoured seat part: FRMASMB16310. HCPCS: E2300
 Comfort contoured seat part: MEC133630. HCPCS: E2300
 Synergy Seat part: FRMASMB16312. HCPCS: E2300

Power Elevating Seat with 6" Stroke \$1,895
 Contoured seat part: FRMASMB16309. HCPCS: E2300
 Comfort contoured seat part: MEC133629. HCPCS: E2300
 Synergy Seat part: FRMASMB16311. HCPCS: E2300

8. CONTOURED SEAT, SEAT TO FLOOR HEIGHT**REQUIRED** with contoured seat selection later in form.

Lowest Seat to Floor
 Part: PTOINDV3558

Middle Seat to Floor
 Part: PTOINDV3559

Highest Seat to Floor
 Part: PTOINDV3560

9. FENDER LIGHTS**Optional.** Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

☐ Switch Fender Lights No Charge
 Part: DGN154117

Fender Light Switch Mount Type


Side Rail/Track Mount Bracket Standard
 Part: MEC1800866
 Allows the egg switch to be mounted to seating side rails using a t-nut assembly. This may also be mounted to other armrest, seat or accessory bar rails. Not available with contoured seating.

Under Arm Mounting Bracket \$25*
 Part: MEC1802726
 *No charge with contoured seating. Allows the egg switch to mount to a joystick tube-style bracket mounted to any armrest. Includes joystick mounting bracket.

10. FULL LIGHT PACKAGES**Optional.** Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package. Select one light package below.

☐ Base Mounted, LED, Full Light Package \$995
 Part without iLevel: ELE1705999. Part with iLevel: ELE168643
 Fender lights will be omitted with Base Mounted, LED, Full Light Package.

☐ Seat Mounted, LED, Full Light Package \$995
 Flex part: ELE131477. Synergy part: ELEASMB5719
 Not available with contoured seating.

**11. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**11.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

☐ NE 4-Key Joystick Standard
 Part: CTLDC1554
 Non available with power seating.

☐ NE+ 6-Key Joystick (2 Actuators) \$200
 Part: CTLDC1560
 Control up to 2 actuators through joystick.

☐ NE+ 8-Key Joystick (2 Actuators) \$400
 Part: CTLDC1561
 Control up to 2 actuators through joystick and full light package.

☐ No Joystick, Setup for NE+ Joystick No Charge
 Part: Configuration selection

**11.1b. Expandable Controls - Quantum Intelligent Systems™**Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. **Q-Logic 3 controls are not available on the TRU-Balance Flex seating.****1. Q-Logic 3 Required Components**

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

☐ Q-Logic 3 EX Controller Module \$695
 Part: See XRef. HCPCS: E2377

☐ Q-Logic 3 Expandable Controller Harness \$425
 Part: HAR144005. HCPCS: E2313

2. Q-Logic 3 Joystick

☐ Q-Logic 3 EX Joystick \$700
 Part: CTL167724

☐ Q-Logic 3 EX Joystick for Full Lights \$700
 Part: CTL167725

☐ No Joystick, Setup for Q-Logic
 Part: Configuration selection.



Electronics section continued on following page.

11. ELECTRONICS CONTINUED**11.2. Q-Logic 3 Joystick Armor****REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2.**Black, Soft Touch Standard
Part: DGN124307Non-matching Armor \$50
Part: PTOINDV3422
Select color belowMatch Armor to Base Color \$50
Part: PTOINDV3421
Part below will be selected with system configurator.Red Corvette
Part: DGN124336-PRD1White Satin
Part: DGN124336-PWH2Electric Blue
Part: DGN124336-PBU9Sugar Plum
Part: DGN124336-PPR7Cherry Bomb
Part: DGN124336-PRD4Deep Purple
Part: DGN124336-PPR5Blue Jean
Part: DGN124336-PBU11Military Digital Camo
Part: DGN124336-WCM7Back in Black
Part: DGN124336-PBK1Raspberry Beret
Part: DGN124336-PPK4Lemon Crush
Part: DGN124336-PYL4Faded Glory
Part: DGN124336-WFL2Night Train
Part: DGN124336-PBK3Evolution Orange
Part: DGN124336-POR7Green Machine
Part: DGN124336-PGN11Carbon Fiber
Part: DGN124336-WCB1Grey Street
Part: DGN124336-PSL3Midnight Blue
Part: DGN124336-PBU10Military Green
Part: DGN124336-PGN8To Be Determined
Follow-up required.**11.3. SCIM Module****Optional.** Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.SCIM Module for Q-Logic 3 \$525
Part: CTL168055**12. JOYSTICK HANDLES****Optional.** The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.☐ Mushroom Handle (B) \$85
Part: ACCASMB1860. HCPCS: E2323☐ Extended I-Handle (G) \$85
Part: ACCASMB1857. HCPCS: E2323☐ Chin Cup (C) \$85
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)☐ Bodypoint 3" Goal Post Handle \$99
Part: PLSKNCP1057. HCPCS: E2323☐ T-Bar (D) \$120
Part: ACCASMB1858. HCPCS: E2323☐ Bodypoint 4" Goal Post Handle \$99
Part: PLSKNCP1058. HCPCS: E2323☐ Medium Ball Handle (E) \$65
Part: ACCASMB1859. HCPCS: E2323☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308☐ Soft Ball Handle (F) \$65
Part: ACCASMB1856. HCPCS: E2323**13. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.


Omit Joystick Mounting Bracket

13.1a. Fixed, In-line Joystick Bracket☐ Fixed, In-line, Bracket Standard
Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744. SAJ Part: FRMASMB10454☐ Left Side Mount
Part: PTOINDV1225☐ Right Side Mount
Part: PTOINDV1226**13.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

☐ Swing-away, Inline, Left Side Mount \$245
Part: See XRef. HCPCS: E1028☐ Swing-away, 1" Offset, Left Side Mount \$245
Part: See XRef. HCPCS: E1028☐ Swing-away, Inline, Right Side Mount \$245
Part: See XRef. HCPCS: E1028☐ Swing-away, 1" Offset, Right Side Mount \$245
Part: See XRef. HCPCS: E1028**13.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

☐ Retract4 Mount \$350
NE/NE+ Part: MEC138876. HCPCS: E1028
Q-Logic Part: MEC138860. HCPCS: E1028
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.☐ Stealth Swing-Away/Flip Down Bracket \$294.25
NE/NE+ Part: ST-TWBM480QL-Q. HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028
SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.☐ Multi-axis Mount \$100
NE/NE+ Part: ACC142528. HCPCS: E1028
Q-Logic Part: ACC131216. HCPCS: E1028
SAJ Part: ACC142740. HCPCS: E1028
Allows height adjustment of joystick.**Specialty Mount Side**☐ Left Side Mount
Part: PTOINDV1225☐ Right Side Mount
Part: PTOINDV1226

Account Number: _____ Date: _____ Marked for: _____

14. ELECTRONICS ACCESSORIES

Optional.

☐ Q-Logic EX/NE Programmer \$545
Part: ELEASMB6976

☐ Q-Logic PC Programming Disk with Cable \$750
Part: ELEASMB5216

15. CONTOURED SEATING

Optional. Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

15.1. Occupied Transit, WC-19

WC-19 Occupied Transit Option for Contoured Seat \$350
Part: KIT122609

Only available with Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

15.2a. Captain Style Seating

1a. Depth Adjustable High Back with Limited Recline

Full length armrests must be selected with a depth adjustable high back.

18"Wx18"-20"D with 60" Lap Belt No Charge
Part: SETASMB6895

22"Wx20"-22"D with 70" Lap Belt \$150
Part: SETASMB6897

20"Wx18"-20"D with 60" Lap Belt \$150
Part: SETASMB6896

1b. Power Lift, Fold Flat Captain Seat

Only available with iLevel. Includes armrests and lap belt.

16"Wx16"-18"D with 60" Lap Belt \$150
Part: SET1701577

20"Wx20"-22"D with 60" Lap Belt \$150
Part: SET1701579

18"Wx18"-20"D with 60" Lap Belt \$150
Part: SET1701578

1c. High Back with Limited Recline

16"Wx16"D \$150
Part: SETASMB2770

16"Wx22"D \$150
Part: SETASMB4887

20"Wx16"D \$150
Part: SETASMB2780

16"Wx18"D \$150
Part: SETASMB2771

18"Wx16"D \$150
Part: SETASMB2775

20"Wx22"D \$150
Part: SETASMB4889

16"Wx20"D \$150
Part: SETASMB2772

18"Wx22"D \$150
Part: SETASMB4888

1d. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline

Full length armrests must be selected with a depth adjustable high back.

16"Wx16"-18"D with 60" Lap Belt \$150
Part: SET131843

20"Wx18"-20"D with 60" Lap Belt \$150
Part: SET131845

18"Wx18"-20"D with 60" Lap Belt \$150
Part: SET131844

15.2b. Contoured Seat with Solid Seat Pan

1a. High Back with Limited Recline and Solid Seat Pan

16"Wx16"D \$150
Part: SETASMB3178

18"Wx16"D \$150
Part: SETASMB3181

20"Wx16"D \$150
Part: SETASMB3183

16"Wx18"D \$150
Part: SETASMB3179

18"Wx18"D \$150
Part: SETASMB2848

20"Wx18"D \$150
Part: SETASMB3184

16"Wx20"D \$150
Part: SETASMB3180

18"Wx20"D \$150
Part: SETASMB3182

20"Wx20"D \$150
Part: SETASMB3185

16"Wx22"D \$150
Part: SETASMB5006

18"Wx22"D \$150
Part: SETASMB5007

20"Wx22"D \$150
Part: SETASMB5008

1b. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline and Solid Seat Pan

Full length armrests must be selected with a depth adjustable high back.

16"Wx14"-16"D with 60" Lap Belt \$150
Part: SET131846

18"Wx16"-18"D with 60" Lap Belt \$150
Part: SET131847

16"Wx16"-18"D with 60" Lap Belt \$150
Part: SET145814

18"Wx18"-20"D with 60" Lap Belt \$150
Part: SET145815

15.3. Depth Selection for Comfort Seat

Required with Comfort Seat selection above.

14" Depth 16" Depth 18" Depth 20" Depth

Part: SETDPTH1004 Part: SETDPTH1006 Part: SETDPTH1008 Part: SETDPTH1010

15.4. Height Selection for Comfort Seat

Required with Comfort Seat selection above.

Set Back Height to Lowest Position Set Back Height to Highest Position

Part: SETHGHT1100 Part: SETHGHT1101

15.5. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: _____ Width: _____ Depth: _____

Contoured seating section continued on following page.

15. CONTOURED SEATING CONTINUED**15.6. Contoured Seat Flip-Up, Height Adjustable, Removable Armrests****REQUIRED except with fold flat, power lift captain seat.** Select a left and right armrest or to omit them.Full Length, Left Armrest \$85
Part: FRMASMB8614. HCPCS: E0973Desk Length, Left Armrest \$85
Part: FRMASMB8612. HCPCS: E0973Full Length, Right Armrest \$85
Part: FRMASMB8615. HCPCS: E0973Desk Length, Right Armrest \$85
Part: FRMASMB8613. HCPCS: E0973Flip-up, Height Adjustable, Comfort Right Armrest \$85
Part: FRMASMB16315. HCPCS: E0973Flip-up, Height Adjustable, Comfort Left Armrest \$85
Part: FRMASMB16316. HCPCS: E0973No Armrest Standard
Part: N/A**16. SYNERGY SEATING SYSTEM****Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 32"D.**16.1. Synergy Frame Selection**Synergy Seating System \$434
Part: PTOINDV2995**16.2. Seat Dimension - Width****REQUIRED with Synergy seating.** Part number is SETWDTH plus number following # below.

10" W (#1000)	12" W (#1002)	14" W (#1004)	16" W (#1006)	18" W (#1008)	20" W (#1010)	22" W (#1012)
11" W (#1001)	13" W (#1003)	15" W (#1005)	17" W (#1007)	19" W (#1009)	21" W (#1011)	

16.3. Seat Dimension - Width Package**REQUIRED with Synergy seating.** Maximum seat width is 22".

10"-14" Width Package Part: FRMASMB18011	14"-18" Width Package Part: FRMASMB18012	16"-20" Width Package Part: FRMASMB18013	20"-24" Width Package Part: FRMASMB18014 (Maximum width is 22")
---	---	---	---

16.4. Seat Dimension - Depth**REQUIRED with Synergy seating.** Part number is SETDPTH plus number following # below.

10" D (#1000)	12" D (#1002)	14" D (#1004)	16" D (#1006)	18" D (#1008)	20" D (#1010)	22" D (#1012)
11" D (#1001)	13" D (#1003)	15" D (#1005)	17" D (#1007)	19" D (#1009)	21" D (#1011)	

16.5. Seat Dimension - Depth Package**REQUIRED with Synergy seating.** Maximum seat depth is 22".

10"-14" Depth Package Part: FRMASMB6981	14"-18" Depth Package Part: FRMASMB17823	15"-20" Depth Package Part: FRMASMB7732	20"-24" Depth Package Part: FRMASMB7733 (Maximum depth is 22")
--	---	--	--

15.6. Synergy Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.☐ WC-19 Occupied Transit Option for Synergy Seat \$275
Part: KITASMB1830To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>**16.7. Synergy Back Options****REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Selection**☐ Synergy Back \$275
Part: See XRef

Select a back cane type and height and back height.

☐ TRU-Comfort, 20" Width or less \$675
Part: See XRef. HCPCS: E2620

Select back height.

☐ TRU-Comfort, 22" Width or greater \$775
Part: See XRef. HCPCS: E2621

Select back height.

☐ Sling Back \$98
Part: See XRef

Select a back height. Available 18"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.

☐ Omit Back No Charge
Part: PTOINDV1287

Select a back cane type and height in following section.

		Width									
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	15"	S	S	S	-	-	-	-	-	-	-
	16"	S	S	S	-	-	-	-	-	-	-
	17"	-	S	S	S	S	S	S	S	S	S
	18"	-	S	S	S	S	S	S	S	S	S
	19"	-	S	S	S	S	S	S	S	S	S
	20"	-	S	S	S TC	S	S TC	S	S TC	S	S
	21"	-	S	S	S TC	S	S TC	S	S TC	S	S
	22"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	23"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	24"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC
	25"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC
S Synergy Back		S TC Synergy or TRU-Comfort Back									

Synergy and TRU-Comfort Back Availability Chart

Synergy seating section continued on following pages.

16. SYNERGY SEATING SYSTEM CONTINUED**2. Back Cane Selection**

- ☐ Synergy Back Cane No Charge
Part: FRMASMB15845
Select a back cane height.
- ☐ Flip Forward Back Cane \$150
Part: FRMASMB15846
Select a back cane height. Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.
- ☐ 2" Mesomorph (1" per side) Back Cane \$195
Part: FRMASMB12067
Select a cane height.
- ☐ 4" Mesomorph (2" per side) Back Cane \$195
Part: FRMASMB12091
Select a cane height.
- ☐ 2" Endomorph (1" per side) Back Cane \$195
Part: FRMASMB12115
Select a cane height.
- ☐ 4" Endomorph (2" per side) Back Cane \$195
Part: FRMASMB12139
Select a cane height.
- ☐ Manual Recline Back \$750
Part: FRMASMB12910. HCPCS: E1226
Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.

3. Synergy Push Handles

- ☐ Push Handles Standard
Part: Configuration selection with back cane height
- ☐ Omit Push Handles \$125
Part: Configuration selection with back cane height
Not available with Sling Back. Minimum cane height is 16".

4. Synergy Accessory Bar

- ☐ Accessory Bar \$100
12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975. 20"W-22"W part: MEC146018

5. Synergy Overall Back Height

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Parts: See XRef, configuration selection with back selection.

- | | | | | | | |
|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|
| ☐ 12" High | ☐ 14" High | ☐ 16" High | ☐ 18" High | ☐ 20" High | ☐ 22" High | ☐ 24" High |
| ☐ 13" High | ☐ 15" High | ☐ 17" High | ☐ 19" High | ☐ 21" High | ☐ 23" High | ☐ 25" High |

6. Synergy Back Cane Height

Omit this selection if a TRU-Comfort back is selected. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles. Parts: See XRef, configuration selection with push handles selection.

- | | | | | | |
|---|---|---|---|---|---|
| ☐ 14" High Canes | ☐ 16" High Canes | ☐ 18" High Canes | ☐ 20" High Canes | ☐ 22" High Canes | ☐ 24" High Canes |
| ☐ 15" High Canes | ☐ 17" High Canes | ☐ 19" High Canes | ☐ 21" High Canes | ☐ 23" High Canes | ☐ 25" High Canes |

16.8. Synergy Armrests

REQUIRED with Synergy seating. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

- | | |
|---|---|
| ☐ ☐ Full Length, Flip Back Armrests Standard
Left part: FRMASMB7038. Right part: FRMASMB7039 | ☐ ☐ Desk Length, Flip Back Armrests Standard
Left part: FRMASMB7040. Right part: FRMASMB7041 |
|---|---|

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

- | | |
|---|---|
| ☐ ☐ Full Length, Two-Post Armrests \$150 ea.
Left Part: FRMASMB7062. Right part: FRMASMB7063. HCPCS: E0973 | ☐ ☐ Desk Length, Two-Post Armrests \$150 ea.
Left part: FRMASMB7060. Right part: FRMASMB7061. HCPCS: E0973 |
|---|---|

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

- | | |
|--|---|
| ☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$75 ea.
Left Part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973 | ☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H \$75 ea.
Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973 |
|--|---|

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

- | | |
|--|--|
| ☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H \$75 ea.
Left Part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973 | ☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H \$75 ea.
Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973 |
|--|--|

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length arm pads.

- | |
|---|
| ☐ ☐ Pediatric, Removable Armrests \$75 ea.
Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973 |
|---|

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

- | | |
|---|---|
| ☐ ☐ Full Length, Cantilever Armrests \$97.50 ea.
Left Part: FRMASMB16188. Right part: FRMASMB16189 | ☐ ☐ Desk Length, Cantilever Armrests \$97.50 ea.
Left part: FRMASMB16190. Right part: FRMASMB16191 |
|---|---|

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

- | | |
|---|---|
| ☐ ☐ Full Length, Single Post Armrests \$75 ea.
Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973 | ☐ ☐ Desk Length, Single Post Armrests \$75 ea.
Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973 |
|---|---|

Synergy seating section continued on following page.

16. SYNERGY SEATING SYSTEM CONTINUED**16.9. Synergy Armpads****REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard
 Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.
 Desk/Full part: UPSASMB1023

Waterfall Armpads..... No Charge
 Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.
 Desk/Full part: UPSASMB1026

No Armpad Selection No Charge
 Part: PTOINDV1373

16.10. Joystick Bracket Receiver**REQUIRED.**

☐ **Standard Fixed Receiver** Standard
 Part: FRMASMB16062

☐ **Stealth Height Adjustable Receiver** \$214
 Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108

**16.11. Synergy Seat to Floor****REQUIRED with Synergy seating.**

15.875"
 Part: PTOINDV3652

17.875"
 Part: PTOINDV3654

19.875"
 Part: PTOINDV3762

16.875"
 Part: PTOINDV3653

18.875"
 Part: PTOINDV3761

20.875"
 Part: PTOINDV3763

17. TRU-BALANCE FLEX SEATING SYSTEM**Optional.** Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D. **Not available with Q-Logic 3 electronics.****17.1. TRU-Balance Flex Frame Selection**

TRU-Balance Flex Seating System No Charge
 Part: PTOINDV3853

17.2. TRU-Balance Flex Seat Dimension - Width**REQUIRED with TRU-Balance Flex seating.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

17.3. TRU-Balance Flex Seat Dimension - Depth**REQUIRED with TRU-Balance Flex seating.** Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------

17.4. TRU-Balance Flex Back Options**REQUIRED with TRU-Balance Flex seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Selection**

Synergy Back Cushion..... \$275
 Part: See XRef

Select a back cane type and height and back height.

TRU-Comfort, 20" Width or less..... \$675
 Part: See XRef. HCPCS: E2620

Select back height.

TRU-Comfort, 22" Width or greater..... \$775
 Part: See XRef. HCPCS: E2621

Select back height.

Sling Back..... \$98
 Part: See XRef

Select a back height. Back cane height must match back height. Not available with flip up cantilever arms.

Omit Back..... No Charge
 Part: PTOINDV1287

Select a back cane type and height in following section.

		Width									
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	15"	S	S	S	-	-	-	-	-	-	-
	16"	S	S	S	-	-	-	-	-	-	-
	17"	-	S	S	S	S	S	S	S	S	S
	18"	-	S	S	S	S	S	S	S	S	S
	19"	-	S	S	S	S	S	S	S	S	S
	20"	-	S	S	S TC	S	S TC	S	S TC	S	S
	21"	-	S	S	S TC	S	S TC	S	S TC	S	S
	22"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	23"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	24"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC
	25"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC
S		Synergy Back									
S TC		Synergy or TRU-Comfort Back									
Synergy and TRU-Comfort Back Availability Chart											

TRU-Balance Flex seating section continued on following pages.

17. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**2. Back Cane Selection**

Synergy Back Cane No Charge

Part: FRMASMB15845

Select a back cane height.

Flip Forward Back Cane \$150

Part: FRMASMB15846

Select a back cane height. Not available with oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane \$195

Part: FRMASMB12067

Select a cane height.

4" Mesomorph (2" per side) Back Cane \$195

Part: FRMASMB12091

Select a cane height.

2" Endomorph (1" per side) Back Cane \$195

Part: FRMASMB12115

Select a cane height.

4" Endomorph (2" per side) Back Cane \$195

Part: FRMASMB12139

Select a cane height.

Manual Recline Back \$750

Part: FRMASMB12910. HCPCS: E1226

Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, or articulating foot platform.

3. TRU-Balance Flex Push Handles

Push Handles Standard

Part: Configuration selection with back cane height

Omit Push Handles \$125

Part: Configuration selection with back cane height

Not available with Sling Back. Minimum cane height is 16".

4. TRU-Balance Flex Accessory Bar

Accessory Bar \$100

12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975.

5. TRU-Balance Flex Overall Back Height

Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Parts: See XRef, configuration selection with back selection.

12" High Back	15" High Back	18" High Back	21" High Back	24" High Back
13" High Back	16" High Back	19" High Back	22" High Back	25" High Back
14" High Back	17" High Back	20" High Back	23" High Back	

6. TRU-Balance Flex Back Cane Height

Omit this selection if a TRU-Comfort back is selected. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles. Parts: See XRef, configuration selection with push handles selection.

14" High Canes	17" High Canes	20" High Canes	23" High Canes
15" High Canes	18" High Canes	21" High Canes	24" High Canes
16" High Canes	19" High Canes	22" High Canes	25" High Canes

17.5. TRU-Balance Flex Armrests**REQUIRED with TRU-Balance Flex seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. TRU-Balance Flex Seat, Drop-In Armrests**
☐ ☐ Full Length, Drop-In Armrests No Charge
☐ ☐ Left part: ARM1803253. Right part: ARM156742

☐ ☐ Desk Length, Drop-In Armrests No Charge
☐ ☐ Left part: ARM1803254. Right part: ARM156743
1b. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests (pair) \$60
☐ ☐ Left part: FRMASMB7038. Right part: FRMASMB7039. HCPCS: E0973

☐ ☐ Desk Length, Flip Back Armrests (pair) \$60
☐ ☐ Left part: FRMASMB7040. Right part: FRMASMB7041. HCPCS: E0973
1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests \$150 ea.
☐ ☐ Part: See XRef. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrests \$150 ea.
☐ ☐ Part: See XRef. HCPCS: E0973
1d. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$75 ea.
☐ ☐ Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests \$75 ea.
☐ ☐ Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973
17.6. TRU-Balance Flex Armpads**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.**1a. Straight Armpads**Straight, Full Length Left and Right No Charge
Part: ARM134275Straight, Full Length Left, Desk Length Right No Charge
Part: ARM134281Straight, Desk Length Left, Full Length Right No Charge
Part: ARM134275Straight, Desk Length Left and Right No Charge
Part: ARM134254**1b. Waterfall Armpads**Waterfall, Full Length Left and Right No Charge
Part: ARM134275Waterfall, Full Length Left, Desk Length Right No Charge
Part: ARM134281Waterfall, Desk Length Left, Full Length Right No Charge
Part: ARM134275Waterfall, Desk Length Left and Right No Charge
Part: ARM134254

TRU-Balance Flex seating section continued on following page.

17. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**17.6. TRU-Balance Flex Armpads continued****1c. Falcon Gel Armpads**

Falcon Gel, Full Length Left and Right \$296.40
 Part: ARM136438

Falcon Gel, Full Length Left, Desk Length Right..... \$274.22
 Part: ARM136442

Falcon Gel, Desk Length Left, Full Length Right..... \$274.22
 Part: ARM136439

Falcon Gel, Desk Length Left and Right \$252.40
 Part: ARM136437

18. LAP BELTS**REQUIRED except with select contoured seats with included lap belt.** Select one lap belt below.**18.1a. Standard Lap Belts**

☐ 38" Lap Belt Part: ACCBELT1006 ☐ 50" Lap Belt Part: ACCBELT1000 ☐ 60" Lap Belt Part: ACCBELT1002 ☐ 70" Lap Belt Part: ACCBELT1003 ☐ 80" Lap Belt Part: ACCBELT1004 ☐ 90" Lap Belt Part: ACCBELT1001

18.1b. Optional Lap Belts

Not available with contoured or TRU-Balance Flex seating.

☐ i-Fit 60" Single Pull, Push Button, Medium \$86
 (7.5") Padded, 1.5" Belt (IFPB5213)
 Part: ST-IFPB5213-Q

☐ i-Fit 60" Single Pull, Push Button, Large \$86
 Padded (9.5"), 2" Belt (IFPB5314)
 Part: ST-IFPB5314-Q

☐ i-Fit 60" Single Pull, Push Button, Extra Large \$86
 Padded (11.5"), 2" Belt (IFPB5315)
 Part: ST-IFPB5315-Q

☐ 60" Sure-Lok Lap Belt (Non-retractable)..... \$120
 Required with WC-19 Transit Option, Synergy Seat only.
 Part: ACCBELT1016

19. HEADRESTS**Optional.** Select a headrest pad and bracket. **Not available with contoured seating.****19.1. Headrest Pads**

☐ Stealth 8" Comfort Plus (CP450)..... \$151.94
 Part: POS148484. HCPCS: E0955


☐ Stealth 10" Comfort Plus (CP250) (A) \$151.94
 Part: SETHEAD1057. HCPCS: E0955

☐ Stealth 14" Comfort Plus (CP150) \$151.94
 Part: INDPART2327. HCPCS: E0955

☐ Pediatric Comfort Plus (6"Wx4.5"H) \$151.94
 (CPN350)
 Part: SETHEAD1069. HCPCS: E0955

☐ Stealth Small Contoured (9.75"Wx3"H) \$64.20
 (STL-OTSM250)
 Part: SETHEAD1052. HCPCS: E0955

☐ Stealth Large Contoured (11"Wx6"H) (B) \$80.25
 (STL-OTLG150)
 Part: SETHEAD1053. HCPCS: E0955

**19.2. Headrest Brackets**

☐ Stealth's Worlds Best Removable.....\$220.50
 Headrest Bracket (TWB485-1) (C)
 Part: FRMASMB15999. HCPCS: E1028

☐ Pediatric Offset (HMO475P-17) (E) \$160.50
 Part: HDWASMB1166. HCPCS: E1028

☐ Stealth Multi-Axis Removable (D).....\$160.50
 Headrest Bracket (HMO475P)
 Part: HDWASMB1084. HCPCS: E1028
 Not available with TRU-Balance Flex seating.

20. SYNERGY POWER OPTION ELECTRONICS**REQUIRED for Synergy Seat with Power Legrests or Elevating Seat.** Select the desired control type for the seating system.

☐ One Actuator Control through NE+ Joystick \$1,895
 Part: See XRef. HCPCS: E2310

☐ Two Actuator Control through NE+ Joystick \$2,375
 Part: See XRef. HCPCS: E2311

☐ One Actuator Control through Q-Logic 3 \$1,895
 Electronics
 Part: See XRef. HCPCS: E2310

☐ Multiple Actuator Control through Q-Logic 3 \$2,895
 Electronics
 Part: See XRef. HCPCS: E2311

21. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**21.1. Center Mount Foot Platform**If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section. **Not available with contoured or TRU-Balance Flex seating.****1. Foot Platform Assembly**

Power Articulating Foot Platform \$2,995
 Part: PTOINDV3754. HCPCS: E1012

Manually Adjustable Center Mount Foot Platform... \$250
 8"-12" extension part: RIG1802138.
 Other extensions: RIG1710345

**2. Lower Extension and Footplates**

8"-12" Extension, 5"x6" Footplates Part: RIG152577 Not available with power articulating platform.	10"-16" Extension, 5"x6" Footplates Part: RIG150206	13"-19" Extension, 5"x6" Footplates Part: RIG151505
8"-12" Extension, 5"x8" Footplates Part: RIG152578 Not available with power articulating platform.	10"-16" Extension, 5"x8" Footplates Part: RIG150203	13"-19" Extension, 5"x8" Footplates Part: RIG151506
	10"-16" Extension, 5"x11" Footplates Part: RIG150212	13"-19" Extension, 5"x11" Footplates Part: RIG151508

Legrests section continued on following page.

21. LEGRESTS CONTINUED

Center Mount Foot Platform Compatibility Matrix	Synergy Seat			☒ No interference issues at this seat-to-floor height. ☐ May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
	Low	Medium	High	
10"-16" Length Extension	!	✓	✓	
13"-19" Length Extension	!	!	!	

3. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H \$150*
Part: RIGASMB7120018

Small Calf Panel 18"Wx8"H \$150
Part: FRMASMB7110012

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$150*
Part: RIGASMB7120019

Large Calf Panel 24"Wx8"H \$200
Part: RIG123824

Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$50
Part: RIG152517

4" Outside Edge (pair) \$50
Part: RIG152510

Width Extension for 11"D Footplates \$50
Part: RIG152518

Only available with 6"D footplates and 8"D footplate extension.

2" Bell-Shaped Length Extension \$50
Part: RIG152519

6" Outside Edge (pair) \$50
Part: RIG152516

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

Available with 8"D and 11"D footplates and 11"D footplate extension.

5. Additional Accessories

Reversed Dartex Calf Pad Cover (pair) \$150
Small part: RIG133717. Large part: RIG133719

Wheel Bumpers (pair) \$12.50 ea.
Part: WHLASMB1578 x2

Reversed Dartex Panel Cover (pair) \$150
Small part: RIG134673. Large part: RIG134744

Footplate Connector \$15
Part: RIGASMB7120023

21.2. 70° Power Articulating, Swing Away Legrests

Not available with contoured or TRU-Balance Flex seating. 4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Articulating Legrest Upper Assembly

☐ ☐ ☐ 70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 14.5"-18" Extension
Pair Part: FRMASMB9073. Left Part: FRMASMB9074.
Right Part: FRMASMB9075. HCPCS: E1010

☐ ☐ ☐ 70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 11"-13" Extension
Pair Part: FRMASMB11444. Left Part: FRMASMB11366.
Right Part: FRMASMB11367. HCPCS: E1010

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	Synergy Seat			☒ This rigging configuration causes no interference issues at this seat-to-floor height. ☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
	Low	Medium	High	
Angle Adjustable, 11" - 13"	✓	✓	✓	
Angle Adjustable, 14.5" - 18"	!	!	!	

2. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 4"Wx6"D
Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566.

☐ ☐ ☐ 5"Wx8"D
Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.

☐ ☐ ☐ 6"Wx11"D¹
Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570.

☐ ☐ ☐ 5"Wx6"D
Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567.

☐ ☐ ☐ 6"Wx8"D¹
Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569.

☐ ☐ ☐ 8"Wx11"D²
Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571.

3. Power Articulating, Swing-Away Legrest Accessories

Calf Pad (pair) 6"Wx6"H No Charge
Part: MEC153735

Falcon Gel Padded Knee Adductor Buttons (KBE)(pair) \$125
Part: INDPART2242. HCPCS: E0953[^]

[^]Please verify the appropriate billing code to be used for this item with all third-party payors.

Legrests section continued on following page.

21. LEGRESTS CONTINUED**21.3. Pediatric, Swing-Away Legrests**

Not available with TRU-Balance Flex seating. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Pediatric Legrests Upper

☐ ☐ ☐ 90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension)\$81 ea.
 P L R Pair part: FRMASMB8808. Left: FRMASMB8809
 Right: FRMASMB8810

2. Angle Adjustable Footplate(s) with Heel Loops

☐ ☐ ☐ 4"Wx6"D Angle Adjustable \$85 ea. ☐ ☐ ☐ 5"Wx6"D Angle Adjustable \$85 ea.
 P L R Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560 P L R Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561
 HCPCS: K0040 HCPCS: K0040

3. Accessories

☐ ☐ ☐ Wheel Bumper(s)\$12.50 ea.
 P L R Pair part: WHL169508. Left: WHL169510. Right: WHL169511

21.4. Swing-Away Legrests

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

☐ ☐ ☐ 60°-90° Swing-Away Legrest \$81 ea. ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short \$81 ea.
 P L R Pair part: RIG120062. Left: FRMASMB17558. P L R Pair part: RIG136481. Left: RIG132607.
 Right: FRMASMB17559 Right: RIG132609.
 Not available with 11"-13" extensions. For use only with 11"-13" lower extensions.

1b. 70° Swing-Away Legrest Uppers

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

☐ ☐ ☐ 70° Swing-Away Legrest (pair)\$162 ☐ ☐ ☐ 70° Tapered, Swing-Away Legrest, Short (pair).....\$162
 Pair part: FRMASMB7843 Pair part: FRMASMB7841
 Not available with 11"-13" extensions. For use only with 11"-13" lower extensions.
☐ ☐ ☐ 70° Tapered, Swing-Away Legrest (pair)\$162
 Pair part: FRMASMB8052
 Not available with 11"-13" extensions.

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat				
	Low	Medium	High		
Extruded, 13" - 16.5"	!	✓	✓	✓	This rigging configuration causes no interference issues at this seat-to-floor height.
Extruded, 15.5" - 19.5"	!	!	!	!	This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Angle Adjustable, 13" - 16.5"	!	✓	✓		
Angle Adjustable, 15.5" - 19.5"	!	!	!		
Multi-Axis, 15.5" - 19.5"	!	!	!	X	This rigging configuration does not work with this seat-to-floor height.

2. Lower Extension

☐ ☐ ☐ 11"-13" Length ☐ ☐ ☐ 11"-16.5" Length ☐ ☐ ☐ 13"-16.5" Length ☐ ☐ ☐ 15.5"-19.5" Length
 P L R Part: See XRef P L R Part: See XRef P L R Part: See XRef P L R Part: See XRef
 N/A with Multi-axis. N/A with Angle Adjustable. N/A with Multi-Axis.

3a. Extruded Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Extruded Footplate Standard ☐ ☐ ☐ 8.5"Wx6"D Extruded Footplate Standard
 P L R Part pair: FRMASMB7847. Left: FRMASMB7861 P L R Part pair: FRMASMB7848. Left: FRMASMB7863
 Right: FRMASMB7862 Right: FRMASMB7864
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 4"Wx6"D \$85 ea. ☐ ☐ ☐ 6"Wx8"D¹ \$85 ea.
 P L R Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 P L R Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
 HCPCS: K0040 HCPCS: K0040
☐ ☐ ☐ 5"Wx6"D \$85 ea. ☐ ☐ ☐ 6"Wx11"D¹ \$85 ea.
 P L R Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 P L R Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
 HCPCS: K0040 HCPCS: K0040
☐ ☐ ☐ 5"Wx8"D \$85 ea. ☐ ☐ ☐ 8"Wx11"D² \$85 ea.
 P L R Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562 P L R Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565
 HCPCS: K0040 HCPCS: K0040

Legrests section continued on following page.

21. LEGRESTS CONTINUED**3c. Multi-Axis Footplate(s)**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐	4"Wx6"D \$247.50 ea. Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552 HCPCS: K0040	☐ ☐ ☐	6"Wx8"D¹ \$247.50 ea. Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555 HCPCS: K0040
☐ ☐ ☐	5"Wx6"D \$247.50 ea. Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553 HCPCS: K0040	☐ ☐ ☐	6"Wx11"D¹ \$247.50 ea. Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558 HCPCS: K0040
☐ ☐ ☐	5"Wx8"D \$247.50 ea. Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554 HCPCS: K0040	☐ ☐ ☐	8"Wx11"D² \$247.50 ea. Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559 HCPCS: K0040

4. Swing-Away Legrest Accessories

Heel Loops (pair) \$50
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable and multi-axis footplates.

Wheel Bumper \$12.50 ea.
Part: See XRef

Gel Neoprene Footrest Sleeve (pair) \$80
Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
QTY Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
QTY Part: INDPART2823. HCPCS: K0108

21.5. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ **Elevating Legrest Assembly** \$162.50 ea.
Pair part: FRMASMB6483. Left: RIG138012
Right: RIG138013. HCPCS: E0990
Select lower extension and footplates

2. Lower Extension

☐ ☐ ☐	11"-14" Length Pair: FRMASMB4755 Left: RIG138026 Right: RIG138030 For reinforced only	☐ ☐ ☐	14"-17" Length Pair: FRMASMB4756 Left: RIG138033 Right: RIG138035 For reinforced only
☐ ☐ ☐	12"-14.5" Length Part: See XRef For angle adjustable only	☐ ☐ ☐	14.5"-18" Length Part: See XRef For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12"-14.5"	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!

 This rigging configuration causes no interference issues at this seat-to-floor height.
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐	7.25"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	☐ ☐ ☐	8.5"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
--	--	--	---

3b. Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐	4"Wx6"D \$85 ea. Pair parts: FRMASMB8812 x2. HCPCS: K0040	☐ ☐ ☐	6"Wx8"D¹ \$85 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040
☐ ☐ ☐	5"Wx6"D \$85 ea. Pair parts: FRMASMB8811 x2. HCPCS: K0040	☐ ☐ ☐	6"Wx11"D¹ \$85 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040
☐ ☐ ☐	5"Wx8"D \$85 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040	☐ ☐ ☐	8"Wx11"D² \$85 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040

4. Elevating Legrest Accessories

Heel Loops (pair) \$50
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable footplates.

Wheel Bumper \$12.50 ea.
Part: See XRef
Available on angle adjustable footplates.

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$125
Part: INDPART2242. HCPCS: E0953^

^Please verify the appropriate billing code to be used for this item with all third-party payors.

22. SET LEGREST LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

23. QUICK STEALTH CUSHION CONFIGURATION

Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the complete [Stealth Cushion Order Form](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

23.1. Cushion Dimensions**Cushion Width**

14" W 16" W 17" W 18" W 19" W 20" W 22" W

Cushion Depth


14" D 16" D 17" D 18" D 19" D 20" D 22" D

23.2. Cushion Model Selection**SIMPLICITY G CUSHION**

General use with 2-way stretch, breathable cover.

Part series: ST-SIMGXXYY


E2601, 16"W to 20" W \$110
E2602, 22"W \$165

**SOLUTION SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SOLSPXXYY

E2607, 14"W to 20" W \$375
E2608, 22"W \$430

**TRU-COMFORT 2 SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover.

Part series: ST-TC2SPPXXYY


E2607, 16"W to 20" W \$400
E2608, 22"W \$500

**SPECTRUM GEL SPP CUSHION**

Positioning with 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCGSPPXXYY

E2607, 14"W to 20" W \$460
E2608, 22"W \$530

**SPECTRUM FOAM SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover.

Part series: ST-SPCFSPXXYY

E2607, 16"W to 20" W \$460
E2608, 22"W \$530

**23.3. Additional Fluid Proof Liners**

Fluid Proof Urethane Liner \$50-\$70

Part: URL1709070-XXYY

Not available with Simplicity and TRU-Comfort 2.

Fluid Proof Polyurethane Liner Coming soon

Not available with Simplicity.

Liner Installation

Install Liner on Cushion

24. PROVIDER INSTALLED ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

Medium, Black, Weather Cover \$168
Part: ACCCOVR1022

Large, Black, Weather Cover \$202
Part: ACCCOVR1021

Oxygen Holder \$118
Contoured seat part: ACCASMB2292. Synergy part: FRMASMB7030.
Flex part: ACC125006. HCPCS: E2208
Not available with Synergy Manual Recline.

Contoured Seat Cup Holder \$20
Part: ACCASMB2287

Synergy/Flex Seat Cup Holder No Charge
Part: ACCASMB2634
Not available on TRU-Balance Flex seating without two-post armrests.

Cell Phone Holder for 12" Lift Captain Seat \$70
Part: ACC1611300

Rear Basket \$118
Contoured seat part: ACCASMB1914
Synergy/Flex part: FRMASMB4821

Walker Holder \$103
Contoured seat part: ACCASMB1912. Synergy/Flex part: FRMASMB7031

Saddle Bag \$35
Part: DWR1020D010

Flag Holder \$25
Contoured seat part: ACCASMB1930
Synergy/Flex part: ACCASMB2415.

Cane/Crutch Holder \$110
Contoured seat part: ACCASMB2294. HCPCS: E2207
Synergy part: FRMASMB4822. Flex part: ACC152515. HCPCS: E2207

Stealth Large Tablet Holder \$375
Part: ACC7120005

Stealth Cell Phone/Media Device Holder \$195
Part: ACC7120006

25. BATTERIES**REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.

- ☐ NF-22, Gel Cell, Interceptor Battery..... \$330 ea.
Part: BAT1704613. HCPCS: E2361
- ☐ NF-22, AGM, MK Battery..... \$360 ea.
Part: BAT121346. HCPCS: E2361
- ☐ No batteries. Provider will assume responsibility of batteries.
Part: Configuration selection
- ☐ NF-22, Gel Cell, MK Battery \$360 ea.
Part: BATGEL1008. HCPCS: E2361

26. BRANCH PATCHES**Optional.**Air Force Patch\$20
Part: ACC140240Coast Guard Patch\$20
Part: ACC140239Navy Patch\$20
Part: ACC140244Army Patch\$20
Part: ACC140241Marines Patch.....\$20
Part: ACC140243Purple Heart Patch\$20
Part: ACC140238**27. MANUFACTURER INSTALLED POSITIONING COMPONENTS****Optional. Not available with contoured seating.** More options are available on the [Positioning Components Order Form](#).**27.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, Rail Mounts, pair (TWBTLTC) \$472.50
Part: INDPART2677. HCPCS: E1028**2. Pads**Curved Pads, 3"x5", pair (TWBL-3X5)..... \$210
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6)..... \$210
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**27.2. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$236.25 ea.
QTY Part: INDPART2804. HCPCS: E1028

Removable Bracket (LP430)..... \$262.50 ea.
QTY Part: INDPART2805. HCPCS: E1028

Removable Bracket with 1" Offset \$400 ea.
QTY (LP430 & LHW-121)
Part: INDPART2854. HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$236.25 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact..... \$236.25 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122).....\$26.25
Part: INDPART2821. HCPCS: K0108
Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129)\$53.50
Part: POS140430. HCPCS: K0108
Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4)..... \$75 ea.
QTY Part: INDPART2801. HCPCS: E0953^

Medium Pad, 4"x8" (TWBADD-4X8)..... \$75 ea.
QTY Part: INDPART2802. HCPCS: E0953^

Long Pad, 4"x12" (TWBADD-4X12)..... \$75 ea.
QTY Part: INDPART2803. HCPCS: E0953^

4. Pelvic/Thigh Guide InstallationInstall on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719

28. PROVIDER INSTALLED POSITIONING COMPONENTSOptional. Not available with contoured seating. More options are available on the [Positioning Components Order Form](#).**28.1. Stealth Foot Boxes**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$498.75
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954^Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$498.75
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954^Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954^**28.2. Stealth Padded Shoe Holder**Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
(SPLE111) (pair)
Part: INDPART2783.Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266
(SPLE114) (pair)
Part: INDPART2785.Stealth Medium Padded Shoe Holder, 4"Wx8.25"D..... \$266
(SPLE112) (pair)
Part: INDPART2784.**28.3. Stealth Elbow Support****1. Elbow Support Bracket**Left, Swing-Away, Elbow Support Bracket \$183.75 ea.
Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028Right, Swing-Away, Elbow Support Bracket \$183.75 ea.
Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028Left, Fixed Elbow Support Bracket..... \$183.75 ea.
TRU-Balance 3 Arms (ARMS-ESFTB3-L)
Part: POS133465Right, Fixed Elbow Support Bracket \$183.75 ea.
TRU-Balance 3 Arms (ARMS-ESFTB3-R)
Part: POS133466**2. Elbow Support Pad**3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$52.50 ea.
Part: POS124065. HCPCS: K01084"Wx4"H Elbow Pad (ARMS-ES-4X4) \$52.50 ea.
Part: POS124067. HCPCS: K0108**28.4. Stealth i-Fit Chest Harness**Upper Classic, Medium (IFPS1103) \$142 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960Upper Contour, Medium (IFPS2103)..... \$142 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960Upper Classic, Large (IFPS1104) \$142 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960Upper Contour, Large (IFPS2104) \$142 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960

Contoured and Classic Chest Harness

28.5. Stealth i-Fit Chest StrapsChest Strap, Medium w/ Belt Anchor \$90 ea.
(IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960Chest Strap, Large w/ Belt Anchor \$90 ea.
(IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960

Chest Strap

28.6. Stealth Swing-Away Headrest Mounts


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

Facial Lateral for Comfort Series..... \$294.25 ea.
Right Side (CPS9-R-HW)
Part: INDPART2647. HCPCS: E1028Facial Laterals for Comfort Series \$347.75 ea.
Left and Right Side (CPS9-B-HW)
Part: INDPART2649. HCPCS: E1028Facial Lateral for Comfort Series..... \$294.25 ea.
Left Side (CPS9-L-HW)
Part: INDPART2648. HCPCS: E1028

Left and Right Facial Lateral with Pads

Facial Lateral Pads and Accessories1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$80.25
Part: POS150532
Required with facial laterals to mount to comfort headrest brackets.Facial Pad, 4"x2" (900-Q)..... \$53.50 ea.
Part: ST-900-Q. HCPCS: K0108Facial Pad, 6"x2" (920-Q)..... \$53.50 ea.
Part: ST-920-Q. HCPCS: K0108Egg Switch Mount (ESM) (A) \$52.50 ea.
Part: INDPART2815Swing-Away Joystick Mount,..... \$250 ea.
18" Rod (C)
Part: ST-SUS18-QFigure 8 Bracket with 13" Rod..... \$267.50 ea.
and Egg Switch Mount (B)
Part: ST-F8675-13-Q

Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

29. SPECIALTY CONTROLS**Optional. Only available with the Synergy Seating system.** Select a type of specialty controls below or to refer to the specialty controls order form.**29.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,720
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard **Match Joystick Shroud to Base Color** \$50
Part: DGN137728 Part: See color below**Non-matching Joystick Shroud**..... \$50
Select color below**Red Corvette**
Part: DGN137729-PRD1**White Satin**
Part: DGN137729-PWH2**Electric Blue**
Part: DGN137729-PBU9**Sugar Plum**
Part: DGN137729-PPR7**Cherry Bomb**
Part: DGN137729-PRD4**Deep Purple**
Part: DGN137729-PPR5**Blue Jean**
Part: DGN137729-PBU11**Military Digital Camo**
Part: DGN137729-WCM7**Back in Black**
Part: DGN137729-PBK1**Raspberry Beret**
Part: DGN137729-PPK4**Lemon Crush**
Part: DGN137729-PYL4**Faded Glory**
Part: DGN137729-WFL2**Night Train**
Part: DGN137729-PBK3**Evolution Orange**
Part: DGN137729-POR7**Green Machine**
Part: DGN137729-PGN11**Carbon Fiber**
Part: DGN137729-WCB1**Grey Street**
Part: DGN137729-PSL3**Midnight Blue**
Part: DGN137729-PBU10**Military Green**
Part: DGN137729-PGN8**To Be Determined**
Follow-up required.

Enhanced Display

29.2. Q-Logic 3 SAJ☐ **Q-Logic 3 EX SAJ**.....\$700
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

29.3. Head Array**Adult, i-Drive Tri-Array Series, 3-Switch** \$4,414
Head Array (IDHBT200-1)
Part: ST-IDHBT200-1-Q, HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461
(IDHBT300-1)
Part: ST-IDHBT300-1-Q, HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

29.4. Sip and Puff Drive**Module****Sip and Puff Module for Q-Logic 3**.....\$1,895
Part: CTL168056

Includes module and mounting hardware.



Sip & Puff Module

Straw and Mount**Stealth Sip and Puff for Comfort Series, Left Mount (SP625-CL)**\$294.25
Part: ST-SP625-CL-Q

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U)\$294.25
Part: ST-SP625-U-Q

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, Right Mount (SP622-CR)\$294.25
Part: ST-SP625-CR-Q

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip & Puff

Specialty electronics section continued on following page.

Account Number: _____ Date: _____ Marked for: _____

29. SPECIALTY CONTROLS CONTINUED

29.5. Egg Switches

Red Egg Switch (ESRED) \$80.85 ea.
 Part: SWTMCRO1077
 QTY _____
 Blue Egg Switch (ESBLU) \$80.85 ea.
 Part: SWTMCRO1079
 QTY _____
 Black Egg Switch (ESBLK) \$80.85 ea.
 Part: SWTMCRO1081
 QTY _____

Green Egg Switch (ESGRN) \$80.85 ea.
 Part: SWTMCRO1078
 QTY _____
 Yellow Egg Switch (ESYEL) \$80.85 ea.
 Part: SWTMCRO1080
 QTY _____



Egg Switch

29.6. Switch Adapter

2 to 1 Switch Adapter \$13.39
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

30. TRU-BALANCE 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA) No Charge
 Part: Configuration selection
 Secondary Multiplier (SM) \$221.21
 Part: HAR144004
 Special Order Configuration (SOF)
 Part: Configuration selection

Arm Mounts	Back Mounts			
		1 Arm	2 Arm	3+ Arm
Joystick (Hand Control) Enhanced Display iAccess USB Charger SAJ	Attendant Control SCIM Sip n Puff ECU	0 Back	NA	SM
		1 Back	NA	SM
		2 Back	SM	SM
		3+ Back	SOF	SOF

Q-Logic 3 Multiplier with Synergy Seat Chart

31. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2018 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.

