

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com**Q6 Edge® 3 Order Form with TRU-Balance® 3 Positioning**

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

**INTRODUCTION**

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


Thank you for choosing Quantum!**PATIENT INFORMATION**

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

**1. BASE MODEL**

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3S-SS \$7,895

Part: Q6 EDGE 3 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Quantum Q6 Edge 3 3SP-SS \$7,895

Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856

Single power option.

Quantum Q6 Edge 3 3MP-SS \$7,895

Part: Q6 EDGE 3 3MP-SS. HCPCS: K0861

Multiple actuator power option.

2. PATIENT'S WEIGHT**REQUIRED.** The patient's weight is required to quote and order the power chair.

Patient's Weight: _____

3. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**3.1a. Standard Shrouds**

 Red Corvette Part: DGN1801101-PRD1	 White Satin Part: DGN1801101-PWH2	 Electric Blue Part: DGN1801101-PBU9	 Military Green Part: DGN1801101-PGN8
 Cherry Bomb Part: DGN1801101-PRD4	 Deep Purple Part: DGN1801101-PPR5	 Blue Jean Part: DGN1801101-PBU11	 Sugar Plum Part: DGN1801101-PPR7
 Back in Black Part: DGN1801101-PBK1	 Raspberry Beret Part: DGN1801101-PPK4	 Lemon Crush Part: DGN1801101-PYL4	To Be Determined Follow-up required.
 Night Train Part: DGN1801101-PBK3	 Evolution Orange Part: DGN1801101-POR7	 Green Machine Part: DGN1801101-PGN11	
 Grey Street Part: DGN1801101-PSL3	 Midnight Blue Part: DGN1801101-PBU10		

3.1b. Hydrographic Shrouds

☐  Military Digital Camo \$250 Part: DGN153201-WCM7	☐  Faded Glory \$250 Part: DGN153201-WFL2	☐  Carbon Fiber \$250 Part: DGN153201-WCB1
---	---	--

4. DRIVE MOTORS**REQUIRED.**

☐ 6 mph Drive Motors Standard Part: MOT147433	☐ Accu-Trac Motor Technology \$1,295 Part: MOT147881. HCPCS: K0108 Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.
--	--

5. BASE MOUNTED FOOT PLATFORM**REQUIRED.** May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard Part: FRMASMB16077	Omit base mounted foot platform No Charge Part: PTOINDV3674
Bell Shaped, Base Mounted Foot Platform \$235 Part: RIGASMB7120041	

6. DRIVE WHEELS**REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".**6.1. Drive Wheels**

3"x14" Flat Free Black Tires Standard Part: PTOINDV3525 (pair)	3"x14" Pneumatic Black Tires No Charge Part: WHL1713925 x2
4"x14" Knobby, Flat Free Black Tires \$275.50 ea. Part: WHL147601 x2	4"x14" Knobby, Pneumatic Black Tires \$275.50 ea. Part: WHL147599 x2

6.2. Drive Wheel Accent Color

Not available with knobby tires.

Back in BlackStandard Part: DGN164245-PBK1 x2	Raspberry BeretNo Charge Part: DGN164245-PPK4 x2	Military GreenNo Charge Part: DGN164245-PGN8 x2
Red CorvetteNo Charge Part: DGN164245-PRD1 x2	Evolution OrangeNo Charge Part: DGN164245-POR7 x2	Sugar PlumNo Charge Part: DGN164245-PPR7 x2
Cherry BombNo Charge Part: DGN164245-PRD4 x2	Midnight BlueNo Charge Part: DGN164245-PBU10 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
Night TrainNo Charge Part: DGN164245-PBK3 x2	Electric BlueNo Charge Part: DGN164245-PBU9 x2	Faded Glory \$50 ea. Part: DGN164245-WFL2 x2
Grey StreetNo Charge Part: DGN164245-PSL3 x2	Blue JeanNo Charge Part: DGN164245-PBU11 x2	Carbon Fiber \$50 ea. Part: DGN164245-WCB1 x2
White SatinNo Charge Part: DGN164245-PWH2 x2	Lemon CrushNo Charge Part: DGN164245-PYL4 x2	To Be Determined Follow-up required.
Deep PurpleNo Charge Part: DGN164245-PPR5 x2	Green MachineNo Charge Part: DGN164245-PGN11 x2	

7. FENDER LIGHTS**Optional.** Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.
☐ Switch Fender Lights No Charge
 Part: DGN154117
Fender Light Switch Mount Type
 Side Rail/Track Mount Bracket Standard
 Part: MEC1800866

Allows the egg switch to be mounted to seating side rails using a t-nut assembly. This may also be mounted to other armrest, seat or accessory bar rails. Not available with contoured seating.

 Under Arm Mounting Bracket \$25
 Part: MEC1802726

Allows the egg switch to mount to a joystick tube-style bracket mounted to any armrest. Includes joystick mounting bracket.

8. FULL LIGHT PACKAGES**Optional.** Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package. Select one available package below.
☐ Base Mounted, LED, Full Light Package \$995
 Part without iLevel: ELE1705999. Part with iLevel: ELE168643

Fender lights will be omitted with Base Mounted, LED, Full Light Package.

☐ Seat Mounted, LED, Full Light Package \$995
 Part: ELE131477


Seat Mounted LED Lights (One Shown)

9. ELECTRONICS**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**9.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

 NE 4-Key Joystick Standard
 Part: CTLDC1554

Available with Tilt Only through Toggle.

 NE+ 6-Key Joystick (1 Actuator) No Charge
 Part: CTLDC1560

Control of one actuator through Joystick on single power base model (3SP-SS) only.

 NE+ 6-Key Joystick (2 Actuators) \$200
 Part: CTLDC1560

Control up to 2 actuators through joystick.

 NE+ 8-Key Joystick (2 Actuators) \$400
 Part: CTLDC1561

Control up to 2 actuators through joystick and full light package.

 No Joystick, Setup for NE+ Joystick No Charge
 Part: Configuration selection


NE+ 6 Key Joystick

9.1b. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

 Q-Logic 3 EX Controller Module \$695
 Part: See XRef. HCPCS: E2377

 Q-Logic 3 Expandable Controller Harness \$425
 Part: HAR144005. HCPCS: E2313
2. Q-Logic 3 Joystick
 Q-Logic 3 EX Joystick \$700
 Part: CTL167724

 No Joystick, Setup for Q-Logic
 Part: Configuration selection.

 Q-Logic 3 EX Joystick for Full Lights \$700
 Part: CTL167725


Q-Logic 3 EX Joystick

9.2. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 Joystick selections above.** Reference colors shown on page 2.
 Black, Soft Touch Standard
 Part: DGN124307

 Non-matching Armor \$50
 Part: PTOINDV3422

Select color below.

 Match Armor to Base Color \$50
 Part: PTOINDV3421

 Red Corvette
 Part: DGN124336-PRD1

 White Satin
 Part: DGN124336-PWH2

 Electric Blue
 Part: DGN124336-PBU9

 Sugar Plum
 Part: DGN124336-PPR7

 Cherry Bomb
 Part: DGN124336-PRD4

 Deep Purple
 Part: DGN124336-PPR5

 Blue Jean
 Part: DGN124336-PBU11

 Military Digital Camo
 Part: DGN124336-WCM7

 Back in Black
 Part: DGN124336-PBK1

 Raspberry Beret
 Part: DGN124336-PPK4

 Lemon Crush
 Part: DGN124336-PYL4

 Faded Glory
 Part: DGN124336-WFL2

 Night Train
 Part: DGN124336-PBK3

 Evolution Orange
 Part: DGN124336-POR7

 Green Machine
 Part: DGN124336-PGN11

 Carbon Fiber
 Part: DGN124336-WCB1

 Grey Street
 Part: DGN124336-PSL3

 Midnight Blue
 Part: DGN124336-PBU10

 Military Green
 Part: DGN124336-PGN8

 To Be Determined
 Follow-up required.
9.3. SCIM Module**Optional.**
☐ SCIM Module for Q-Logic 3 \$525
 Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

10. JOYSTICK HANDLES

Optional. The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- | | |
|---|--|
| ☐ Mushroom Handle (B) \$85
Part: ACCASMB1860. HCPCS: E2323 | ☐ Extended I-Handle (G) \$85
Part: ACCASMB1857. HCPCS: E2323 |
| ☐ Chin Cup (C) \$85
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | ☐ Bodypoint 3" Goal Post Handle \$99
Part: PLSKNCP1057. HCPCS: E2323 |
| ☐ T-Bar (D) \$120
Part: ACCASMB1858. HCPCS: E2323 | ☐ Bodypoint 4" Goal Post Handle \$99
Part: PLSKNCP1058. HCPCS: E2323 |
| ☐ Medium Ball Handle (E) \$65
Part: ACCASMB1859. HCPCS: E2323 | ☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| ☐ Soft Ball Handle (F) \$65
Part: ACCASMB1856. HCPCS: E2323 | |


**11. JOYSTICK MOUNTING BRACKETS**

Optional. Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket


11.1a. Fixed, In-line Joystick Bracket

- | | |
|--|--|
| ☐ Fixed, In-line, Bracket Standard
Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744 | |
| ☐ Left Side Mount
Part: PTOINDV1225 | ☐ Right Side Mount
Part: PTOINDV1226 |

**11.1b. Swing-Away (Retractable) Joystick Bracket**


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

- | | |
|---|--|
| Swing-away, Inline, Left Side Mount \$245
NE/NE+ Part: ACC142526. Q-Logic: ACC142524.
HCPCS: E1028 | Swing-away, 1" Offset, Left Side Mount \$245
NE/NE+ Part: ACC142747. Q-Logic: ACC142742.
HCPCS: E1028 |
| Swing-away, Inline, Right Side Mount \$245
NE/NE+ Part: ACC142527. Q-Logic: ACC142525.
HCPCS: E1028 | Swing-away, 1" Offset, Right Side Mount \$245
NE/NE+ Part: ACC142748. Q-Logic: ACC142743.
HCPCS: E1028 |

**11.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

- | | |
|--|---|
| Retract4 Mount \$350
NE/NE+ Part: MEC138876. HCPCS: E1028
Q-Logic Part: MEC138860. HCPCS: E1028
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted. | Stealth Swing-Away/Flip Down Bracket \$294.25
NE/NE+ Part: ST-TWBM480QL-Q. HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028
SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted. |
| Multi-Axis Mount \$100
NE/NE+ Part: ACC142528. HCPCS: E1028
Q-Logic Part: ACC131216. HCPCS: E1028
SAJ Part: ACC142740. HCPCS: E1028
Allows height adjustment of joystick. | |

**Specialty Mount Side**

- | | |
|---|--|
| ☐ Left Side Mount
Part: PTOINDV1225 | ☐ Right Side Mount
Part: PTOINDV1226 |
|---|--|

12. ELECTRONICS ACCESSORIES

Optional.

- | | |
|---|---|
| USB Mobile Device Charger \$175
Part: ELE168069
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier. | Q-Logic EX/NE Programmer \$545
Part: ELEASMB6976 |
| | Q-Logic PC Programming Disk with Cable \$750
Part: ELEASMB5216 |



13. TRU-BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™**REQUIRED.** Maximum seat size is 22"W x 22"D.**13.1. TRU-Balance 3 Configuration**

TRU-Balance 3 Static Seat\$2,295
 Part: PTOINDV3442
 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Power Tilt Only\$6,995
 Part: PTOINDV3443. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Power Recline Only\$7,595
 Part: PTOINDV3444. HCPCS: E1004
 175° of power recline.

TRU-Balance 3 Power Tilt & Recline.....\$11,895
 Part: PTOINDV3446. HCPCS: E1007
 50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

13.2. 12" Power Lift or Manual Recline**Optional.** Select one option below.

TRU-Balance 3, 12" Power Lift with iLevel®\$2,600
 Part: PTOINDV3773. HCPCS: E2300
 12" of scissor lift with iLevel® System. Not available with high speed motors or manual recline.



Attendant Manual Recline*\$500
 Part: PTOINDV3501. HCPCS: E1226

*Available through Special Order when ordered with fender lights. Not available with NE+ joystick with Tilt through toggle. Available with Static Seat or Power Tilt. Not available with Power Lift.

13.3. Seat Dimension - Width

REQUIRED with seat selection. Maximum seat width is 22". Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.4. Seat Dimension - Width Package**REQUIRED with seat selection.** Width and depth package part: See XRef: TB3SPACK

12"-16" Width Package	14"-16" Width Package	16"-20" Width Package	20"-24" Width Package (Maximum Width is 22")
-----------------------	-----------------------	-----------------------	---

13.5. Seat Dimension - Depth**REQUIRED with seat selection.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.6. Seat Dimension - Depth Package**REQUIRED with seat selection.** Width and depth package part: See XRef: TB3SPACK

12"-16" Depth Package	14"-16" Depth Package	16"-20" Depth Package	20"-24" Depth Package (Maximum Depth is 22")
-----------------------	-----------------------	-----------------------	---

13.7. Occupied Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option.....\$275
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

13.8. Seat to Floor Height**REQUIRED.** Select a seat to floor height compatible with the seat functions previously selected.

16.5" Seat to Floor Part: PTOINDV3459 Not available with Power Tilt or Power Lift Configurations.	17.5" Seat to Floor Part: PTOINDV3460	18.5" Seat to Floor Part: PTOINDV3463	19.5" Seat to Floor Part: PTOINDV3505 Not available with Static Seat or Power Recline only.
---	--	--	---

13.9. Back Type Selection**REQUIRED.** Match choice in this section with the options in the following back cushion configuration section.

TRU-Balance 3 Back Assembly (Static Back)
 Part: Configuration selection
 Select a back cushion configuration in the subsequent section.

Back Cane Option
 Part: Configuration selection
 Available only with TRU-Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Flip Forward Back.....\$150
 Part without Lift: MEC131957, Part with Lift: MEC148117
 Available with Static Seat and Power Tilt (1 actuator). Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width or with attendant manual recline. Select a back cushion type in the subsequent section.

13.10. Back Cane Selection**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**

Angle Adjustable No Charge Part: MEC140266	1" Endo.....\$70 Part: MEC141455
Flip Forward\$25 Part: MEC140272	2" Endo\$70 Part: MEC141456

2. Back Cane Height

16"H Cane Part: FRMASMB15834	20"H Cane Part: FRMASMB15838	24"H Cane Part: FRMASMB15842
17"H Cane Part: FRMASMB15835	21"H Cane Part: FRMASMB15839	25"H Cane Part: FRMASMB15843

TRU-Balance 3 Seating section continued on following pages.

13. TRU-BALANCE 3 SEATING CONTINUED**13.11. Back Cushion Configuration - Type**

REQUIRED. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

Sport Back No Charge
Part: See XRef

Complete following width package and size.

ADI Q-Back, Active Contour \$525

Part: See XRef. HCPCS: E2613

Complete following width package and size.

ADI Q-Back, Deep Contour \$595

Part: See XRef. HCPCS: E2615

Complete following width package and size.

ADI Q-Back, Extra Deep Contour \$650

Part: See XRef. HCPCS: E2620

Complete following width package and size.

20" Wide Extra Deep Contour is not compatible with Gel and Waterfall Armpads, Endo/Meso Downposts, 2 Post Arms, or Ottobock arm blocks unless wider arm bars are installed on the chair.

TRU-Comfort 2 Back, 20" Width or less \$675

Part: See XRef. HCPCS: E2620

Complete following width package and size.

TRU-Comfort 2 Back, 22" Width or greater \$775

Part: See XRef. HCPCS: E2621

Complete following width package and size.

Omit Back Pan/Cushion No Charge

Part: PTOINDV1287

Select a following width package.

After-market Back Interface \$200

Part 16"-20"W: SET131920. Part 20"-24"W: SET131921

Select a set width in following section, range of 16"-24". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs. A Special Order is required on back widths less than 16".

13.12. Back Cushion Configuration - Back Width Package

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections

Part: PTOINDV3468

10"W-14"W Package

Part: Configuration selection, back cushions

16"W-20"W Package

Part: Configuration selection, back cushions

20"W-24"W Package

Part: Configuration selection, back cushions

13.13. Back Cushion Configuration - Back Width

REQUIRED except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

12" W (#3562)	14" W (#3563)	16" W (#3564)	17" W (#3565)	18" W (#3566)	19" W (#3567)	20" W (#3568)	21" W (#3569)	22" W (#3570)	24" W (#3571)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.14. Back Cushion Configuration - Overall Back Height

Static Seat and Power Tilt available height range: 14"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height Part: SETHGHT1027	22" Height Part: SETHGHT1035
15" Height Part: SETHGHT1028	23" Height Part: SETHGHT1036
16" Height Part: SETHGHT1029	24" Height Part: SETHGHT1037
17" Height Part: SETHGHT1030	25" Height Part: SETHGHT1038
18" Height Part: SETHGHT1031	26" Height Part: SETHGHT1039
19" Height Part: SETHGHT1032	27" Height Part: SETHGHT1040
20" Height Part: SETHGHT1033	28" Height Part: SETHGHT1041
21" Height Part: SETHGHT1034	

	Tilt or Static Width												Recline Width											
	12"	14"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	16"	17"	18"	19"	20"	21"	22"	24"				
Overall Height	14"	ST	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
	15"	ST	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
	16"	S	ST	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-				
	17"	S	ST	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-				
	18"	-	-	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-				
	19"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	S	ST	ST	ST	ST	ST	S	ST	S				
	20"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	S	ST	ST	ST	ST	ST	S	ST	S				
	21"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	-	-	ST	ST	ST	ST	S	ST	S				
	22"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	23"	-	-	ST	ST	S	ST	ST	S	ST	ST	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	24"	-	-	S	S	ST	S	ST	S	S	S	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	25"	-	-	S	S	ST	S	ST	S	S	S	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST			
	26"	-	-	-	-	-	-	-	-	-	-	-	-	ST	ST	ST	ST	ST	S	ST	ST			
	27"	-	-	-	-	-	-	-	-	-	-	-	-	S	S	ST	S	ST	S	S	S			
	28"	-	-	-	-	-	-	-	-	-	-	-	-	S	S	ST	S	ST	S	S	S			

ST

 Sport Back

ST

 Sport or TRU-Comfort 2 Back

SQT

 Sport, Q-Back, or TRU-Comfort 2 Back

Sport, Q-Back, and TRU-Comfort 2 Back Availability Chart

S Sport Back **T** Sport or TRU-Comfort 2 Back
ST Sport, Q-Back, or TRU-Comfort 2 Back

Sport, Q-Back, and TRU-Comfort 2 Back Availability Chart

13.15. TRU-Balance 3 Back Shroud

REQUIRED. The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette
Part: DGN148034

Cherry Bomb
Part: DGN1801259-PRD4

Back in Black
Part: DGN148032

Night Train
Part: DGN1706558

Grey Street
Part: DGN148074

White Satin
Part: DGN148075

Deep Purple
Part: DGN148046

Raspberry Beret
Part: DGN1801259-PPK4

Evolution Orange
Part: DGN1801259-POR7

Midnight Blue
Part: DGN148035

Electric Blue
Part: DGN148076

Blue Jean
Part: DGN1801259-PBU11

Lemon Crush
Part: DGN1801259-PYL4

Green Machine
Part: DGN1801259-PGN11

Military Green
Part: DGN148037

Sugar Plum
Part: DGN1801259-PPR7

Military Digital Camo
Part: DGN1707171

Faded Glory
Part: DGN1707173

Carbon Fiber
Part: DGN1707169

To Be Determined
Follow-up required.

TRU-Balance 3 Seating section continued on following page.

13. TRU-BALANCE 3 SEATING CONTINUED**13.16. TRU-Balance 3 Power Base Interface Shroud**

REQUIRED. The TRU-Balance 3 power base interface shrouds are located on the left and right interface plates of the TRU-Balance 3 system, which mount to the power base, and are available in most matching power base shroud colors. Reference colors shown on page 2.

Red Corvette

Lift part: DGN1802631-PRD1 x2
Non-lift part: DGN1801547-PRD1

Cherry Bomb

Lift part: DGN1802631-PRD4 x2
Non-lift part: DGN1801547-PRD4

Back in Black

Lift part: DGN1802631-PBK1 x2
Non-lift part: DGN1801547-PBK1

Night Train

Lift part: DGN1802631-PBK3 x2
Non-lift part: DGN1801547-PBK3

Grey Street

Lift part: DGN1802631-PSL3 x2
Non-lift part: DGN1801547-PSL3

White Satin

Lift part: DGN1802631-PWH2 x2
Non-lift part: DGN1801547-PWH2

Deep Purple

Lift part: DGN1802631-PPR5 x2
Non-lift part: DGN1801547-PPR5

Raspberry Beret

Lift part: DGN1802631-PPK4 x2
Non-lift part: DGN1801547-PPK4

Evolution Orange

Lift part: DGN1802631-POR7 x2
Non-lift part: DGN1801547-POR7

Midnight Blue

Lift part: DGN1802631-PBU10 x2
Non-lift part: DGN1801547-PBU10

Electric Blue

Lift part: DGN1802631-PBU9 x2
Non-lift part: DGN1801547-PBU9

Blue Jean

Lift part: DGN1802631-PBU11 x2
Non-lift part: DGN1801547-PBU11

Lemon Crush

Lift part: DGN1802631-PYL4 x2
Non-lift part: DGN1801547-PYL4

Green Machine

Lift part: DGN1802631-PGN11 x2
Non-lift part: DGN1801547-PGN11

Military Green

Lift part: DGN1802631-PGN8 x2
Non-lift part: DGN1801547-PGN8

Sugar Plum

Lift part: DGN1802631-PPR7 x2
Non-lift part: DGN1801547-PPR7

To Be Determined

Follow-up required.

13.17. Additional TRU-Balance 3 Back Options**Optional.**

TRU-Comfort Lateral Positioning Wedge, pair..... \$30
Part: SET123465

Available with TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

See Comfort Company Back order form

Privacy Flap, 9.75" \$35
Part: DWR1044E456B

Standard on recline seating systems.

Privacy Flap, 12.75" \$35

Part: DWR1044E456A

Standard on recline seating systems.

14. TRU-BALANCE 3 ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

14.1a. TRU-Balance 3 Flip-Back Armrests

Choose a left and right TRU-Balance 3 armrest and set height in the following section. An optional downpost is also available. Different left (L) and right (R) armrest height packages are allowed. Not available with Back Cane Option.

☐ L ☐ R TRU-Balance 3 Flip-Back, 8"-12"H Standard
Part: See XRef

☐ L ☐ R TRU-Balance 3 Flip-Back, 12"-16"H Standard
Part: See XRef

Flip-Back Armrest Support

TRU-Balance 3 Armrest Downpost (Pair)..... \$87.50 ea.
Part: See XRef

Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available on 14" Width seat or back.

14.1b. Two-Post Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

2-Post, Flip Back, Height Adjustable, and Removable Armrests (9"-13" Height) \$150 ea. or No Charge*
Part: See XRef. HCPCS: E0973

Not available with a depth less than 16". Required with flip forward back and lap trays. *No charge with power seating.

14.1c. Heavy Duty, Removable Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$75 ea.
Heavy Duty (Sold in pairs)

Part: See XRef. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable, \$75 ea.
Heavy Duty (Sold in pairs)

Part: See XRef. HCPCS w/ Static Seat: E0973

14.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H

PTOINDV3535

10" H

PTOINDV3537

12" H

PTOINDV3539

14" H

PTOINDV3541

16" H

PTOINDV3543

9" H

PTOINDV3536

11" H

PTOINDV3538

13" H

PTOINDV3540

15" H

PTOINDV3542

2. Right Armrest Height

8" H

PTOINDV3450

10" H

PTOINDV3452

12" H

PTOINDV3454

14" H

PTOINDV3456

16" H

PTOINDV3458

9" H

PTOINDV3451

11" H

PTOINDV3453

13" H

PTOINDV3455

15" H

PTOINDV3457

14.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ L ☐ R 14" Full Length
Part: Configuration selection

☐ L ☐ R 10" Desk Length
Part: Configuration selection

Account Number: _____ Date: _____ Marked for: _____

15. TRU-BALANCE 3 ARMPADS

REQUIRED. Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

No Armpad Selection and/or See Positioning Components Page or Order Form No Charge
Part: Configuration selection

15.1a. Foam Armpads

☐ L ☐ R Straight, Armpad Standard
Part: See XRef

☐ L ☐ R Waterfall, Armpad No Charge
Part: See XRef

15.1b. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for TRU-Balance 3 Armrests

Falcon Gel, Full Length Left and Right \$185
Part: ARM136417

Falcon Gel, Desk Length Left, Full Length Right \$162.50
Part: ARM136418

Falcon Gel, Full Length Left, Desk Length Right \$162.50
Part: ARM136426

Falcon Gel, Desk Length Left and Right \$140
Part: ARM136416

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right \$296.40
Part: ARM136438

Falcon Gel, Desk Length Left, Full Length Right \$274.22
Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right \$274.22
Part: ARM136442

Falcon Gel, Desk Length Left and Right \$252.40
Part: ARM136437

16. TRU-BALANCE 3 JOYSTICK BRACKET RECEIVER

REQUIRED.

Standard Fixed Receiver Standard
Part: ACC123525

Stealth Height Adjustable Receiver \$214
Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108
TRU-Balance 3 Arm Part: ST-ARM250T-Q. HCPCS: K0108



Height Adjustable
Receiver (Tube Arm)



Height Adjustable
Receiver (Track Arm)

17. TRU-BALANCE 3 HEADRESTS

Optional. Select a headrest pad and bracket.

17.1. Headrest Pads

☐ Stealth 8" Comfort Plus (CP450) \$151.94
Part: POS148484. HCPCS: E0955

☐ Stealth Small Contoured (9.75"Wx3"H) \$64.20
(STL-OTS250)
Part: SETHEAD1052. HCPCS: E0955

☐ Stealth 10" Comfort Plus (CP250)(A) \$151.94
Part: SETHEAD1057. HCPCS: E0955

☐ Stealth Large Contoured (11"Wx6"H) (B) \$80.25
(STL-OTLG150)
Part: SETHEAD1053. HCPCS: E0955


☐ Stealth 14" Comfort Plus (CP150) \$151.94
Part: INDPART2327. HCPCS: E0955

☐ Pediatric Comfort Plus (6"Wx4.5"H) \$151.94
(CPN350)
Part: SETHEAD1069. HCPCS: E0955

17.2. Headrest Brackets

☐ Stealth's Worlds Best Removable \$220.50
Headrest Bracket (TWB485-1) (C)
Part: FRMASMB15999. HCPCS: E1028

☐ Stealth Multi-Axis Removable Headrest \$160.50
Bracket (HMO475P) (D)
Part: HDWASMB1084. HCPCS: E1028



Headrests

18. LAP BELTS

REQUIRED. Select a lap belt below.

18.1a. Standard Lap Belts

☐ 38" Lap Belt
Part: ACCBELT1006

☐ 50" Lap Belt
Part: ACCBELT1000

☐ 60" Lap Belt
Part: ACCBELT1002

☐ 70" Lap Belt
Part: ACCBELT1003

☐ 80" Lap Belt
Part: ACCBELT1004

☐ 90" Lap Belt
Part: ACCBELT1001

18.1b. Optional Lap Belts

☐ i-Fit 60" Single Pull, Push Button, Medium \$86
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

☐ i-Fit 60" Single Pull, Push Button, Extra Large \$86
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

☐ i-Fit 60" Single Pull, Push Button, Large \$86
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

60" Sure-Lok Lap Belt (Non-retractable) \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016

19. TRU-BALANCE 3 MANUFACTURER INSTALLED ACCESSORIES**Optional.**

Trilogy Vent Tray\$400
 Part: ACC136791. HCPCS: E1029

Push Handles (pair) \$175
 Part: ACC147772

Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.

Compact Vent Tray\$1,395
 Part: ACC133889. HCPCS: E1029

Accessory Bar \$100
 Part: MEC125098

Standard with the Trilogy vent tray, and oxygen holder. May be added to mount additional accessory to back of unit. Not available with back cane option.

**20. TRU-BALANCE 3 LEGREST CONTROLS**

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

☐ Independent Legrests (adds 2 actuator control requirement)
 Part: Configuration selection

☐ Articulating Foot Platform (adds 1 actuator control requirement)
 Part: Configuration selection

☐ Combined Legrests (adds 1 actuator control requirement)
 Part: Configuration selection

21. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS

REQUIRED with power seating or legrest options.

21.1a. Tilt Only through Toggle

Single Actuator Control, Through Toggle (Power Tilt Only) No Charge
 Part: See XRef

Only available with Power Tilt only. Select toggle side. Will require extended lead time when selected with Fender Lights.

Left Side Toggle
 Part: PTOINDV1838

Right Side Toggle
 Part: PTOINDV1839

21.1b. NE+ 6 Key Controls through Joystick

Single Actuator Control Through NE+\$1,895
 Joystick
 Part: See XRef. HCPCS: E2310

Two Actuator Control Through NE+ \$2,895
 Joystick
 Part: See XRef. HCPCS: E2311

21.1c. Q-Logic 3 Controls through Joystick

Single Actuator Through Q-Logic 3\$1,895
 Controller and Interface (AM1)
 Part: See XRef. HCPCS: E2310

Multiple Actuator Through Q-Logic 3..... \$2,895
 Controller and Interface (AAM)
 Part: See XRef. HCPCS: E2311

Single Actuator Through Q-Logic 3 \$2,895
 Controller and Interface (AAM)
 Part: See XRef. HCPCS: E2310

Optional Switch Control

Not available with AM1 electronics. Switch will be mounted to same side of joystick.

4 Way Toggle
 Part: ELEASMB5925

4 Way Button Toggle
 Part: ELEASMB5926

**21.1d. Q-Logic 3 Controls with iAccess Switch Control**

iAccess control is not available with combined or independent power legrests.

Single Actuator Control, Through Flat Button\$2,295
 iAccess and Q-Logic 3 Controller and Interface (AM1)
COMING SOON. Part: See XRef. HCPCS: E2310

Single Actuator Control, Through Toggle Button\$2,295
 iAccess and Q-Logic 3 Controller and Interface (AM1)
COMING SOON. Part: See XRef. HCPCS: E2310

Single Actuator Control, Through Flat Button\$3,295
 iAccess and Q-Logic 3 Controller and Interface (AAM)
 Part: See XRef. HCPCS: E2310

Single Actuator Control, Through Toggle Button\$3,295
 iAccess and Q-Logic 3 Controller and Interface (AAM)
 Part: See XRef. HCPCS: E2310

Multiple Actuator Control Through Flat Button\$3,295
 iAccess and Q-Logic 3 Controller and Interface (AAM)
 Part: See XRef. HCPCS: E2311

Multiple Actuator Control Through Toggle Button.....\$3,295
 iAccess and Q-Logic 3 Controller and Interface (AAM)
 Part: See XRef. HCPCS: E2311

iAccess Mount Location

Left Side iAccess
 Part: PTOINDV3678

Right Side iAccess
 Part: PTOINDV3679

21.1e. Q-Logic 3 Controls with 4 Function Push Button

Single Actuator Control Through 4 Function\$2,895
 Push Button and Q-Logic 3 Controller and Interface (AAM)
 Only available with TB3 armrests. Select a side to mount.
COMING SOON. Part: See XRef. HCPCS: E2310

Multiple Actuator Control Through 4 Function.....\$2,895
 Push Button and Q-Logic 3 Controller and Interface (AAM)
 Only available with TB3 armrests. Select a side to mount.
COMING SOON. Part: See XRef. HCPCS: E2311

4 Function Push Button Mount Location

Left Side Mount
 Part: PTOINDV3678

Right Side Mount
 Part: PTOINDV3679

22. TRU-BALANCE 3 ATTENDANT CONTROL

Optional. *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$500*
 Part: ELE167862. HCPCS: E2331

23. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

23.1. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Power Articulating Foot Platform \$2,995
Part: PTOINDV3754. HCPCS: E1012

Manually Adjustable Center Mount Foot Platform... \$250
8"-12" extension part: RIG1802138.
Other extensions: RIG1710345

**2. Lower Extension and Footplates**

8"-12" Extension, 5"x6" Footplates
Part: RIG152577

Not available with power articulating platform.

8"-12" Extension, 5"x8" Footplates
Part: RIG152578

Not available with power articulating platform.

10"-16" Extension, 5"x6" Footplates
Part: RIG150206

10"-16" Extension, 5"x8" Footplates
Part: RIG150203

10"-16" Extension, 5"x11" Footplates
Part: RIG150212

13"-19" Extension, 5"x6" Footplates
Part: RIG151505

13"-19" Extension, 5"x8" Footplates
Part: RIG151506

13"-19" Extension, 5"x11" Footplates
Part: RIG151508

Center Mount Foot Platform Compatibility Matrix for TRU-Balance® 3	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
10"-16" Length Extension	!	✓	✓	✓	✓	✓	✓	✓	✓
13"-19" Length Extension	!	!	!	!	!	!	!	!	✓

✓ No interference issues at this seat-to-floor height.

! May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H \$150*
Part: RIGASMB7120018

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$150*
Part: RIGASMB7120019

Small Calf Panel 18"Wx8"H \$150
Part: FRMASMB7110012

Large Calf Panel 24"Wx8"H \$200
Part: RIG123824
Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$50
Part: RIG152517

Width Extension for 11"D Footplates \$50
Part: RIG152518

2" Bell-Shaped Length Extension \$50
Part: RIG152519

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

4" Outside Edge (pair) \$50
Part: RIG152510

Only available with 6"D footplates and 8"D footplate extension.

6" Outside Edge (pair) \$50
Part: RIG152516

Available with 8"D and 11"D footplates and 11"D footplate extension.

5. Additional Accessories

Reversed Dartex Calf Pad Cover (pair) \$150
Small part: RIG133717. Large part: RIG133719

Reversed Dartex Panel Cover (pair) \$150
Small part: RIG134673. Large part: RIG134744

Wheel Bumpers (pair) \$12.50 ea.
Part: WHLASMB1578 x2

Footplate Connector \$15
Part: RIGASMB7120023

23.2. 70° Power Articulating, Swing Away Legrests

4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Articulating Legrest Upper Assembly

P L R 70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 14.5"-18" Extension
Pair Part: FRMASMB9073. Left Part: FRMASMB9074.
Right Part: FRMASMB9075. HCPCS: E1010

P L R 70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 11"-13" Extension
Pair Part: FRMASMB11444. Left Part: FRMASMB11366.
Right Part: FRMASMB11367. HCPCS: E1010

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Angle Adjustable, 11"-13"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	!	!	!	✓

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following page.

23. LEGRESTS CONTINUED**2. Angle Adjustable Footplate with Heel Loop and Wheel Bumper**¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P	L	R	4"Wx6"D Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566.	P	L	R	5"Wx8"D Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.	P	L	R	6"Wx11"D¹ Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570.
P	L	R	5"Wx6"D Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567.	P	L	R	6"Wx8"D¹ Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569.	P	L	R	8"Wx11"D² Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571.

3. Power Articulating, Swing-Away Legrest AccessoriesCalf Pad (pair) 6"Wx6"H..... No Charge
Part: MEC153735Falcon Gel Padded Knee Adductor Buttons (KBE)(pair)\$125
Part: INDPART2242. HCPCS: E0953[^][^]Please verify the appropriate billing code with all third-party payors.**23.3. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

P	L	R	60°-90° Swing-Away Legrest\$81 ea. or N/C* Pair part: RIG120062. Left: FRMASMB17558. Right: FRMASMB17559 *No charge with Power Positioning. Not available with 11"-13" extensions.
P	L	R	60°-90° Swing-Away Legrest, Short\$81 ea. or N/C* Pair part: RIG136481. Left: RIG132607. Right: RIG132609. *No charge with Power Positioning. For use only with 11"-13" lower extensions.

1b. 70° Swing-Away Legrest Upper**Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.**

P	L	R	70° Swing-Away Legrest (pair)\$162 Pair part: FRMASMB7843 Not available with 11"-13" extensions.
P	L	R	70° Tapered, Swing-Away Legrest (pair)\$162 Pair part: FRMASMB8052 Not available with 11"-13" extensions.
P	L	R	70° Tapered, Swing-Away Legrest, Short (pair)\$162 Pair part: FRMASMB7841 For use only with 11"-13" lower extensions.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
11"-13" extension causes no interference at all seat-to-floor heights.									

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2. Lower Extension

P	L	R	11"-13" Length Part: See XRef N/A with Multi-axis.	P	L	R	11"-16.5" Length Part: See XRef N/A with Angle Adjustable.	P	L	R	13"-16.5" Length Part: See XRef N/A with Multi-Axis.	P	L	R	15.5"-19.5" Length Part: See XRef
----------	----------	----------	---	----------	----------	----------	---	----------	----------	----------	---	----------	----------	----------	---

3a. Extruded Footplate(s)

P	L	R	7.25"Wx6"D Extruded Footplate Standard Pair part: FRMASMB7847. Left: FRMASMB7861 Right: FRMASMB7862	P	L	R	8.5"Wx6"D Extruded Footplate Standard Pair part: FRMASMB7848. Left: FRMASMB7863 Right: FRMASMB7864 May cause interference on seat widths less than 18".
----------	----------	----------	--	----------	----------	----------	---

3b. Angle Adjustable Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P	L	R	4"Wx6"D\$85 ea. Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 HCPCS: K0040	P	L	R	6"Wx8"D¹\$85 ea. Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563 HCPCS: K0040
P	L	R	5"Wx6"D\$85 ea. Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 HCPCS: K0040	P	L	R	6"Wx11"D¹\$85 ea. Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564 HCPCS: K0040
P	L	R	5"Wx8"D\$85 ea. Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562 HCPCS: K0040	P	L	R	8"Wx11"D²\$85 ea. Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565 HCPCS: K0040

3c. Multi-Axis Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P	L	R	4"Wx6"D\$247.50 ea. Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552 HCPCS: K0040	P	L	R	6"Wx8"D¹\$247.50 ea. Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555 HCPCS: K0040
P	L	R	5"Wx6"D\$247.50 ea. Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553 HCPCS: K0040	P	L	R	6"Wx11"D¹\$247.50 ea. Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558 HCPCS: K0040
P	L	R	5"Wx8"D\$247.50 ea. Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554 HCPCS: K0040	P	L	R	8"Wx11"D²\$247.50 ea. Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559 HCPCS: K0040

Legrests section continued on following page.

23. LEGRESTS CONTINUED**4. Swing-Away Legrest Accessories**

Heel Loops (pair) \$50
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi-axis footplates.

Wheel Bumper \$12.50 ea.
 Part: See XRef

Gel Neoprene Footrest Sleeve (pair) \$80
 Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 Part: INDPART2823. HCPCS: K0108

23.4. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ P ☐ L ☐ R Elevating Legrest Assembly \$162.50 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2. Lower Extension

☐ P ☐ L ☐ R 11"-14" Length
 Pair: FRMASMB4755
 Left: RIG138026
 Right: RIG138030
 For reinforced only

☐ P ☐ L ☐ R 14"-17" Length
 Pair: FRMASMB4756
 Left: RIG138033
 Right: RIG138035
 For reinforced only

☐ P ☐ L ☐ R 12"-14.5" Length
 Part: See XRef
 For angle adjustable only

☐ P ☐ L ☐ R 14.5"-18" Length
 Part: See XRef
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-14.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	✓	!	✓	✓

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

☐ P ☐ L ☐ R 7.25"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ P ☐ L ☐ R 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s) with Heel Loops

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". *4"x6" and 5"x6" footplates are only available in pairs.

☐ P ☐ L ☐ R 4"Wx6"D \$85 ea.
 Pair parts: FRMASMB8812 x2.
 HCPCS: K0040

☐ P ☐ L ☐ R 6"Wx8"D1 \$85 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040

☐ P ☐ L ☐ R 5"Wx6"D \$85 ea.
 Pair parts: FRMASMB8811 x2.
 HCPCS: K0040

☐ P ☐ L ☐ R 6"Wx11"D1 \$85 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040

☐ P ☐ L ☐ R 5"Wx8"D \$85 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040

☐ P ☐ L ☐ R 8"Wx11"D2 \$85 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040

4. Elevating Legrest Accessories

Heel Loops (pair) \$50
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable footplates.

Wheel Bumper \$12.50 ea.
 Part: See XRef
 Available on angle adjustable footplates.

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$125
 Part: INDPART2242. HCPCS: E0953^

^Please verify the appropriate billing code to be used for this item with all third-party payors.

24. SET LEGREST LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

☐ 7" L (#3470)	☐ 8.5" L (#3473)	☐ 10" L (#3476)	☐ 11.5" L (#3479)	☐ 13" L (#3482)	☐ 14.5" L (#3485)	☐ 16" L (#3488)	☐ 17.5" L (#3491)	☐ 19" L (#3494)	☐ 20.5" L (#3497)	☐ 22" L (#3669)
☐ 7.5" L (#3471)	☐ 9" L (#3474)	☐ 10.5" L (#3477)	☐ 12" L (#3480)	☐ 13.5" L (#3483)	☐ 15" L (#3486)	☐ 16.5" L (#3489)	☐ 18" L (#3492)	☐ 19.5" L (#3495)	☐ 21" L (#3498)	☐ 22.5" L (#3670)
☐ 8" L (#3472)	☐ 9.5" L (#3475)	☐ 11" L (#3478)	☐ 12.5" L (#3481)	☐ 14" L (#3484)	☐ 15.5" L (#3487)	☐ 17" L (#3490)	☐ 18.5" L (#3493)	☐ 20" L (#3496)	☐ 21.5" L (#3600)	

25. QUICK STEALTH CUSHION CONFIGURATION

Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the complete [Stealth Cushion Order Form](#). This section allows a quick selection of a single cushion. The XYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

25.1. Cushion Dimensions**Cushion Width**

14" W 16" W 17" W 18" W 19" W 20" W 22" W

Cushion Depth

14" D 16" D 17" D 18" D 19" D 20" D 22" D


25.2. Cushion Model Selection**SIMPLICITY G CUSHION**

General use with 2-way stretch, breathable cover.

Part series: ST-SIMGXYY

E2601, 16"W to 20" W \$110

E2602, 22"W \$165


**SOLUTION SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SOLSPXYY

E2607, 14"W to 20" W \$375

E2608, 22"W \$430

**TRU-COMFORT 2 SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover.

Part series: ST-TC2SPPXYY

E2607, 16"W to 20" W \$400

E2608, 22"W \$500


**SPECTRUM GEL SPP CUSHION**

Positioning with 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCGSPPXYY

E2607, 14"W to 20" W \$460

E2608, 22"W \$530

**SPECTRUM FOAM SPP CUSHION**

Skin protection and positioning with 4-way stretch, breathable cover.

Part series: ST-SPCFSPXYY

E2607, 16"W to 20" W \$460

E2608, 22"W \$530

**25.3. Additional Fluid Proof Liners**

Fluid Proof Urethane Liner \$50-\$70

Part: URL1709070-XXYY

Not available with Simplicity and TRU-Comfort 2.

Fluid Proof Polyurethane Liner Coming soon

Not available with Simplicity.

Liner Installation

Install Liner on Cushion

26. PROVIDER INSTALLED ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.
☐ TRU-Balance 3 Cup Holder No Charge
Part: ACC125003

For use with TRU-Balance 3 (track style) armrests.

☐ Cup Holder (Tube-style Armrest) No Charge
Part: ACCASMB2634

For use with tube-style armrests.

☐ TRU-Balance 3 Backpack Holder (Pair) (A) \$50
Part: ACC147774

Not available with back cane option.

☐ TRU-Balance 3 Hydration System (B) \$200
Part: ACC125010

☐ TRU-Balance 3 Glove Box (C) \$125
Part: ACC125011

☐ TRU-Balance 3 Transfer Bar, \$200
Quick Release (Pair) (D)
Part: ACC147773. HCPCS: E0705

☐ TRU-Balance 3 Front Mounted \$50
Cell Phone/ Media Holder (E)
Part: ACC155835

☐ Large, Black, Weather Cover \$202
Part: ACCCOVR1021

☐ TRU-Balance 3 Reflector Kit \$50
Part: ACC124670

☐ TRU-Balance 3 Personal Item Hook \$50
Part: ACC125008

Not available with back cane option.


☐ TRU-Balance 3 Rear View Mirror \$50
Part: ACC155833

☐ TRU-Balance 3 Oxygen Holder \$200
Part: ACC125006. HCPCS: E2208

☐ TRU-Balance 3 Clothing Guard (Pair) \$200
Part: ACC125007

☐ Quantum Backpack \$75
Part: ACCINDV1102

 Stealth Large Tablet Holder \$375
Part: ACC7120005

 Stealth Cell Phone/Media Device Holder \$195
Part: ACC7120006
TRU-Balance 3
Accessories

27. BATTERIES**REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.☐ NF-22, Gel Cell, Interceptor Battery..... \$330 ea.
Part: BAT1704613. HCPCS: E2361☐ No batteries. Provider will assume responsibility of batteries.
Part: Configuration selection☐ NF-22, AGM, MK Battery..... \$360 ea.
Part: BAT121346. HCPCS: E2361☐ NF-22, Gel Cell, MK Battery \$360 ea.
Part: BATGEL1008. HCPCS: E2361**28. BRANCH PATCHES****Optional.**Air Force Patch \$20
Part: ACC140240Army Patch \$20
Part: ACC140241Coast Guard Patch \$20
Part: ACC140239Marines Patch..... \$20
Part: ACC140243Navy Patch \$20
Part: ACC140244Purple Heart Patch \$20
Part: ACC140238**29. MANUFACTURER INSTALLED POSITIONING COMPONENTS****Optional. Not available with contoured seating.** More options are available on the [Positioning Components Order Form](#).**29.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, Rail Mounts, pair (TWBTLTC) \$472.50
Part: INDPART2677. HCPCS: E1028**2. Pads**Curved Pads, 3"x5", pair (TWBL-3X5)..... \$210
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6)..... \$210
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**29.2. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide BracketFlip Down Bracket (TWBTLADD-HW) \$236.25 ea.
Part: INDPART2804. HCPCS: E1028Removable Bracket (LP430)..... \$262.50 ea.
Part: INDPART2805. HCPCS: E1028Removable Bracket with 1" Offset \$400 ea.
(LP430 & LHW-121)
Part: INDPART2854. HCPCS: E1028Flip Down with LEFT Full Surface Contact \$236.25 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028Flip Down with RIGHT Full Surface Contact..... \$236.25 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028**2. Pelvic/Thigh Guide Adapter**2" Height Adapter Plate (LHW-122)..... \$26.25
Part: INDPART2821. HCPCS: K0108

Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) \$53.50
Part: POS140430. HCPCS: K0108

Only available with the flip down full surface bracket.


3. Pelvic/Thigh Guide Pad

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4)..... \$75 ea.
Part: INDPART2801. HCPCS: E0953^Medium Pad, 4"x8" (TWBADD-4X8)..... \$75 ea.
Part: INDPART2802. HCPCS: E0953^Long Pad, 4"x12" (TWBADD-4X12)..... \$75 ea.
Part: INDPART2803. HCPCS: E0953^**4. Pelvic/Thigh Guide Installation**Install on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719

30. PROVIDER INSTALLED POSITIONING COMPONENTSOptional. Not available with contoured seating. More options are available on the [Positioning Components Order Form](#).**30.1. Stealth Foot Boxes**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$498.75
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954^Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$498.75
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954^Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954^**30.2. Stealth Padded Shoe Holder**Stealth Small Padded Shoe Holder 3.75"Wx6.25"D..... \$266
(SPLE111) (pair)
Part: INDPART2783Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D..... \$266
(SPLE114) (pair)
Part: INDPART2785Stealth Medium Padded Shoe Holder, 4"Wx8.25"D..... \$266
(SPLE112) (pair)
Part: INDPART2784**30.3. Stealth Elbow Support****1. Elbow Support Bracket**Left, Swing-Away, Elbow Support \$183.75 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028Right, Swing-Away, Elbow Support \$183.75 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028Left, Fixed Elbow Support Bracket..... \$183.75 ea.
TRU-Balance 3 Arms (ARMS-ESFTB3-L)
Part: POS133465Right, Fixed Elbow Support Bracket \$183.75 ea.
TRU-Balance 3 Arms (ARMS-ESFTB3-R)
Part: POS133466**2. Elbow Support Pad**3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$52.50 ea.
Part: POS124065. HCPCS: K01084"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$52.50 ea.
Part: POS124067. HCPCS: K0108**30.4. Stealth i-Fit Chest Harness**Upper Classic, Medium (IFPS1103) \$142 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960Upper Contour, Medium (IFPS2103)..... \$142 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960Upper Classic, Large (IFPS1104) \$142 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960Upper Contour, Large (IFPS2104) \$142 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960**30.5. Stealth i-Fit Chest Straps**Chest Strap, Medium w/ Belt Anchor \$90 ea.
(IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960Chest Strap, Large w/ Belt Anchor \$90 ea.
(IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960**30.6. Stealth Swing-Away Headrest Mounts**


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

Facial Lateral for Comfort Series.....\$294.25 ea.
Right Side (CPS9-R-HW)
Part: INDPART2647. HCPCS: E1028Facial Laterals for Comfort Series \$347.75 ea.
Left and Right Side (CPS9-B-HW)
Part: INDPART2649. HCPCS: E1028Facial Lateral for Comfort Series.....\$294.25 ea.
Left Side (CPS9-L-HW)
Part: INDPART2648. HCPCS: E1028**Facial Lateral Pads and Accessories**1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$80.25
Part: POS150532Facial Pad, 4"x2" (900-Q)..... \$53.50 ea.
Part: ST-900-Q. HCPCS: K0108

Required with facial laterals to mount to comfort headrest brackets.

Facial Pad, 6"x2" (920-Q)..... \$53.50 ea.
Part: ST-920-Q. HCPCS: K0108Egg Switch Mount (ESM) (A) \$52.50 ea.
Part: INDPART2815Swing-Away Joystick Mount, 18" Rod (C) \$250 ea.
Part: ST-SUS18-QFigure 8 Bracket with 13" Rod..... \$267.50 ea.
and Egg Switch Mount (B)
Part: ST-F8675-13-Q

Allows mounting of facial pads and attaches to the headrest



31. SPECIALTY CONTROLS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**31.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,720
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard **Match Joystick Shroud to Base Color** \$50
Part: DGN137728 Part: See color below**Non-matching Joystick Shroud**..... \$50
Select color below**Red Corvette**
Part: DGN137729-PRD1**White Satin**
Part: DGN137729-PWH2**Electric Blue**
Part: DGN137729-PBU9**Sugar Plum**
Part: DGN137729-PPR7**Cherry Bomb**
Part: DGN137729-PRD4**Deep Purple**
Part: DGN137729-PPR5**Blue Jean**
Part: DGN137729-PBU11**Military Digital Camo**
Part: DGN137729-WCM7**Back in Black**
Part: DGN137729-PBK1**Raspberry Beret**
Part: DGN137729-PPK4**Lemon Crush**
Part: DGN137729-PYL4**Faded Glory**
Part: DGN137729-WFL2**Night Train**
Part: DGN137729-PBK3**Evolution Orange**
Part: DGN137729-POR7**Green Machine**
Part: DGN137729-PGN11**Carbon Fiber**
Part: DGN137729-WCB1**Grey Street**
Part: DGN137729-PSL3**Midnight Blue**
Part: DGN137729-PBU10**Military Green**
Part: DGN137729-PGN8**To Be Determined**
Follow-up required.

Enhanced Display

31.2. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ**.....\$700
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick..



EX SAJ

31.3. Head Array**Adult, i-Drive Tri-Array Series, 3-Switch** \$4,414
Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q, HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461
(IDHBT300-1)

Part: ST-IDHBT300-1-Q, HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

31.4. Sip and Puff Drive**Module****Sip and Puff Module for Q-Logic 3**.....\$1,895
Part: CTL168056

Includes module and mounting hardware.



Sip & Puff Module

Straw and Mount**Stealth Sip and Puff for Comfort Series, Left Mount (SP625-CL)**\$294.25

Part: ST-SP625-CL-Q

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U)\$294.25

Part: ST-SP625-U-Q

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, Right Mount (SP622-CR)\$294.25

Part: ST-SP625-CR-Q

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip & Puff

Specialty controls section continued on following page.

31. SPECIALTY CONTROLS CONTINUED**31.5. Egg Switches**

Red Egg Switch (ESRED)..... \$80.85
 QTY Part: SWTMCRO1077
 Blue Egg Switch (ESBLU)..... \$80.85
 QTY Part: SWTMCRO1079
 Black Egg Switch (ESBLK)..... \$80.85
 QTY Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$80.85
 QTY Part: SWTMCRO1078
 Yellow Egg Switch (ESYEL) \$80.85
 QTY Part: SWTMCRO1080

**31.6. Switch Adapter**

2 to 1 Switch Adapter..... \$13.39
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

32. TRU-BALANCE 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA)..... No Charge
 Part: Configuration selection

Secondary Multiplier (SM)..... \$221.21
 Part: HAR144004

Special Order Configuration (SOF)
 Part: Configuration selection

		Power Positioning Chart			
		1 Arm	2 Arm	3 Arm	4+ Arm
Arm Mounts	Joystick (Hand Control)	NA	SM	SM	SOF
	Enhanced Display	NA	SM	SM	SOF
	iAccess	NA	SM	SM	SOF
	USB Charger	SM	SM	SM	SOF
Back Mounts	SAJ	SOF	SOF	SOF	SOF
	Attendant Control	NA	SM	SOF	SOF
	SCIM	NA	SM	SOF	SOF
	Sip n Puff	SM	SM	SOF	SOF
	ECU	SM	SM	SOF	SOF
	NA	SOF	SOF	SOF	SOF
	SM	SOF	SOF	SOF	SOF
	SOF	SOF	SOF	SOF	SOF

Static Seat Chart

Q-Logic 3 Multiplier Charts

33. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2018 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.

