

New!

Get Out There!

# GO-CHAIR

GET OUT THERE IN STYLE WITH THE ALL-NEW GO-CHAIR



NEW, dual under-seat compartments for convenient storage

Shown in: Ruby Red



The all-new Go-Chair is re-engineered from the ground up, offering a sleek, bold new style available in an array of contemporary colors. Enhanced performance and comfort, along with feather-touch disassembly, allows you to enjoy light-weight travel and independence on the go. With an increased weight capacity of 300 lbs. and a sleek, bold look, the new Go-Chair makes travel easy.



One-hand "Feather Touch" disassembly

Sapphire Blue    Rose Quartz    Citrine Yellow    Amber Orange    Pearl White



Variety of color options

# SPECIFICATIONS

## GO CHAIR

Drive Wheels	9" (solid)
Front Casters	6" (solid)
Rear Anti-tips	3" (solid)
Maximum Speed <sup>1</sup>	Up to 3.7 mph
Ground Clearance <sup>2</sup>	1.6" at anti-tip bracket 1.9" at frame
Turning Radius <sup>2</sup>	25.4" with foot platform
Length	33" with front rigging
Width <sup>6</sup>	22.6"
Seat Type	Foldable, compact, lightweight
Batteries <sup>3,5</sup>	(2) 18 Ah, 12 V, deep-cycle
Range <sup>1,4</sup>	Up to 8.74 miles
Charger	Off-board, 2 amp
Electronics	40 amp, PG GC 2 (non-programmable)
Weight Capacity	300 lbs.
Base Weight	67 lbs. (front & rear)
Standard Seat Weight	30 lbs.
Battery Pack Weight <sup>5</sup>	31 lbs.
Suspension	Limited
Heaviest Piece	36 lbs. (rear)

- <sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition tire type and condition. These specifications can be subject to a variance of +10%, -5%.  
<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.  
<sup>3</sup> AGM or Gel-Cell type required. \*See "Batteries and Charging."  
<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.  
<sup>5</sup> Battery weight may vary based on manufacturer.  
<sup>6</sup> Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

**NOTE: This product conforms to all applicable ANSI/RESNA and ISO 7176 series standards. This product is not coded by DME PDAC or does not meet criteria: K0899. All specifications subject to change without notice.**

## COMPATIBLE LIFTS

- Outlander Full Platform
- Olympian
- Backpacker Series
- Commander

To find a compatible lift for your vehicle/mobility product, use our new tool:  
[pridemobility.com/lift-finder/](http://pridemobility.com/lift-finder/)

## 4-wheel

### FEATURES

- Choice of 6 colors
- 18" x 17", 5-point swivel seat with latch release
- Increased 300 lb. weight capacity
- Faster speed - up to 3.7 mph
- Larger, 18 amp batteries
- Improved agility - longer wheel base; wider stance; compact turning radius 25.4"
- Dual swing-away, under seat storage bins
- Larger, height & width adjustable armrests
- Larger footplate
- One hand feather-touch disassembly

### OPTIONS

- Cup Holder
- Lap Belt- 50", 60", 70"
- Safety Flag
- Oxygen Holder
- Dual Saddle Bag
- Weather Cover, Micro
- Additional Battery Pack
- Rear Mount Bracket Cane/Crutch Holder Single (2", 3"), Double (2", 3")
- Armrest Mount Cane Clip

\*Military Patches Available



## COMPATIBLE LIFTS

- Outlander Full Platform
- Olympian
- Backpacker Series
- Commander

To find a compatible lift for your vehicle/mobility product, use our new tool:  
[pridemobility.com/lift-finder/](http://pridemobility.com/lift-finder/)

**Pride**  
 Mobility Products Corp.  
 — LIVE YOUR BEST —

**GO GO**  
 TRAVEL MOBILITY

US: 800-800-8586 · CAN: 888-570-1113

▶ /PrideMobility    📌 /PrideMobilityProductsCorp

[pridemobility.com](http://pridemobility.com)